

## Product datasheet for **SC318880**

### Endothelin B Receptor (EDNRB) (NM\_001122659) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Endothelin B Receptor (EDNRB) (NM_001122659) Human Untagged Clone
Tag:	Tag Free
Symbol:	EDNRB
Synonyms:	ABCD5; ET-B; ET-BR; ETB; ETB1; ETBR; ETRB; HSCR; HSCR2; WS4A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC318880 representing NM_001122659. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGCCGCTCCAAGTCTGTGCGGACGCGCCCTGGTTGCGCTGGTTCTTGCCTGCGGCCTGTGCGGG
ATCTGGGGAGAGGAGAGAGGCTTCCCGCTGACAGGGCCACTCCGCTTTTGCAAACCGCAGAGATAATG
ACGCCACCCACTAAGACCTTATGGCCCAAGGGTCCAACGCCAGTCTGGCGGGTCTGTTGGCACCTGCG
GAGGTGCCTAAGGAGACAGGACGGCAGGATCTCCGCCACGCACCATCTCCCTCCCCCGTCCAAGGA
CCCATCGAGATCAAGGAGACTTTCAAATACATCAACACGGTTGTGCTCCTGCCTTGTGTTCTGCTGGGG
ATCATCGGGAACCCACACTTCTGAGAATTATCTACAAGAACAAGTGCATGCGAAACGGTCCCAATATC
TTGATCGCCAGCTTGGCTCTGGGAGACCTGCTGCACATCGTCATTGACATCCCTATCAATGTCTACAAG
CTGCTGGCAGAGGACTGGCCATTTGGAGCTGAGATGTGTAAGCTGGTGCCTTTCATACAGAAAGCCTCC
GTGGGAATCACTGTGCTGAGTCTATGTGCTCTGAGTATTGACAGATATCGAGCTGTTGCTTCTTGGAGT
AGAATTAAGGAATTGGGTTCCAAAATGGACAGCAGTAGAAAATTGTTTTGATTTGGGTGGTCTCTGTG
GTTCTGGCTGTCCCTGAAGCCATAGGTTTTGATATAATTACGATGGACTACAAAGGAAGTTATCTGCGA
ATCTGCTTGCCTCATCCCGTTCAGAAGACAGCTTTCATGCAGTTTTACAAGACAGCAAAAGATTGGTGG
CTATTCAGTTTCTATTTCTGCTTGCCATTGGCCATCACTGCATTTTTTTATACACTAATGACCTGTGAA
ATGTTGAGAAAAGAAAAGTGGCATGCAGATTGCTTTAAATGATCACCTAAAGCAGAGACGGGAAGTGGCC
AAAACCGTCTTTTGCCTGGTCTTGTCTTTGCCCTCTGCTGGCTTCCCTTACCTCAGCAGGATTCTG
AAGCTCACTCTTTATAATCAGAATGATCCAATAGATGTGAACCTTTTGGCTTTCTGTTGGTATTGGAC
TATATTGGTATCAACATGGCTTCACTGAATTCCTGCATTAACCCAATTGCTCTGTATTTGGTGAGCAAA
AGATTCAAAACTGCTTTAAGTCATGCTTATGCTGCTGGTGGCAGTCATTTGAAGAAAAACAGTCTTG
GAGGAAAAGCAGTCGTGCTTAAAGTTCAAAGCTAATGATCACGGATATGACAACCTCCGTTCCAGTAAT
AAATACAGCTCATCTGA
ACGGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001122659

Insert Size: 1329 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001122659.2](#)

RefSeq Size: 4320 bp

RefSeq ORF: 1329 bp

Locus ID: 1910

UniProt ID: [P24530](#)

<b>Cytogenetics:</b>	13q22.3
<b>Protein Families:</b>	Druggable Genome, GPCR, Transmembrane
<b>Protein Pathways:</b>	Calcium signaling pathway, Melanogenesis, Neuroactive ligand-receptor interaction
<b>MW:</b>	49.6 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a G protein-coupled receptor which activates a phosphatidylinositol-calcium second messenger system. Its ligand, endothelin, consists of a family of three potent vasoactive peptides: ET1, ET2, and ET3. Studies suggest that the multigenic disorder, Hirschsprung disease type 2, is due to mutations in the endothelin receptor type B gene. Alternative splicing and the use of alternative promoters results in multiple transcript variants. [provided by RefSeq, Oct 2016]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1 and 3 both encode isoform 1.</p>