

Product datasheet for **SC318888**

alpha 2 Glycine Receptor (GLRA2) (NM_001118886) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	alpha 2 Glycine Receptor (GLRA2) (NM_001118886) Human Untagged Clone
Tag:	Tag Free
Symbol:	GLRA2
Synonyms:	GLR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC318888 representing NM_001118886. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACCGGCAGCTAGTGAACATTTTGACAGCCTTGTTTGCATTTTCTTAGAGACAAACCACTTCAGG
ACGGCTTTCTGCAAAGACCATGACTCCAGGTCTGGAAAACAACCTTCACAGACCCATCTCCTTCAGAT
TTCTTGGACAAGTTAATGGGAAGGACATCAGGATATGATGCAAGAATCAGGCCAAATTTAAAGGGCCT
CCTGTAATGTTACCTGCAACATATTTATCAACAGCTTTGGGTCAATAGCAGAACTACAATGGACTAC
CGAGTGAATATTTTCTGAGACAACAGTGAATGATTCACGGCTGGCGTACAGTGAGTACCCAGATGAC
TCCCTGGACTTGGACCCATCCATGCTAGACTCCATTTGAAAACCCAGATTTGTTCTTTGCCAATGAGAAG
GGTGCCAACTTCCACGATGTACCACCTGACAACAAATTGCTACGGATTTGAAAAATGGCAAAGTGCTC
TACAGTATCAGACTCACCTTGACCTTATCCTGTCCCATGGACTTGAAGAATTTCCGATGGATGTCCAG
ACCTGTACAATGCAGCTGGAGAGTTTGGGTACACGATGAATGACCTGATATTTGAGTGGTTAAGTGAT
GGTCCAGTGCAAGTTGCTGAAGGATTGACCCTGCCCCAGTTTATTTTGAAGAAGAGAAGGAACTTGGC
TACTGTACAAAGCACTACAACACTGAAAGTTTACCTGCATTGAGGTCAAGTTTCATCTGGAACGCCAA
ATGGGATATTATTTGATCCAGATGTACATCCCAAGCCTGCTTATAGTAATTTTGTCTGGGTTTCCTTT
TGGATAAATATGGATGCAGCCCCTGCCAGGGTCCGACTGGGCATCACCACAGTCTTAACGATGACCACC
CAGAGTTCAGGCTCCAGGGCATCTGCCAAAGTCTCCTATGTAAGGCGATTGACATCTGGATGGCG
GTGTGCCTTCTGTTTGTGTTGCTGCCTTACTGGAATACGACGGTGAACCTTCGCTCCAGGCAACAC
AAGGAGTTCTGCGCCTCCGAAGAAGACAGAAGAGGCAGAATAAGGAAGAAGACGTTACTCGTGAAGT
CGTTTTAATTTAGCGTTATGGGATGGGTCACTGCCTCCAAGTGAAGATGGAACAGCTGTCAAGGCC
ACACCTGCCAACCCACTCCACAACCGCCAAAAGATGGAGATGCTATCAAGAAGAAGTTTGTGGACCGG
GCAAAAAGGATTGACACGATATCTCGAGCTGCCTTCCATTGGCCTTCTCATTTTCAACATCTTTTAC
TGGATCACATAAAGATCATTCCGCATGAAGATGTCCACAAGAAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001118886

Insert Size: 1359 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001118886.1](#)

RefSeq Size: 3458 bp

RefSeq ORF: 1359 bp

Locus ID: 2742

UniProt ID: [P23416](#)

Cytogenetics:	Xp22.2
Protein Families:	Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction
MW:	52 kDa
Gene Summary:	<p>The glycine receptor consists of two subunits, alpha and beta, and acts as a pentamer. The protein encoded by this gene is an alpha subunit and can bind strychnine. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2010]</p> <p>Transcript Variant: This variant (3) contains an alternate in-frame exon and lacks an exon, compared to variant 1, resulting in a protein (isoform B) of identical length but differing in 2 amino acids, compared to isoform A.</p>