

## Product datasheet for **SC318890**

### MANEL (MANEAL) (NM\_001113482) Human Untagged Clone


#### Product data:

Product Type:	Expression Plasmids
Product Name:	MANEL (MANEAL) (NM_001113482) Human Untagged Clone
Tag:	Tag Free
Symbol:	MANEAL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC318890 representing NM_001113482. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCGGCGGGCGCGCCGCGCTGCATCGCTCTGTTCTGGTCTGCTCTTCGCCTTCGGCACCTC
ATGGGTCTGCGCACGCTCAAGGCTCCGGACGACTCCCGCGCTGGGCCCGGGCTGGAGCTGGCGCCC
TTTGAGCGACGCCAGAGGGGGCCCGCGCCGCTGCCAGGGCCCGGACGCCCTGCCGCGCCGCC
CCGCGCCCGCGCCGCCCGCCGACCGCGGACCCTGGCGGCTCCCTGGGCCGGCACCCGCGGAGGCCGAG
CCCGCCCCGTGCAGAGCCTGCGCGTCTACTCGGACCTGCACGCCTTCTACTACTCGTGGTACGGGAGC
CCGCGCGCGAGGGCCACTACATTCACCTGGGACCAGTCACTGGTCCGCACTGGGACCCCAAGATCTCG
GCCAGTACCCCGCGGCCGCCACAGCCCCAGACACTTGGGCTCCAGCTTCTACCCGGAGCTGGGG
CCCTACAGCTCCCGGACCCCGAAGTGCTGCGGGAGCATATGACCCAGCTCAAGGAAGCCGCCATCGGC
GTCCTGGTCTGTCTGTACCCACCTGGCATGGCTGATGATAACGGGGAGCCCTCAGATGACCTGGTG
CCCGCCATTCTGGACACCGCCATCAGTACAGCATCCAGGTGGCCTTCCACATCCAACCTACAAGGGC
CGGGATGACACTGTACATGACAACATCAAGTACATCATTGACACGTATGGCTCCCATGGTGCATTT
TACCGCTATAAAGACAGCATGGCAAGAGCCTCCCACTCTTTATATCTACGACTCATACCTGACGTCC
CCTGAGGCTGGGCCACCTCCTGACACCAACGGGCCCATTCGATCCGCAACACGCCCTACGATGGG
GTCTTCATAGCGTGTGGTGGAGAGGGCCACCCACGATATCCTGGCCGCGGATTTGACGGCATG
TACACCTACTTTGCCTCCAATGGTTTCTCTTTGGTTCTCCCATCAGAAGCTGGAAGAGCTGAAGAAC
TTTTGTGATGCCAACAACCTCATGTTATCCCCAGTGTGGGGCCTGGCTACATAGACACCAGCATTCCG
CCCTGGAACAACCACAATACGCGCAACAGGGTCAATGGCAAGTACTATGAGACGGCCCTGCAGGCGGCC
CTGACAGTGGAGCCCGAGATCGTTTCCATTACCTCCTTCAATGAGTGGCAGGAGGCCACCCAGATTGAG
AAGGCCATTCCCAAGAAGACACCCACCGCCTGTATTTGGACTACCTGCCTACCCAGCCAGCCTGTAC
CTGGAGCTGACACGCGCTGGCGGAGCACTTCATCAAAGAGAAGGAGCAGTGGCTCATGTA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



**Plasmid Map:**


**ACCN:** NM\_001113482

**Insert Size:** 1374 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001113482.1](#)

**RefSeq Size:** 2858 bp

**RefSeq ORF:** 1374 bp

**Locus ID:** 149175

**UniProt ID:** [Q5VSG8](#)

**Cytogenetics:** 1p34.3

**Protein Families:** Transmembrane

**MW:** 51.3 kDa