

Product datasheet for **SC318962**

SLC44A3 (NM_001114106) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC44A3 (NM_001114106) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC44A3
Synonyms:	CTL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318962 representing NM_001114106.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCACTGCCTGGGCGCCGAGTACCTGGTTTCTGCAGAAGGAGCCCTAGGCAAAGGGAGTGGCGACCC
CAGATTTATAGGAAATGCACAGATACGGCATGGTTATTCTGTTCTTTCTTTTGGACTGGTTGGTG
TTTATCATGGGCTACTCGGTGGTGGCTGGAGCCGCGGAAGACTCCTCTTTGGCTATGACAGCTTTGGC
AACATGTGTGGCAAGAAGAAGTCCCCCGTGAAGGGGCCCTCTTTCAGGGCAGGACATGACCCTAAAA
AAACACGTGTTCTTTATGAATCCTGCAACCTGGAAGTCAAAGGTACGCAGCTCAACCGCATGGCCCTC
TGTGTATCCAACCTGCCCTGAAGAGCAGCTTGACTCCCTGGAAGAGGTCCAGTTCTTTGCAAAACCCAGT
GGGTCTTCTGTGTGTTATAGTTTGAATTCCTCAACTATACCCACAGTCCAAAAGCAGACTCACTG
TGTCCCAGGCTACCAGTTCCTCCAAGCAAGTCATTTCCCTATTTAACCGATGTGTCCTCAAACACCT
GAGTGCTACTCCCTATTTGCATCTGTTTTGATAAATGATGTTGACACCCTCCACCGAATTCTAAGTGA
ATCATGTGCGGAAGAGATACAATCCTTGGCCTGTGTATCCTCGCATTAGCCTTGCTTTGGCCATGATG
TTTACCTTCAGATTCATACCACCTTCTGGTTACATTTTCATTTTCATTGGTTATTTGGGATTGTTG
TTTGTCTGCGGTGTTTTATGGTGGCTGATTATGACTATACCAACGACCTCAGCATAGAATTGGACACA
GAAAGGGAAAATGAAGTGCCTGCTGGGGTTGCTATCGTATCCACAGGCATCACGGCAGTGTGCTC
GTCTTGATTTTTGTTCTCAGAAAGAGAATAAAATTGACAGTTGAGCTTTTCCAAATCACAAATAAAGCC
ATCAGCAGTGTCCCTTCTGCTGTTCCAGCCACTGTGGACATTTGCCATCCTCATTCTTCTGGGTG
CTCTGGGTGGCTGTGCTGCTGAGCCTGGAACTGCAGGAGCTGCCAGGTTATGGAAGGCGGCCAAGTG
GAATATAAGCCCCCTTCGGGCATTCGGTACATGTGGTACGTTACCATTTAATTGGCCTCATCTGGACTAGT
GAATTCATCCTTGCCTGCCAGCAATGACTATAGCTGGGCAGTGGTTACTTGTTATTTCAACAGAAGT
AAAAATGATCCTCTGATCATCCCATCCTTTCTGCTCTCTCCATTCTTCTTCTACCATCAAGGAACC
GTTGTGAAAGGGTCATTTTTAATCTCTGTGGTGAAGATTCCGAGAATCATTGTCATGTACATGCAAAAC
GCACTGAAAGAACAGCAGCATGGTGCATTGTCCAGGTACCTGTTCCGATGCTGCTACTGCTGTTTCTGG
TGTCTTGACAAATACCTGCTCCATCTCAACCAGAAATGCATATACTACAACCTGCTATTAATGGGACAGAT
TTCTGTACATCAGCAAAGATGCATTCAAATCTTGTCCAAGAACTCAAGTCACTTTACATCTATTAAC
TGCTTTGGAGACTTCATAATTTTTCTAGGAAAGGTGTTAGTGGTGTGTTTCACTGTTTTTGGAGGACTC
ATGGCTTTTAACTACAATCGGGCATTCCAGGTGTGGGCAGTCCCTCTGTTATTGGTAGCTTTTTTGGCC
TACTTAGTAGCCATAGTTTTTATCTGTGTTTGAAGTGTGCTGGATGCACTTTTCTGTGTTTTGCT
GTTGATCTGGAACAAATGATGGATCGTCAGAAAAGCCCTACTTTATGGATCAAGAAATTTCTGAGTTTC
GTAAAAAGGAGCAACAAATTAACAATGCAAGGGCACAGCAGGACAAGCACTCATTAAAGGAATGAGGAG
GGAACAGAACTCCAGGCCATTGTGAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_001114106

Insert Size: 1962 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001114106.2</u>
RefSeq Size:	2397 bp
RefSeq ORF:	1962 bp
Locus ID:	126969
UniProt ID:	<u>Q8N4M1</u>
Cytogenetics:	1p21.3
Protein Families:	Transmembrane
MW:	73.8 kDa