

Product datasheet for **SC318970**

Lebercilin (LCA5) (NM_001122769) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lebercilin (LCA5) (NM_001122769) Human Untagged Clone
Tag:	Tag Free
Symbol:	LCA5
Synonyms:	C6orf152
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318970 representing NM_001122769.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGGAAAGAGCAGGAAGTCCAGTACTGATCAAGAAAGAAAGGCAGGCAAACACCATTATTCTTAC
TTATCTGATTTTGAACGCCACAGTCTTCTGGCCGATCATCGCTGGTCAGTTCTTACCTGCAAGTGTT
AGGAGAAAAAATCCTAAAAGACAAACTTCAGATGGCCAAGTACATCACCAGCCCTCGAAACCAAGC
CCTAAGGGTCTACCAAACAGAAAGGGAGTCCGAGTGGGATTTTCGCTCCCAGAGCCTCAATAGAGAGCCA
CTTCGAAAGATACTGATCTTGTACAAAACGGATTCTGTCTGCAAGACTGCTAAAAATCAATGAGTTG
CAGAATGAAGTATCTGAACTCCAGGTCAAGTTAGCTGAGCTGCTAAAAGAAAATAAATCTTTGAAAAGG
CTTCAGTACAGACAGGAGAAAGCCCTGAATAAGTTTGAAGATGCCGAAAATGAAATCTCACAACTATA
TTTCGTCATAACAATGAGATTACAGCACTCAAAGAACGCTTAAGAAAATCTCAAGAGAAAGAACGGGCA
ACTGAGAAAAGGGTAAAAGATACAGAAAGTGAATTTTAGGACAAAATTTTCTTACAGAACTGAAA
GAGATCTCTGAAGCTAGACACCTACCTGAACGAGATGATTTGGCAAAGAACTAGTTTCAGCAGAGTTA
AAGTTAGATGACACCGAGAGAAGAAATTAAGGAGCTATCGAAAAACCTTGAAGTACTAACAGTTTC
CAACGACAGTTGCTTGTGAAAGGAAAAGGGCATATGAGGCTCATGATGAAAATAAAGTTCTTCAAAG
GAGGTACAGCGACTATATCACAATTAAGGAAAAGGAGAGAGAACTGGATATAAAAAATATATATTCT
AATCGTCTGCCAAAGTCTCTCCAAATAAAGAGAAAGAACTTGCATTAAGAAAAATGCTGCATGCCAG
AGTGATTTTGCAGACCTGTGTACAAAAGGAGTACAACCATGGAAGACTCAAGCCAGAAGAATATCCT
TAACTCCAGAAACAATTATGTGTACGAAAACAAATGGGAAGAACCAGGACATCTACTTTGGACTTG
CAATCTCAAAGCAAGACAGGCATGGAGAAGCAGGGATTCTAAACCAATTATGGAAAGAGAAGAAAAA
TTTGTACAGATGAAGAAGTCCATGTCGTAACAGGAGGTTGAAAAGCTGGAGGATGAATGGGAAAGA
GAAGAATTGATAAAAAGCAAAAAGAAAAGGCATCTTTACTGGAAAGAGAAGAAAAGCCAGAGTGGGAA
ACTGGAAGGTACCAACTAGGAATGTATCCAATTCAGAATATGGATAAATTGCAAGGAGAGGAAGAAGAA
AGACTGAAGAGAGAAATGCTACTTGCTAAACTGAATGAAATTGACAGAGAACTCCAAGTTCTCGAAAT
CTAAAATACCCTGTTTTGCCATTGTTACCTGATTTTGAATCAAACTACACTCCCAGAGAGAAGCCCC
AAACATACAGTTCTCTGAATCCTCAGAGAGATTATTTAATGGGCATCATTTGCAAGACATCAGTTTC
TCAACTCCAAAAGGAGAAGGTGAGAATTCAGGAAATGTTAGAAGTCCAGCCTCCCTAATGAGTTCGCA
TTTGGTAGCTACGTGCCTTCGTTTGCAAAACATCAGAGAGGTCAAATCCATTTAGTCAAAAAGTAGT
TTTTTGGATTTCAAAGAACAGTATGGAAAACCTTAGTAAAGATGGTGTAGATTTAATACAAGAAAA
GAGAAAAAAGCTAATTTGATGGAACAGTTATTTGGTGCCAGTGGTAGCAGCACCATTTCTCCAAAAGC
AGTGACCCAAATCTGTGGCTTCCAGTAAAGGAGACATTGACCCTCTAAATTTTCTCCCTGGGAATAAA
GGCAGCAGAGATCAAGAACATGATGAAGATGAAGGCTTTTTCTCAGTGAAGGAAGAAGTTTAAATCCA
AATAGGCCACCGATTAACATGCAGACGATAAACCAGCAGTAAAAGCAGCTGATTCTGTAGAAGATGAA
ATTGAAGAAGTAGCACTGAGATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:

ACCN: NM_001122769

Insert Size: 2094 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001122769.2</u>
RefSeq Size:	4574 bp
RefSeq ORF:	2094 bp
Locus ID:	167691
UniProt ID:	<u>Q86VQ0</u>
Cytogenetics:	6q14.1
MW:	80.6 kDa
Gene Summary:	<p>This gene encodes a protein that is thought to be involved in centrosomal or ciliary functions. Mutations in this gene cause Leber congenital amaurosis type V. Alternatively spliced transcript variants are described. [provided by RefSeq, Oct 2009]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein.</p>