

## Product datasheet for **SC318978**

### **EPLIN (LIMA1) (NM\_001113546) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	EPLIN (LIMA1) (NM_001113546) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPLIN
Synonyms:	EPLIN; LDLCQ8; SREBP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC318978 representing NM\_001113546.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAATCATCTCCATTTAATAGACGGCAATGGACCTCACTATCATTGAGGGTAACAGCCAAAGAAGCTT
TCTCTTGTCAACAAGAACAAGTCACTCGGCTATTGTGAAATATTCTCCAAGTACCAGAAAGCAGCTGAA
GAAACAAACATGGAGAAGAAGAGAAGTAACACCGAAAATCTCTCCAGCACTTTAGAAAGGGGACCCCTG
ACTGTGTTAAAGAAGAAGTGGGAGAACCAGGGCTGGGAGCAGAGTCTCACACAGACTCTCTACGGAAC
AGCAGCACTGAGATTAGGCACAGAGCAGACCATCCTCCTGCTGAAGTGACAAGCCACGCTGCTTCTGGA
GCCAAAGCTGACCAAGAAGAACAATCCACCCAGATCTAGACTCAGGTCACCTCCTGAAGCCCTCGTT
CAGGGTCGATATCCACATCAAGGACGGTGAAGTCTTAAAGACCACTCAACAGAAAGTAAAAAATG
GAAAATGTCTAGGAGAATCCAGGCATGAAGTAGAAAAATCAGAAATCAGTGAAAACACAGATGCTTCG
GGCAAAATAGAGAATATAATGTTCCGCTGAACAGGCTTAAGATGATGTTTGAGAAAGGTGAACCAACT
CAAATAAGATTCTCCGGCCCCAAGCCGAAGTGAAGTGAAGGAAGATCTCTGAAAACAGCTATTCT
CTAGATGACCTGGAAATAGGCCACAGTCACTTGTATCTTCTACATTTGACTCGGAGAAAAATGAGAGT
AGACGAAATCTGGAACCTCCACGCCTCTCAGAAACCTCTATAAAGGATCGAATGGCCAAGTACCAGGCA
GCTGTGTCAAACAAAGCAGCTCAACCAACTATACAAATGAGCTGAAAGCCAGTGGTGGCGAAATCAAA
ATTCATAAAATGGAGCAAAGGAGAATGTGCCCCAGGCTCCTGAGGTCTGCATCACCCATCAGGAAGGG
GAAAAGATTTCTGCAATGAGAAATAGCCTGGCAGTCCGTTCCACCCCTGCCGAAGATGACTCCCCAGGT
GACTCCCAGGTTAAGAGTGAAGTCAACAGCCTGTCCATCCCAAGCCACTAAGTCCAGATTCCAGAGCC
TCCAGTCTTTCTGAAAGTCTCCTCCAAAGCAATGAAGAAGTTTCAGGCACCTGCAAGAGAGACCTGC
GTGGAATGTGAGAAGACAGTCTATCCAATGGAGCGTCTCTGGCCAACCAGCAGGTGTTTACATCAGC
TGCTTCCGTTGCTCCTATTGCAACAACAACACTCAGTCTAGGAACATATGCATCTTTACATGGAGAATC
TATTGTAAGCCTCACTTCAACTCTTTAAATCTAAGGGCAACTATGATGAAGCTTTGGGCACAGA
CCACACAAGGATCTATGGGCAAGCAAAAATGAAAACGAAGAGATTTGGAGAGACCAGCCAGCTTGCA
AATGCAAGGGAGACCCCTCACAGCCAGGGGTAGAAGTGCCCTATTGCTAAGGTGGGTGCTCTGGCT
GCAAGTATGGAAGCCAAGGCTCCTCTCAGCAGGAGAAGGAAGACAAGCCAGCTGAAACCAAGAAGCTG
AGGATCGCTGGCACCCTCACTGAAGTTCAGGAAGTGCCTTGAGGAAAGGGATCAAAATG
TCAAAGCCAAAATGGCTCCTGAAGACGAAATCAGCAAGCCGAAGTTCCTGAGGATGTCGATCTAGAT
CTGAAGAAGCTAAGACGATCTTCTTCACTGAAGGAAAGAAGCCGCCATTCACTGTAGCAGCTTCATTT
CAAAGCACCTCTGTCAAGAGCCAAAACCTGTGTCCCACTATCAGGAAAGGCTGGAGCATGTCAGAG
CAGAGTGAAGAGTCTGTGGTGGAAAGATTGCAGAAAGGAAACAAGTGGAAAATGCCAAGGCTTCTAAG
AAGAATGGGAATGTGGGAAAAACAACCTGGCAAAACAAGAATCTAAAGGAGAGACAGGGAAGAGAAGT
AAGGAAGTCAATAGTTTGGAGATGGAGAATGAGAATCTTGTAGAAAATGGTGCAGACTCCGATGAAGAT
GATAACAGCTTCTCAACAACAATCTCCACAAGAACCAAGTCTCTGAATTGGTCGAGTTTTGTAGAC
AACACCTTTGCTGAAGAATTCACTACTCAGAATCAGAAATCCCAGGATGTGGAACCTCTGGGAGGGAGAA
GTGGTCAAAGAGCTCTCTGTGGAAGAACAGATAAAGAGAAATCGGTATTATGATGAGGATGAGGATGAA
GAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001113546

**Insert Size:** 2283 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001113546.1</a></u>
<b>RefSeq Size:</b>	3748 bp
<b>RefSeq ORF:</b>	2283 bp
<b>Locus ID:</b>	51474
<b>UniProt ID:</b>	<u><a href="#">Q9UHB6</a></u>
<b>Cytogenetics:</b>	12q13.12
<b>MW:</b>	85.2 kDa
<b>Gene Summary:</b>	<p>This gene encodes a cytoskeleton-associated protein that inhibits actin filament depolymerization and cross-links filaments in bundles. It is downregulated in some cancer cell lines. Alternatively spliced transcript variants encoding different isoforms have been described for this gene, and expression of some of the variants maybe independently regulated. [provided by RefSeq, Aug 2011]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1).</p>