

Product datasheet for **SC319007**

NCX1 (SLC8A1) (NM_001112802) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NCX1 (SLC8A1) (NM_001112802) Human Untagged Clone
Tag:	Tag Free
Symbol:	NCX1
Synonyms:	NCX1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001112802, the custom clone sequence may differ by one or more nucleotides

```

ATGTACAACATGCGGCGATTAAGTCTTTCACCCACCTTTCAATGGGATTCATCTGTTA
GTTACTGTGAGTCTTTATTTTCCCATGTGGACCATGTAATTGCTGAGACAGAAATGGAA
GGAGAAGGAAATGAACTGGTGAATGTACTGGATCATATTACTGTAAGAAAGGGGTGATT
TTGCCCATTTGGGAACCCCAAGACCTTCTTTTGGGACAAAAATTGCTAGAGCTACTGTG
TATTTTGTGGCCATGGTCTACATGTTTCTTGGAGTCTCTATCATAGCTGATCGGTTTCATG
TCCTCTATAGAAGTCATCACATCTCAAGAAAAAGAAATAACCATAAAGAAACCCCAATGGA
GAGACCACCAAGACAACCTGTGAGGATCTGGAATGAAACAGTTTCTAACCTGACCTTGATG
GCCCTGGGATCTTCTGCTCCTGAGATTCTCCTTTCAGTAATTGAAGTGTGTGGCCATAAC
TTCCTGACGAGGACCTCGGTCCTAGCACCATCGTGGGAAGTGTGCATTCAATATGTTT
ATCATTATTGCACTCTGTGTTTATGTGGTGCCTGACGGAGAGACAAGGAAGATTAAGCAT
TTGCGTGTCTTCTTGTGACAGCAGCCTGGAGCATCTTTCCTACACCTGGCTTTACATT
ATTTTGTCTGTCATATCTCCTGGTGTGTGGAGGTCTGGGAAGGTTTGCTTACTTTCTTC
TTCTTCCCATCTGTGTTGTGTTGCTTGGGTAGCGGATAGGAGACTTCTGTTTTACAAG
TATGTCTACAAGAGGTATCGAGCTGGCAAGCAGAGGGGGATGATTATTGAACATGAAGGA
GACAGGCCATCTTCTAAGACTGAAATTGAAATGGACGGGAAAGTGGTCAATTCTCATGTT
GAAAAATTTCTTAGATGGTCTCTGGTTCTGGAGGTGGATGAGAGGGACCAAGATGATGAA
GAAGCTAGGCGAGAAATGGCTAGGATTCTGAAGGAACTTAAGCAGAAGCATCCAGATAAA
GAAATAGAGCAATTAATAGAATTAGCTAACTACCAAGTCTAAGTCAGCAGCAAAAAAGT
AGAGCATTTTATCGCATTCAAGCTACTCGCCTCATGACTGGAGCTGGCAACATTTTAAAG
AGGCATGCAGCTGACCAAGCAAGGAAGGCTGTCAGCATGCACGAGGTCAACACTGAAGTG
ACTGAAAATGACCCTGTTAGTAAGATCTTCTTTGAACAAGGGACATATCAGTGTCTGGAG
AATGTGTGACTGTGGCCCTTACCATTATCCGCAGAGGTGGTGGTGGTACTAACACTGTG
TTTGTGACTTCAAGACAGAGGATGGCACAGCAATGCTGGGTCTGATTATGAATTTACT
GAAGGAACTGTGGTGTAAAGCCTGGTGTATCCCAAGGAAATCAGAGTGGGTATCATA
GATGATGATATCTTTGAGGAGGATGAAAAATTTCTTGTGCATCTCAGCAATGTCAAAGTA
TCTTCTGAAGCTTCAAGATGGCATACTGGAAGCCAATCATGTTTCTACACTTGCTTGC
CTCGGATCTCCCTCCACTGCCACTGTAACATTTTTGATGATGACCACGCAGGCATTTTT
ACTTTTGAGGAACCTGTGACTCATGTGAGTGAAGCATTGGCATCATGGAGGTGAAAGTA
TTGAGAACATCTGGAGCTCGAGGAAATGTTATCGTTCCATATAAAACCATCGAAGGGACT
GCCAGAGGTGGAGGGGAGGATTTGAGGACACTTGTGGAGAGCTCGAATTCAGAATGAT
GAAATTGTGAAGATCATTACCATTAGAATATTTGACCGTGAGGAATATGAGAAAGAGTGC
AGTTTTCTCCCTGTGCTTGGGAACCAAAATGGATAAGAAGAGGAATGAAAGGTGGCTTC
ACAATAACAGACGAATATGATGACAAGCAGCCACTGACCAGCAAAGAGGAAGAGGAGAGG
CGCATTGCAGAAATGGGGCGCCCATCTGGGAGAGCACACCAAGTTGGAAGTGATCATT
GAAGAATCCTATGAATCAAGAGTACTGTGGACAACTCATTAAAGACAAACCTGGCC
CTTGTGGTTGGGACTAACAGCTGGAGAGAACAGTTCATTGAAGCTATCACTGTCACTGCT
GGGGAAGATGATGACGACGATGAATGTGGGAAGAGAAGCTGCCCTCCTGTTTCGATTAC
GTGATGCACTTTCTGACTGTGTTCTGGAAGGTCTGTTTGCCTTCGTCCTCCCTACTGAA
TACTGGAATGGCTGGGCGTGTTCATTGTCTCCATCCTCATGATTGGCCTACTGACAGCT
TTCATTGGAGACCTGGCTTCCCCTTTGGCTGCACCATTTGGCCTGAAAGATTCTGTGACT
GCAGTCGTGTTTCGTCGCACTTGGAAATCAGTGCCAGACACATTTGCCAGCAAAGTGGCA
GCCACCCAGGACAGTATGCAGACGCCTCCATAGGTAACGTACGGGCAGCAACGCGGTG
AATGTCTTCTGGGAATCGGTGTGGCCTGGTCCATCGCTGCCATCTACCACGCAGCCAAT
GGGGAACAGTTCAAAGTGTCCCCTGGCACACTAGCTTTCTGTCACTCTTTCACCATT
TTTGCTTTTCATCAATGTGGGGTGTGCTGTATCGGCGGAGGCCAGAAATCGGAGGTGAG
CTGGGTGGGCCCCGACTGCCAAGCTCCTCACATCCTGCCTCTTTGTGCTCCTATGGCTC
TTGTACATTTTCTTCTCCTCCCTGGAGGCTACTGCCACATAAAAGGCTTC

```

Restriction Sites: Please inquire

ACCN:	NM_001112802
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001112802.1</u> , <u>NP_001106273.1</u>
RefSeq Size:	6032 bp
RefSeq ORF:	2814 bp
Locus ID:	6546
UniProt ID:	<u>P32418</u>
Cytogenetics:	2p22.1
Protein Families:	Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Calcium signaling pathway, Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)

Gene Summary:

In cardiac myocytes, Ca(2+) concentrations alternate between high levels during contraction and low levels during relaxation. The increase in Ca(2+) concentration during contraction is primarily due to release of Ca(2+) from intracellular stores. However, some Ca(2+) also enters the cell through the sarcolemma (plasma membrane). During relaxation, Ca(2+) is sequestered within the intracellular stores. To prevent overloading of intracellular stores, the Ca(2+) that entered across the sarcolemma must be extruded from the cell. The Na(+)-Ca(2+) exchanger is the primary mechanism by which the Ca(2+) is extruded from the cell during relaxation. In the heart, the exchanger may play a key role in digitalis action. The exchanger is the dominant mechanism in returning the cardiac myocyte to its resting state following excitation.[supplied by OMIM, Apr 2004]

Transcript Variant: This variant (D) has multiple differences in the coding region but maintains the reading frame, compared to variant A. This variant encodes isoform D, also known as NaCa3, NCX1.3, or the kidney isoform, which is 36 aa shorter than isoform A. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.