

## Product datasheet for **SC319076**

### PTBP1 (NM\_002819) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PTBP1 (NM_002819) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTBP1
Synonyms:	HNRNP-I; HNRNPI; HNRPI; pPTB; PTB; PTB-1; PTB-T; PTB2; PTB3; PTB4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

**Fully Sequenced ORF:** >OriGene sequence for NM\_002819.3  
 CTCGGAGCCGTTGGGTCGGTTCCTGCTATTCCGGCGCCTCCACTCCGTCCTCCCGCGGGTC  
 TGCTCTGTGTGCCATGGACGGCATTGTCCAGATATAGCCGTTGGTACAAAGCGGGGATC  
 TGACGAGCTTTTCTACTTGTGTCACTAACGGACCGTTTATCATGAGCAGCAACTCGGC  
 TTCTGCAGCAAACGAAATGACAGCAAGAAGTTCAAAGGTGACAGCCGAAGTGCAGGCGT  
 CCCCTCTAGAGTGATCCACATCCGGAAGCTCCCCATCGACGTCACGGAGGGGGAAGTCAT  
 CTCCTGGGGCTGCCCTTTGGGAAGGTCACCAACCTCCTGATGCTGAAGGGGAAAAACCA  
 GGCCTTCATCGAGATGAACACGGAGGAGGCTGCCAACCCATGGTGAAGTACTACACCTC  
 GGTGACCCCTGTGCTGCGCGCCAGCCCATCTACATCCAGTTCTCCAACCACAAGGAGCT  
 GAAGACCGACAGCTCTCCAACCAGGCGGGCCAGGCGGCCCTGCAGGCGGTGAAGTCT  
 GGTCCAGTCGGGAACCTGGCCTTGGCTGCCTCGGCGGCGCCGTGGACGCAGGGATGGC  
 GATGGCCGGGCAGAGCCCGTGCTCAGGATCATCGTGGAGAACCTTCTACCCTGTGAC  
 CCTGGATGTGCTGCACCAGATTTTCTCCAAGTTCGGCACAGTGTTGAAGATCATCACCTT  
 CACCAAGAACAACCAGTTCAGGCCCTGCTGCAGTATGCGGACCCCGTGAGCGCCAGCA  
 CGCCAAGCTGTGCTGGACGGGCAGAATCTACAACGCTGCTGCACGCTGCGCATCGA  
 CTTTTCCAAGCTCACCAGCCTCAACGTCAAGTACAACAATGACAAGAGCCGTGACTACAC  
 ACGCCAGACCTGCCTTCCGGGGACAGCCAGCCCTCGTGGACCAGACCATGGCCGCGGC  
 CTTGCGTGACCTGGTATAATCTCAGCCTCTCCGTATGCAGGAGCTGGTTTCCCTCCCAC  
 CTTTGCATTCCCTCAAGCTGCAGGCCTTCCGTTCCGAACCTCACGGCGCCCTGGCCCC  
 CCTGGCCATCCCCTCGGCGGCGGCGCAGCTGCGGCGGACAGTTCGGATCGCCATCCCGGG  
 CCTGGCGGGGACAGAAATCTGTATTGCTGGTACGCAACCTCAACCCAGAGAGAGTAC  
 ACCCCAAAGCCTCTTTATTCTTTTCGGCGTCTACGGTGCAGTGCAGCGCTGAAGATCCT  
 GTTCAATAAGAAGGAGAACGCCCTAGTGCAGATGGCGGACGGCAACCAGGCCAGCTGGC  
 CATGAGCCACCTGAACGGGCACAAGCTGCACGGGAAGCCCATCCGCATCACGCTCTCGAA  
 GCACCAGAACGTGCAGCTGCCCCGAGGGCCAGGAGGACCAGGGCCTGACCAAGGACTA  
 CGGCAACTCACCCCTGCACCGCTTCAAGAAGCCGGGCTCCAAGAAGTCCAGAACATATT  
 CCCGCCCTCGGCCACGCTGCACCTCTCCAACATCCCGCCTCAGTCTCCGAGGAGGATCT



[View online »](#)

```

CAAGGTCCTGTTTTCCAGCAATGGGGGCGTCGTCAAAGGATTCAAGTTCTTCCAGAAGGA
CCGCAAGATGGCACTGATCCAGATGGGCTCCGTGGAGGAGGCGGTCCAGGCCCTCATTGA
CCTGCACAACCACGACCTCGGGGAGAACCACCACTGCGGGTCTCCTTCTCCAAGTCCAC
CATCTAGGGGACAGGCCCCACGGCCGGGCCCCCTGGCACAATTCCATCATTCCAGA
GAAAAGCCACTTTAAAAACAGCTGAAGTGACCTTAGCAGACCAGAGATTTATTTTTTTA
AAGAGAAATCAGTTTACCTGTTTTTAAAAAATTAATCTAGTTCACCTTGCTCACCCCTG
CGGTGACAGGGACAGCTCAGGCTCTGGTGACTGTGGCAGCGGGAGTTCCCGGCCCTCCA
CACCCGGGGCCAGACCCTCGGGGCCATGCCTTGGTGGGGCCTGTGTCGGGCGTGGGGCCT
GCAGGTGGGCGCCCCGACCACGACTTGGCTTCTTGTGCCTTAAAAAACCTGCCTTCTG
CAGCCACACACCCACCCGGGGTGTCTGGGGACCCAAGGGGTGGGGGGTACACCAGAG
AGAGGCAGGGGGCTGGCCGGCTCCTGCAGGATCATGCAGCTGGGGCGCGCGCCGCGG
CTGCGACACCCCAACCCAGCCCTAATCAAGTCACGTGATTCTCCCTTACCCCGCCC
CCAGGGCCTTCCCTTCTGCCCCAGGCGGGCTCCCGCTGCTCCAGCTGCGGAGCTGGTC
GACATAATCTCTGTATTATATACTTTGCAGTTGCAGACGTCTGTGCCTAGCAATATTTCC
AGTTGACCAAATATTCTAATCTTTTTTCAATTTATATGCAAAGAATAGTTTTAAGTAAC
TTTTTATAGCAAGATGATACAATGGTATGAGTGAATCTAAACTCCTTGTGGTATTACC
TTGTATGCTGTTACTTTTATTTTATTCCTTGAATTAAGTCACAGGCAGGACCCAGTTTC
CAGAGAGCAGGGGGGCGCCAGTGGGTGAGGCACAGGGAGCCCGGTCTATCTTAGA
GCCCTGAGCTTCAGGGAAGGGGCGGGCGTGTGCGCCCTCTGGCATCGCTCCGGTTGC
CTTACACCACGCCTTACCTGCAGTCGCTAGAAAAGTCTCAAACTTACAGGTTTTT
TTCTTCTTCAAATTTTGGACCAAAGTCTCATTCTGTGTTTTGCCTGCCTCTGATGCTG
GGACCCGGAAGGGGCGGCTCCTCCTGTCTTTGTGCTCTTCTACCGCCCCCGTCTCTG
TCCCGGGGCTCTCCTAGGATCCCCTTTCCGTAAGAGCGTGAACAAGGGTGAATATT
TATAATTTTTTATACCTGTTGTGAGACCCGAGGGGCGGGCGGTTTTTTATGGTGAC
ACAAATGTATATTTTGTAAACAGCAATCCAGGCTCAGTATTGTGACCGGGAGCCACAG
GGGACCCACGCACATTCGTTGCCTTACCCGATGGCTTGTGACGCGGAGAGAACCATT
AAAACCGTTTGGAACTCCTCCTTGTCTAGCCCTGTGTTGCTGTGGACGCTGTAGAG
GCAGGTTGGCCAGTCTGTACCTGGACTTCGAATAAATCTTCTGTATCCTCAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
    
```

**Restriction Sites:**

Please inquire

**ACCN:**

NM\_002819

**OTI Disclaimer:**

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:**

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_002819.3</a></u> , <u><a href="#">NP_002810.1</a></u>
<b>RefSeq Size:</b>	3321 bp
<b>RefSeq ORF:</b>	1674 bp
<b>Locus ID:</b>	5725
<b>UniProt ID:</b>	<u><a href="#">P26599</a></u>
<b>Cytogenetics:</b>	19p13.3
<b>Domains:</b>	RRM
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	<p>This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA-binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has four repeats of quasi-RNA recognition motif (RRM) domains that bind RNAs. This protein binds to the intronic polypyrimidine tracts that requires pre-mRNA splicing and acts via the protein degradation ubiquitin-proteasome pathway. It may also promote the binding of U2 snRNP to pre-mRNAs. This protein is localized in the nucleoplasm and it is also detected in the perinucleolar structure. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (a).</p>