

Product datasheet for **SC319079**

PPPDE1 (DESI2) (NM_016076) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPPDE1 (DESI2) (NM_016076) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPPDE1
Synonyms:	C1orf121; CGI-146; DESI; DeSI-2; DESI1; FAM152A; PNAS-4; PPPDE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_016076.3 GCCAACGGTCCGGCCCCGTCGACAGACGCTCCTGTCGGCGCGCCCGGGAGCGGCTCG GCTGCCCGATGCTTCCGCCCGGCTGCCGCGGGCCGGGCTGTACGCTTAGTGCCCGGCTC AGGCCCCCTGAAGCGCCCGCGGGGTGAGAGCGCCTCCGGCCCCGCGGAGACGGAGCGG CTTGAGGACGAGCGCGCGCGGGGAGGAGGATGGGGCTAACCAGTTAGTGGTGCTC AACGTGTACGACATGTATTGGATGAACGAATATACCTCATCCATTGGAATTGGAGTTTT CATTGAGAAATTGAAGTCTATGGCAGAGAATTTGCTTATGGTGGCCATCCTTACCCTTT TCTGGAATATTTGAAATTTCCCCAGGAAATGCTTCTGAACTAGGAGAAACATTTAAATTT AAAGAAGCTGTTGTTTTAGGGAGCACGGACTTCTAGAAGATGATATAGAAAAAATTGTA GAAGAAGTGGGAAAAGAATACAAAGGCAATGCTTATCATTTAATGCATAAAAACTGCAAT CATTTTTCTTCAGCTTTATCAGAGATTCTTTGTGGGAAAGAGATTCTCGCTGGATCAAT CGACTTGCTACTTCAGCTCCTGTATACCCTTTCTACAGAGTTGCCTCCGGAAGGAGTGG CTCACGCCCGCAGCCCTGCACTAGTGTGAGCAAGAACTCCAGGATGAACTGGAGGAA GCAGAGGATGCTGCCGATCCGCTTCCGTGGCAAGCACTGCAGCAGGCTCCAGACCCGGG CGCCACACTAACTATAAATGCTCCAAAGTCACACATTCAGAACTGTCTCTGGCAGTCG AATATCACTAGAGAAAAGTAAACAGAGAAGCATCCTTTAGATATTTGTATGCAAAGATG GCTCTCCCCAAATCCCAGTTTTTCAGCTCAGGATTATATTTGTAATCAAAAAAAAAAAAA AAAAAA
Restriction Sites:	Please inquire
ACCN:	NM_016076
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016076.3 , NP_057160.2
RefSeq Size:	4093 bp
RefSeq ORF:	585 bp
Locus ID:	51029
UniProt ID:	Q9BSY9
Cytogenetics:	1q44
Gene Summary:	<p>Has deubiquitinating activity towards 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains. Deubiquitinates 'Lys-48'-linked polyubiquitination of RPS7 leading to its stabilization (PubMed:28483520).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>