

Product datasheet for **SC319271**

CD151 (NM_004357) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD151 (NM_004357) Human Untagged Clone
Tag: Tag Free
Symbol: CD151
Synonyms: GP27; MER2; PETA-3; RAPH; SFA1; TSPAN24
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_004357.4
CGCCCCGCAGCTGCCGCCGCCAGGGCCCCGACTCGGACCGCTGCCTAGAGTCTGG
GGAGCTTCTGTCCACCTGTCCTGCAGAGGAGTCGTTTCCAGCCCCGCTGCCCCAGGATGG
GTGAGTTCAACGAGAAGAAGACAACATGTGGCACCGTTTTGCCTCAAGTACCTGCTGTTTA
CCTACAATTGCTGCTTCTGGCTGGCTGGCTGTCATGGCAGTGGGCATCTGGACGC
TGGCCCTCAAGAGTGACTACATCAGCCTGCTGGCCTCAGGCACCTACCTGGCCACAGCCT
ACATCCTGGTGGTGGCGGGCACTGTCGTCATGGTGACTGGGGTCTTGGGCTGCTGCGCCA
CCTTCAAGGAGCGTCGGAACCTGCTGCGCCTGTACTTCATCCTGCTCCTCATCATCTTTC
TGCTGGAGATCATCGCTGGTATCCTCGCCTACGCCTACTACCAGCAGTGAACACGGAGC
TCAAGGAGAACCCTGAAGGACACCATGACCAAGCGCTACCACCAGCCGGCCATGAGGCTG
TGACCAGCGCTGTGGACCAGCTGCAGCAGGAGTCCACTGCTGTGGCAGCAACAACCTCAC
AGGACTGGCGAGACAGTGAAGTGGATCCGCTCACAGGAGGCCGGTGGCCGTGTTGGTCCCAG
ACAGCTGCTGCAAGACGGTGGTGGCTCTTTGTGGACAGCGAGACCATGCCTCCAACATCT
ACAAGGTGGAGGGCGGGTGCATACCAAGTTGGAGACCTTCATCCAGGAGCACCTGAGGG
TCATTGGGGCTGTGGGGATCGGCATTGCCTGTGTGCAGGTCTTTGGCATGATCTTCACGT
GCTGCCTGTACAGGAGTCTAAGCTGGAGCACTACTGACCCTGCCTTGGGCCTTGTCTGCT
GCTGCACCCAACCTACTGAGCTGAGACCACTGAGTACCAGGGGCTGGGCTCCCTGATGACA
CCCACCCTGTGCCATCACCATAACCTCTGGGGACCCCAACCTCAGAGGCAGCTTCAAGTG
CCTTTTGTGCGCACCAATGCCAGCAGGGGAGGTGAGGGGGCTGGCGGGCGAAGTTT
GGGGGGTGTGTTTGTGGGGCTCCCCGACATACTCTGCTGCTGGTGGTGCAGATGCAGGTTG
GAAGGGGCTTGTGAGTGGCGCAAGGCCGAGCGTTCACAGCAGGGGAGAAACCTTCA
CACCCAGGCCCTCAGGAACTGGGGCTTTCCTTGCAGCCACATGGCCCATCCCAGTT
GGGGAAGCCAGGTGAGCTCTGACCCTTGGGCTGGGCTCTGCCCTCCCAACCCAGCCG
TCGTCTCCCTCGACAGCGCCCTGCTGTCTTCCCACCGCAGTACCACCACCCGAAATG
CCACGTGGTCACTGTGCACTGCCCTGTTTATGTGCTCTGCGGGGAGGGCTTCTGCTG
TTTGTACTGCTGTACCAGATGCCTACAACCATCCCTGCCACATACAGGTGCTCAATA
AACACTTGTAGAGCAAAAAAAAAAAAAAAAAAAAAA



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_004357
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_004357.4 , NP_004348.2
RefSeq Size:	1574 bp
RefSeq ORF:	762 bp
Locus ID:	977
UniProt ID:	P48509
Cytogenetics:	11p15.5
Domains:	transmembrane4
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	<p>The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins and other transmembrane 4 superfamily proteins. It is involved in cellular processes including cell adhesion and may regulate integrin trafficking and/or function. This protein enhances cell motility, invasion and metastasis of cancer cells. Multiple alternatively spliced transcript variants that encode the same protein have been described for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) is the longest transcript. All variants encode the same protein.</p>