

Product datasheet for **SC319373**

CD320 (NM_016579) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD320 (NM_016579) Human Untagged Clone
Tag: Tag Free
Symbol: CD320
Synonyms: 8D6; 8D6A; TCBLR; TCN2R
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_016579.2
 CTGGACAGCGCGTGGCCGGCCGCTGTGGGGACAGCATGAGCGCGGTTGGATGGCGCA
 GGTTGGAGCGTGGCGAACAGGGGCTCTGGGCTGGCGCTGCTGCTGCTCGGCTCGG
 ACTAGGCCTGGAGCGCCGCGAGCCCGCTTCCACCCGACCTTGCCAGGCCGAGG
 CCCAGCTCAGGCTCGTGCCACCCACCAAGTTCCAGTGCCGACCAAGTGGCTTATGCGT
 GCCCTCACCTGGCGCTGCGACAGGACTTGGACTGCAGCGATGGCAGCGATGAGGAGGA
 GTGCAGGATTGAGCCATGTACCCAGAAAGGGCAATGCCACCGCCCTGGCTCCCTG
 CCCCTGCACCGCGCTCAGTACTGCTCTGGGGAACTGACAAGAACTGCGCAACTGCAG
 CCGCTGGCTGCCTAGCAGGCGAGCTCCGTTGCACGCTGAGCGATGACTGCATTCCACT
 CACGTGGCGCTGCGACGGCCACCAGACTGTCCGACTCCAGCGACGAGCTCGGCTGTGG
 AACCAATGAGATCCTCCCGAAGGGGATGCCACAACCATGGGGCCCTGTGACCCTGGA
 GAGTGTCACCTCTCTCAGGAATGCCACAACCATGGGGCCCTGTGACCCTGGAGAGTGT
 CCCCTCTGTGGGAATGCCACATCCTCCTCTGCCGGAGACCAGTCTGGAAGCCAACTGC
 CTATGGGGTTATTGCAGCTGCTGCGGTGCTCAGTGCAAGCCTGGTACCGCCACCCTCCT
 CCTTTTGTCTGGCTCCGAGCCAGGAGCGCTCCGCCACTGGGGTTACTGGTGGCCAT
 GAAGGAGTCCCTGCTGCTGTGAGAACAGAAAGACCTCGCTGCCCTGAGGACAAGCACTTGC
 CACCACCGTCACTCAGCCCTGGGCGTAGCCGGACAGGAGGAGAGCAGTATGCGGATGGG
 TACCCGGGCACACCAGCCCTCAGAGACCTGAGCTCTTCTGGCCACGTGGAACCTCGAACC
 CGAGCTCCTGCAGAAGTGGCCCTGGAGATTGAGGGTCCCTGGACACTCCCTATGGAGATC
 CGGGGAGCTAGGATGGGGAACCTGCCACAGCCAGAAGTGGGGGCTGGCCCGAGGAGCT
 CCCAGGGGTAGAACGGCCCTGTGCTTAAAGACTCCTGCTGCCCGTCTGAGGGTGGCG
 ATTAAGTTGCTTACATCCTCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire
ACCN: NM_016579



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016579.2 , NP_057663.1
RefSeq Size:	1260 bp
RefSeq ORF:	849 bp
Locus ID:	51293
UniProt ID:	Q9NPF0
Cytogenetics:	19p13.2
Domains:	Idl_recept_a
Protein Families:	Transmembrane
Gene Summary:	<p>This gene encodes the transcobalamin receptor that is expressed at the cell surface. It mediates the cellular uptake of transcobalamin bound cobalamin (vitamin B12), and is involved in B-cell proliferation and immunoglobulin secretion. Mutations in this gene are associated with methylmalonic aciduria. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2011]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1).</p>