

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for SC319469

Casein Kinase 1 delta (CSNK1D) (NM_001893) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Casein Kinase 1 delta (CSNK1D) (NM_001893) Human Untagged Clone
Tag:	Tag Free
Symbol:	Casein Kinase 1 delta
Synonyms:	ASPS; CKI-delta; CKId; CKIdelta; FASPS2; HCKID
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Fully Sequenced ORF:

>OriGene sequence for NM_001893.3 GCGGCGGCTCCTTTAGGCAGCTGAAAGGGGATTTAGGCCCGGAAGATCCGAGTCCATCCG CGGCGGGGAGAGGGCAAGCGGGACCGGTAGGGGCCGGAGCAGCGGCGGCGGCGCCCGGAC CGCCGCCTCCCAGAGTTCAGACCCAGGAAGCGGCCGGGAGGGCAGGAGCGAATCGGGCCG CCGCCGCCATGGAGCTGAGAGTCGGGAACAGGTACCGGCTGGGCCGGAAGATCGGCAGCG GCTCCTTCGGAGACATCTATCTCGGTACGGACATTGCTGCAGGAGAAGAGGTTGCCATCA AGCTTGAATGTGTCAAAACCAAACACCCTCAGCTCCACATTGAGAGCAAAATCTACAAGA TGATGCAGGGAGGAGTGGGCATCCCCACCATCAGATGGTGCGGGGCAGAGGGGGGACTACA ACGTCATGGTGATGGAGCTGCTGGGGGCCAAGCCTGGAGGACCTCTTCAACTTCTGCTCCA GGAAATTCAGCCTCAAAACCGTCCTGCTGCTGCTGACCAAATGATCAGTCGCATCGAAT ACATTCATTCAAAGAACTTCATCCACCGGGATGTGAAGCCAGACAACTTCCTCATGGGCC TGGGGAAGAAGGGCAACCTGGTGTACATCATCGACTTCGGGCTGGCCAAGAAGTACCGGG ATGCACGCACCCACCAGCACATCCCCTATCGTGAGAACAAGAACCTCACGGGGACGGCGC GGTACGCCTCCATCAACACGCACCTTGGAATTGAACAATCCCGAAGAGATGACTTGGAGT CTCTGGGCTACGTGCTAATGTACTTCAACCTGGGCTCTCTCCCCTGGCAGGGGCTGAAGG CTGCCACCAAGAGACAGAAATACGAAAGGATTAGCGAGAAGAAAATGTCCACCCCATCG AAGTGTTGTGTAAAGGCTACCCTTCCGAATTTGCCACATACCTGAATTTCTGCCGTTCCT TGCGTTTTGACGACAAGCCTGACTACTCGTACCTGCGGCAGCTTTTCCGGAATCTGTTCC ATCGCCAGGGCTTCTCCTATGACTACGTGTTCGACTGGAACATGCTCAAATTTGGTGCCA GCCGGGCCGCCGATGACGCCGAGCGGGAGCGCAGGGACCGAGAGGAGCGGCTGAGACACT CGCGGAACCCGGCTACCCGCGGCCTCCCTTCCACAGCCTCCGGCCGCCTGCGGGGGGACGC AGGAAGTGGCTCCCCCACACCCCTCACCCCTACCTCACACGGCTAACACCTCCCCCC GGCCTGTCTCCGGCATGGAGAGAGAGGGGAAAGTGAGTATGCGGCTGCACCGCGGGGCCC CCGTCAACATCTCCTCGTCCGACCTCACAGGCCGACAAGATACCTCTCGCATGTCCACCT CACAGATTCCTGGTCGGGTGGCTTCCAGTGGTCTTCAGTCTGTCGTGCACCGATGAGAAC TCTCCTTATTGCTGTGAAGGGCAGACAATGCATGGCTGATCTACTCTGTTACCAATGGCT TTACTAGTGACACGTCCCCCGGTCTAGGATCGAAATGTTAACACCGGGAGCTCTCCAGGC CACTCACCCAGCGACGCTCGTGGGGGGAAACATACTAAACGGACAGACTCCAAGAGCTGCC ACCGCTGGGGCTGCACTGCGGCCCCCCACGTGAACTCGGTTGTAACGGGGCTGGGAAGAA **GGTGGCCGCCCTGTACTCCCTGACGATTCCACTGTAACTACCAATCTTCTACTTGGTTAA**


Restriction Sites: ACCN: Please inquire

NM_001893

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	asein Kinase 1 delta (CSNK1D) (NM_001893) Human Untagged Clone – SC319469
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by calling 301.340.3188 option 3 for pricing and delivery.
	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Meth	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 001893.3, NP 001884.2</u>
RefSeq Size:	2030 bp
RefSeq ORF:	1248 bp
Locus ID:	1453
UniProt ID:	<u>P48730</u>
Cytogenetics:	17q25.3
Domains:	pkinase, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Circadian rhythm - mammal, Gap junction, Hedgehog signaling pathway

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Gene Summary:This gene is a member of the casein kinase I (CKI) gene family whose members have been
implicated in the control of cytoplasmic and nuclear processes, including DNA replication and
repair. The encoded protein may also be involved in the regulation of apoptosis, circadian
rhythm, microtubule dynamics, chromosome segregation, and p53-mediated effects on
growth. The encoded protein is highly similar to the mouse and rat CK1 delta homologs.
Three transcript variants encoding different isoforms have been found for this gene.
[provided by RefSeq, Feb 2014]
Transcript Variant: This variant (1) encodes the longer isoform (1).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US