

Product datasheet for **SC319721**

PLGLB1 (NM_001032392) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PLGLB1 (NM_001032392) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLGLB1
Synonyms:	PLGL; PLGP1; PRGB; PRP-B
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<p>>OriGene sequence for NM_001032392.2</p> <pre> ATGTAAGTCAACAACATCCTGGGATTGGGACACACTTTCTGGGCACTGCTGGCCAGTCCC AAAATGGAACATAAGGAAGTGGTTCTTCTACTTCTTTTATTTCTGAAATCAGGTCAAGGA GAGCCCCCTGGATGACTATGTGAATACCCAGGGCCTTCACTGTTCAAGTGTCACTAAGAAG CAGCTGGGGGCAGGAAGCAGAGAAGAATGTGCAGCAAAATGTGAAGAGGACAAAGAATTC ACCTGCAGGGCATTCCAATATCACAGTAAAGAGCAACAATGTGTGATAATGGCTGAAAAC AGGAAGTCCTCCATAATCATTAGGATGAGAGATGCAGTTTTATTTGAAAAGTAAATGTAT CTTTTCAGAGTGAAGACTGGGAATGAAAGAAGTACAGAGGGACGATGTCCAAAACAAAA AATGGCATCACCTGTCAAAAATGGAGTTCCTCTCCCCGCAGACCTAGGTCAGACTTT TCCTTTCATCTTTGTGTTTCTACTGTAAAGTGTCCGCTGTGTCTGTGAGGGATTGGT TCCAGGACCCCTGTGGCTACCAAAATCCATGCTTCTCAAGTCCCTTATGTAAAATGGTGC AGTATTTGCATATAACCTACATACCTTCTTGTGCAATCCCTAATATAATGTAATGCT ATGTAATCGTTGTTATACTGTATTGTTTTATTTGTATTATGTTTTATTGTCATATTGTT ATTTTCTCTCATTTTTCAAGTCTTTCCATCCACAGTTGGTTGAATTTGTGGATGTGG AACCCATGGATACAGAGGGCCAAGTGTATTTAGGATAAATTCATCACTTTTAATTCAAAC CAGAATATGTGAATAAGCAGACAGAAAGAATCTTTTGTGTCGATGTTCAACTATTTTT GGCACCATAGTAGAACATGGTTGCTTTCTATTTTTCTTGGGATGGAGGTTTCTTGAAG ACCTAGAACATAGAAGAATGCCTAGTAAAAAAAAAATCAATGAAACTACGAGTTTTAG GCCAAATCTGAGAAAAGACCAAGATGACTATGTTTGGGACTGAAGTAAACATATCAAGT TAGAACTCTCATCACATGTTTGACTCAAATTTGTGGAGCAAAATAGTAAATAAAATATAAA AATGAAAATGTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001032392



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001032392.2</u> , <u>NP_001027564.1</u>
RefSeq Size:	3011 bp
RefSeq ORF:	291 bp
Locus ID:	5343
UniProt ID:	<u>Q02325</u>
Cytogenetics:	2p11.2
Gene Summary:	May bind noncovalently to lysine binding sites present in the kringle structures of plasminogen. This may interfere with the binding of fibrin or alpha-2-antiplasmin to plasminogen and may result in the localization of activity at sites necessary for extracellular matrix destruction.[UniProtKB/Swiss-Prot Function]