

Product datasheet for **SC319863**

BOD1 (NM_138369) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: BOD1 (NM_138369) Human Untagged Clone
Tag: Tag Free
Symbol: BOD1
Synonyms: FAM44B
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_138369.1
 GCCGTTTCTGCCGCTCCACCACAACCTCTGCTGTATCAGAGGGGCCTTCGCCGCCGCC
 GCCCCACACCGCCCTTCCCGCCATCGCCACCTCCGCCGCTGCCGCCTGGGGCCCTCC
 TCCTTACCGCCCTTACCACCTCTACACGTTCTAGGCCTTCTGTCCTGGAGAAGAAG
 CTATAGTCGTTTCCCTTGTGGGCCCGGGCGCAGCCATGGCGGACGGCGCGCGCGG
 GGGAACTGGCGCGGTGGCGGCGCGGAACTAGCCAGGCCTTCCGGGGCAGCGACTGG
 CGCTACTGGGGCCAGCGGGGCGGTGGCCCATCAACCCGGCCTCGCTGCCTCCCGCGA
 CCCGCAGCTCATCGTCTCATCGTGGAGCAGCTCAAGAGCCGGGCCTTTTTGACAGCTT
 CCGCCGGGACTGCCTGGCCGACGTGGACACCAAGCCAGCTTACCAAACTGAGGCAGAA
 AGTGGATAATTTTGTGCAACACATCTGGACAAGCAGGAATGGAATCCTACGATGAACAA
 AAACCAGTTGCGAAATGGTCTGAGGCAGAGTGTGGTTCAGTCAGGGATGTTGGAAGCTGG
 AGTAGACAGGATTATTTCTCAGGTGGTGGATCCAAACTTAACCACATCTTCAGGCCACA
 AATAGAACGAGCAATTCATGAGTTCTTCCGCGCCAGAAAAAGCAGCTGTGCCAGCACC
 CCCTCCAGAGCCGAAGGCCAGGACCCTCCAGCTCCATCTCAGGACACTTCTAAGAATA
 CGCCAGACACCTTTTAAAGCTAATTTTGGTGAAGAAATGGATTTCGTTACATAAGAGT
 GCAACTTCAGACTGAAGTAGGCCAAGGTGCTCACTGATCTCAAGATTTCAACCTTGACC
 ATGGGCAGTGACCAGATTGAAAGGGGAGCAAGTTCGGCAGTGGGAGAGTTGACCGTGCA
 CCCCTGCATTGTGCTGCCATTTGGCCAGCCTGTCCAAGGCATGACACCAAGTAGACAC
 TACAGAGAGAGAAACACTACAGCAACCCAGGTTGCTTGAACAGACTTTTATACTTGA
 ACATGGAGACTGCACATGGACTTTAGGGTTTGTGCTGTGGGATAAACGGAAGCTACAGTG
 AGAACATAGCCAGTCCCAAAGACAATTTCAAAGAAAAATGACAGTAAAGATTAGCTGGGA
 GTAGTCTTTGACAGTGCTTATTTGATACTGTCTCTCAGAGTTTGCAAACAGATTGTACA
 AGTCATTAGCGTCAGATAGCTTTAAAGTTGTGACCTTCTGTACATGAATCTTCTAGCCA
 GTTTCCTTTCCTTTGTAACGAAACATGAAATCCTAGAATGTATGAGAAGTTCAGACATTA
 GGCATAAGGAACTCGTTTGCAGGCTCTGTCCAGGGCTGCTTCTGCTGCTGGAGGGG
 CAGTGAGTCTTAGGTATGTTATTTTATTCTCACATTTGTGTTTTTTAGAAAAGTGAAT
 GGTCATAAATGGCTTATCTTTCATAATAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_138369
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_138369.1</u> , <u>NP_612378.1</u>
RefSeq Size:	1554 bp
RefSeq ORF:	558 bp
Locus ID:	91272
UniProt ID:	<u>Q96IK1</u>
Cytogenetics:	5q35.2
Gene Summary:	<p>Required for proper chromosome biorientation through the detection or correction of syntelic attachments in mitotic spindles.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>