

## Product datasheet for **SC319876**

### Protein **FAM195A (MCRIP2) (NM\_138418) Human Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Protein FAM195A (MCRIP2) (NM_138418) Human Untagged Clone
Tag:	Tag Free
Symbol:	Protein FAM195A
Synonyms:	C16orf14; c349E10.1; FAM195A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

**Fully Sequenced ORF:** >OriGene sequence for NM\_138418.2  
CGAGCCCGTGC GGCCCTTTAAGGGCCGGGGCGTGTAGCGGGCCCGCCCTCCCCGCG  
GCGCCCGCAGTCCGTTAAGTGCAGCCCCGGCGCAGGGCCGGATCTGGCCGGGGCCGG  
CGGCGGTGTGGGAGCGGCGCATGTACACCATACCAAGGGGCCAGCAAGCTGGTCCG  
CGCAGCGCCGACAGGTCCCACGCAGCAGCAGGTGGAGGGCCGGCTCGGCGAGCTCCTGA  
AATGCCGGCAGCCCGCGCCGACCTCGCAGCCCCGCGGGCGCAGCCCTTTGCGCAGC  
CGCCGGGACCCTGGCCCTGTGAGTCCAGGGCCAAGGCTTGTGTTCAATCGTGTGAATG  
GCCGGCGGGCCCTCCACGTCCCCATCCTTCGAGGGGACCCAGGAGACCTACACAGTGG  
CCCACGAGGAGAATGTCCGCTTTGTGTCGAAGCCTGGCAGCAGGTGCAACAGCAGCTGG  
ATGGTGGCCAGCCGGTGAAGGGCCGCAAGGCCTGTGCAGTACGTGGAGAGACCCCA  
ATCCCCGCTGCAGAATTTGTGCCATTGACCTAGACGAGTGGTGGGCGCAGCAGTTCC  
TGGCGAGAATCACCAGCTGTTCTAGTGGCTGCTGGGAGGGGGCGCTGCTACACGGCCGA  
CCTGTGCGCAGGAGAGAAGCATGGCGCCCTGCCACCCACTGCGCCTGGCTGGGTGCCGG  
CCACACCTGAAGTGCCAGCATTTGGACTTTTGCACCTTTTTTCCCTTGGCCCGGTGTC  
CCAACCAAGCTGCCATGGCCAAGGGCCGAACCCGTCTGACCTCAGCCCTGCTCACTGTGC  
CCAGGGACCAGCGACCAGCCCTGGGGCTGGCAGGGAGGAGCTCCAGGCTAATAAAGTGG  
AGAAACTGTCAAAAAAAAAAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_138418
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_138418.2</a></u> , <u><a href="#">NP_612427.2</a></u>
<b>RefSeq Size:</b>	928 bp
<b>RefSeq ORF:</b>	483 bp
<b>Locus ID:</b>	84331
<b>UniProt ID:</b>	<u><a href="#">Q9BUT9</a></u>
<b>Cytogenetics:</b>	16p13.3