

Product datasheet for **SC319882**

FBXL6 (NM_012162) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	FBXL6
Synonyms:	FBL6; FBL6A; PP14630
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



Fully Sequenced ORF:	>OriGene sequence for NM_012162.1 GGTGTGCTGCTATAAGCCTCGCCCGTTCCGCTCCCTGGGGCTTCCCCGAGCGCCGTCGGTGGTCATGGCTGCCCCAGCCTCCCGGCAGGTCCGACGACAGCCCGGGCAGCGCCGCGGCCGCTCGGCCGAGGACTGGTGGTGGGACCGCTGGCGCCGAGGGGCTCGGGGTACCACTGCTGCAGTCCGACAGCATGCTGCTGGTGTGTCCGAACCCGGCCCGCCCGGCCCGCGCACAGCGGCGCGCTTCCCGCCCACTCCCGGCAGCCGCCCGGGGCCCCAGCGCCGCGGCCAAGCCCCAAGGCCGGGCTCAGGTCCGAGGCGGCGCGCCCGCACCCGACCCGGCACCCACGCCCCACGCCGAGGAAGGGCCGACGCGGGCTGGGAGACCGCATTCCTTGGAATCCTGTGTCAGATTTTCGGGTTGTTGGTGGCGGCGACGGCCCATGCCCTTCTGGGCAGGGCTGCTCGCTGTGCCCGCTGGCAGGAGGCGCTTCCCAACCCGCGCTCTGGCACA CCGTGACCTGTGCTCCCGCTGGTGGCGCGGCTGCCAAGGGCGGGTCAAGCGGAGAGAAGAGCTCCTTGCTTCCCTGGAGTGGCTTATGCCCAATCGGTTTTACAGCTCCAGAGGCTGACCCTCATCCACTGGAAGTCTCAGGTACACCCGTTTGAAGTGGTAGGTGAGTGTGTCTCGGCTCACTTTCCTCAAGCTCTCCGGTGCCACGGTGTGACTGCTGACGCTCTGGTCATGTAGCCAAAGCCTGCTGCCAGCTCCATAGCCTGGACCTACAGACTCCATGGTGGAGTCCACAGCTGTGGTGTGAGTCTTGGAGGAGGAGGGTCCCGAATGCGCAAGTTGTGGCTGACCTACAGCTCCAGACGACAGCCATCCTGGGCGCACTGCTGGGAGCTGCTGCCCCAGCTCCAGGTCCTGGAGGTGAGCACCAGGATCAACCGTAATAGCATTCCCTTCAGTGTCTGTGAGGCTCTGCAGAAAGGCTGCCCTCAGTCCAGGTGCTGCGGCTGTTGAACCTGATGTGGTGTGCCAAGCCTCCGGACGAGGGTGGTCCCGGACCAGGCTTCCCTAGCCTAGAGGAGCTCTGCTGGCGAGCTCAACCTGCAACTTTGTGAGCAACGAGGTCTGGGCCGCC TACTCCACGGCTCTCCCAACCTGCGCTTACTGGATCTTCGTGGTGTGCGCGCATCACGCGGGCTGGCCTTCAAGATCTGCCATGTCCGGAGCTGGAGCAGTTTCTGGGCTGTATGGCAGCTGACACCGGCTGACTCTAGCAAGGAGGCGAGCCCTTTTGGACCCAGAAAGTGGTGCCATACACTGCGAGAACTGGACTTGAGTGGCCAGGGGTTCAAGTGAAGGACCTGGAGCAGGCCCTGGCTGCTTCTTAAGCACCCCTGGGGCTCACACCCAGCCCTGTGCTCTTAACTCAGGGGACCCGGGTACACCAAGCACTGTGAGTCTGTGATCAGCAGTGCCCGGGCTGCTCTACCTCAACCTGGAGTCTGCGCTGCTTCCCGGGGCTGAAGCGGGCTACCGGGGCTGGAGGAAGTCCAGTGGTGTCTGGAGCAGCTGCTCACCAGCCCTCACCGAGCTAGGACGACAGCTGGGACACCTCAGCCAGCTTGCCCCACCCTCCACCTTTGCCCAATTTGAGATTTTGTTAAAATAAAACATTTTAGGAAAAAAAAAAAAAAAAAA
Restriction Sites:	Please inquire
ACCN:	NM_012162
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_012162.1, NP_036294.1</u>
RefSeq Size:	1799 bp
RefSeq ORF:	1620 bp
Locus ID:	26233
UniProt ID:	<u>Q8N531</u>
Cytogenetics:	8q24.3
Domains:	LRR, F-box, LRR_CC
Gene Summary:	This gene encodes a member of a family of proteins that are characterized by an F-box motif. The encoded protein also contains leucine-rich repeats. F-box-containing proteins comprise one of the subunits of the SCF (SKP1-cullin-F-box) complex, which functions in phosphorylation-dependent ubiquitination. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).