

Product datasheet for **SC319912**

FGFR1 Oncogene Partner (FGFR1OP) (NM_194429) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR1 Oncogene Partner (FGFR1OP) (NM_194429) Human Untagged Clone
Tag:	Tag Free
Symbol:	FGFR1 Oncogene Partner
Synonyms:	FGFR1OP; FOP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_194429.1
GAGAGAGAGATGGCGCGACGGCGCCGAGTGGTGGCCGAGGAGACACGGAGCTGCGG
GACCTGCTGGTGCAGACGCTGGAGAACAGCGGGTCTGAACCGCATCAAGGCTGAAC
CGAGCAGCTGTGTTTTAGCACTAGAGGAGCAAGAAAAAGTAGAGAACAAAACCTTTA
GTTAATGAGAGCCTGAGAAAGTTTTAAATACCAAAGACGGTCGTTTGTGGCTAGTCTT
GTTGCAGAAATTTCTCAGTTTTTAACTTACTTTACTTTGGCTGTTTTCAACCTGAA
ACTAGCACACTGCAAGGTCTCGAAGGTGAGAGAATTTAGCCCGAGATTTAGGTATAATT
GAAGCAGAAGGTACTGTGGGTGGACCTTATTATTAGAAGTGATCAGGCGCTGTCAACAG
AAAGAAAAGGGCCAACCACTGGGGAAGGTGCACTTGATCTATCTGATGTACATTCTCCA
CCAAAGTCACCAGAGGAAAAACAAGTGCACAGACAACCAAGTAAGAAGGCCAATGAT
GAGGCCAATCAGAGTGATACAAGTGTCTCCTTGTGAGAACCAAGAGCAAAGCAGCCTT
CACTTACTGTCCCATGAAACAAAATTGGATCTTTTCTAAGCAACAGAACTTTAGATGGC
AAAGACAAAGCTGGCCTTTGTCCAGATGAAGATGATATGGAAGGAGATTCTTTCTTTGAT
GATCCCATTCTAAGCCAGAGAAAACCTTACGGTTTGAGGAATGAACCTAGGAAGCAAGCA
GGAAGTCTGGCCTCGCTCTCGGATGCACCCCTTAAAAAGTGGACTCAGCTCCCTGGCG
GGAGCCCTTCTTTAAAAGACTCTGAGAGTAAAAGGGGAAATACAGTTTTGAAAGATCTG
AAATTGATCAGTGATAAAATTGGATCACTTGGATTAGGAACTGGAGAAGATGATGACTAT
GATAGAGAAGAAGACCTTTCTGTGGAAATAGATGACATCAATACCAGTGATAAGCTTGAA
GACCTCACACAAGATCTGACTGTATCCCAGCTCAGTGATGTTGCGGATTATCTGGAAGAT
GTTGCATAGACACGAAGAAGGAAGTATTCTAATTAACAAGGACAGAGGACTGACCGGTTT
CATTTTTTTTTTCCAGACAATCACTCAGCTGGAATGTCTGCTCTTATTGGTGCCTTG
CATTTCAAAAACACTGCAGATATTTTTAAAAGTAATTTTCACTTTACTAAACAAAATAC
TTCCTATTTGAGCCCATGTGTGGAAGATTAATATTCTTAATTTAGCTGTACATTTCTTT
ATGGAAATTGATTATCTACATTCAGTTTCATTACAGGGAAGGAACCCATGAAAACATCAG
TGTTAAGAGCATGATGAAAGGTGCAATAAAGCCGTAGGATCACGCAAAAAAAAAAAAAA
AAAAAA



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_194429
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_194429.1 , NP_919410.1
RefSeq Size:	1842 bp
RefSeq ORF:	1140 bp
Locus ID:	11116
UniProt ID:	O95684
Cytogenetics:	6q27
Protein Families:	Druggable Genome
Gene Summary:	<p>This gene encodes a largely hydrophilic centrosomal protein that is required for anchoring microtubules to subcellular structures. A t(6;8)(q27;p11) chromosomal translocation, fusing this gene and the fibroblast growth factor receptor 1 (FGFR1) gene, has been found in cases of myeloproliferative disorder. The resulting chimeric protein contains the N-terminal leucine-rich region of this encoded protein fused to the catalytic domain of FGFR1. Alterations in this gene may also be associated with Crohn's disease, Graves' disease, and vitiligo. Alternatively spliced transcript variants that encode different proteins have been identified. [provided by RefSeq, Jul 2013]</p> <p>Transcript Variant: This variant (2) lacks an in-frame exon in the internal coding region, compared to variant 1. The encoded isoform (b) is shorter, compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>