

Product datasheet for **SC319955**

TCL1B (NM_004918) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TCL1B (NM_004918) Human Untagged Clone
Tag:	Tag Free
Symbol:	TCL1B
Synonyms:	SYN-1; TML1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_004918.2

```

TGGAAAGCTACACGTGTGAGCCTAGAGGCGGGTCCCGTTGCAGACTTGCCATGGCCTCC
GAAGCTTCTGTGCGTCTAGGGGTGCCCCCTGGCCGTCTGTGGATCCAGAGGCCTGGCATC
TACGAAGATGAGGAGGGGAGAACCTGGGTGACTGTGGTCTGCGGTTCAATCCCTCGCGT
AGGGAATGGGCCAGGGCCTCCCAGGGCAGCAGATGAACCCAGCATCACAGTGCAGTTG
TGGCAGATGGCAGTGCATACCCGGGAGCTACTCTCCTCCGGCCAGATGCCCTTCTCCAG
CTGCCCCGCGTGTGGCAGCTCTACCCAGGAGGAAGTACCGAGCAGCGGATTCCAGTTTC
TGGGAAATAGCAGACCATGGCCAGATTGACTCTATGGAGCAGCTGGTCTAACATATCAG
CCGGAGAGGAAAGACTGACACTGGGAGTGGCTGGCCCTGCTGGCCCTGCCTCTTCTGGCC
TGGTGTCTCCTCATGCCCCCTCAGTGAGGATCTTCATGTACCTGCTCTTCTGTTTGCACA
CCCAGCATAGCCTCCTTGCAGGCAGAAGGCAGTAGGGCCCTGCACACTCAGTTTCTCTC
GTTTTCTTAGTTATCAGTCCTGTCTGTCCCACTCAGGTCTGTACTTAGGGCAGCTGGC
CTGGATGGGCTTCACTGGGGCCCTGTCTGTGTGCTGAGCCAGTTTCCCTGCTGGCTGCA
AGCTGTGGGTCTTTCTCCTCTGTGCCCTCATGCTGATCTTCTAGATGCCACTCCCAA
TCCCTTCATACCCACCAGGATGTGTGCCAGCCAGGCCCTCCAGCACCCCAAGTGACGCT
CGTGATTGGAACTCACCATCGGCAGGCAGTGGTTCGGTTTAAGAGATGGCATTAGAGGG
AGCCAGTCTGGATGTGGACTTGGATGCCCTGTGGGTATCAGTTCTGCTGACACTTTGGC
CCGAAATAGATCCAGTGTGAGCAAGCAATGTACACCAGAGCCTCAGTGAGCCCATCTGC
ACAGTGGGGAGCATGGAGGGATGGGTTTGGCCTGTGCTTCTGCTTATTCAGTCTTCAGC
TCACGGAAGGGATGCTAGTCCGTGAAGGTGACCTCACAGTACTGGTTAATTAACTTTAT
TGCTCACTGTCAAAAAAAAAAAAAAAAAA

```

Restriction Sites:	Please inquire
ACCN:	NM_004918



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_004918.2, NP_004909.1</u>
RefSeq Size:	1161 bp
RefSeq ORF:	387 bp
Locus ID:	9623
UniProt ID:	<u>O95988</u>
Cytogenetics:	14q32.13
Gene Summary:	Enhances the phosphorylation and activation of AKT1 and AKT2.[UniProtKB/Swiss-Prot Function]