

Product datasheet for **SC320014**

GM CSF Receptor alpha (CSF2RA) (NM_172245) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GM CSF Receptor alpha (CSF2RA) (NM_172245) Human Untagged Clone
Tag:	Tag Free
Symbol:	GM CSF Receptor alpha
Synonyms:	alphaGMR; CD116; CDw116; CSF2R; CSF2RAX; CSF2RAY; CSF2RX; CSF2RY; GM-CSF-R-alpha; GMCSFR; GMCSFR-alpha; GMR; GMR-alpha; SMDP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_172245.1
 GGCACGAGGCGAGAGAAGAAAAGCAGGTGGAAGGAGAGGAAGCGGATGCCGTGGGGTTTA
 CAGCAGGAAAAATCCGTGGAGACAGCAGATCCGAGAAGCGGCGATGTTTGCGTAGAACCTT
 GTACGTGCTTCCCTTCGGCCTGTCGCTCTCCCTTCTCTGACCAGCACCATGCTTCTCC
 TGGTGACAAGCCTTCTGCTCTGTGAGTTACCACACCCAGCATTCTCCTGATCCCAGAGA
 AATCGGATCTGCGAACAGTGGCACCAGCCTCTAGTCTCAATGTGAGGTTTGACTCCAGGA
 CGATGAATTTAAGCTGGGACTGCCAAGAAAACACAACCTTCAGCAAGTGTTCCTTAACTG
 ACAAGAAGAACAGAGTCGTGGAACCCAGGCTCAGTAACAACGAATGTTTCGTGCACATTTT
 GTGAAATTTGTCTGCATGAAGGAGTACATTTGAGGTTACGTGAATACTAGTCAAAGAG
 GATTTCAACAGAACTGCTTTATCCAAATTCAGGAAGGGAGGGTACCGCTGCTCAGAATT
 TCTCTGTTTCATCTACAATGCGGATTTAATGAACTGTACCTGGGCGAGGGGTCCGACGG
 CCCCCGTGACGTCCAGTATTTTTGTACATACGAACTCAAAGAGAAGGAGGGAGATCC
 GGTGTCCTTATTACATACAAGACTCAGGAACCCATGTGGGATGTCACCTGGATAACCTGT
 CAGGATTAACGTCTCGCAATTACTTTCTGGTTAACGGAACCCAGCCGAGAAATTGGCATCC
 AATTCTTTGATTCACTTTTGGACACAAAGAAAATAGAACGATTCAACCCTCCCAGCAATG
 TACCCTGACGTTGCAACACGACGCACTGCCTCGTACGGTGGAAACAGCCAGGACCTATC
 AGAAGCTGTCGTACCTGGACTTTCAGTACCAGCTGGACGTCCACAGAAAGAAATACCCAGC
 CTGGCACGGAAAACCTACTGATTAATGTTTCTGGTATTGGAAAATAGATACAACCTTTT
 CAAGCTCTGAGCCAGAGCAAAACACAGTGTGAAGATCAGAGCTGCAGACGTCCGCATCT
 TGAATTGGAGCTCCTGGAGTGAAGCCATTGAATTTGGTTCTGACGACGGGAACCTCGGCT
 CTGTGTACATTTATGTGCTCCTAATCGTGGGAACCTTGTCTGTGGCATCGTCTCGGCT
 TCCTCTTTAAAAGGTTCCCTTAGGATACAGCGGCTGTTCCCGCCAGTTCACAGATCAAAG
 ACAAACTGAATGATAACCATGAGGTGGAAGACGAGATCATCTGGGAGGAATTCACCCAG
 AGGAAGGGAAAGGCTACCGCGAAGAGGTCTTGACCGTGAAGGAAATTACCTGAGACCCAG
 AGGGTGTAGGAATGGCATGGACATCTCCGCTCCGCGACACGGGGAACTGTTTTCTTGA
 TGATGCTGTGAACCTTTATATCATTTTCTATGTTTTTATTTAAAAACATGACATTTGGGG
 CCAGGCGCGGTGGCTCACGCTGTAATCCAGCACTTTGGGAGGCCAAGGCAGGCGGATC
 ACCTGAGGTGAGGAGTTCAAGACCAGCCTGCCAACATGGTGAACCCCATCTGGACTAA
 AAATGCAGAAATTTACCCAGGCACGGCGCGGACGCCATCATCCAGCTACTTGGGAGG
 CTGAGGCAGGAGAATTGCTTGAACCCGTGAGGCGGAGGTTGTAGTGAGCCAAGATCGCAC
 CATTGCACACCAACCTGGGTGACAGAGCAAGATTGCATCTCAAAACAAACAATAATAATA
 AATAATAAAAACCTGAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_172245

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172245.1 , NP_758448.1
RefSeq Size:	1820 bp
RefSeq ORF:	1203 bp
Locus ID:	1438
UniProt ID:	P15509
Cytogenetics:	X;Y
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway, Pathways in cancer
Gene Summary:	<p>The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 7 encode the same isoform a, which is membrane-bound.</p>