

## Product datasheet for **SC320118**

### SMRP1 (C9orf24) (NM\_032596) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SMRP1 (C9orf24) (NM\_032596) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** SMRP1  
**Synonyms:** bA573M23.4; CBE1; NYD-SP22; SMRP1  
**Mammalian Cell Selection:** None  
**Vector:** [pCMV6-XL5](#)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**Fully Sequenced ORF:** >OriGene sequence for NM\_032596.3  
AATGGGCAGTGCCCATCCTGGAGGGGCATCAGTGTGGGTGGGAGCAGGAAGGCTGTTG  
GGGAAGGAAGCTGAGAGATCCTCTGCAGGAGGGGATTATAGCAGCTAAGGGCACCTGCCA  
GGCTGGAGGTGTGGCCTGTCTAGGGGAAGGGTGGTCACTGATGTCCAGCAGGAGCCCTTG  
ACAGGTAGAGCCACCACGGCATCTCCAAGACCACCCATGTTCTGTTCTCCCCTAAGA  
CCAGGACCCCTATCAGTACCTACAGTACTCCTACAGGGCTCCCACCTCCATCAAGGAGG  
TCTATAAGGACCCACCCCTGTGTGCCTGGGAAGCCAACAAGTTTCTGACTCCGGGTCTGA  
CTCATAACCATGGAGCGACATGTGGATCCCAGGAGCCCTGCAGAAGATGGCCAAATGTGCTG  
TACAGGACTACACTTATAGGGGGTCCATATCAGGCCACCCCTACTTGCCTGAGAAGTACT  
GGCTTTCTCAAGAGGAAGCAGACAAATGCAGCCAAACTACCTGGGCAGTACTGGTACA  
ACACATGGAGGATGGAACCTTACAACAGCAGCTGCTGCAACAAGTATACCACCTACCTTC  
CTCGGCTGCCTAAGGAGGCCAGGATGGAGACAGCAGTTCGAGGAATGCCCTTGGAATGCC  
CTCCTAGGCCGGAGCGGCTCAATGCCTACGAGCGCGAAGTGATGGTGAACATGCTGAACT  
CACTGTGCGGGAACCAGCAGCTGCCGCGGATCACGCCCGATGCGGGTGCCTGGACCCGC  
TGCCCGGCGCCTGCCCTTCCATGGTTACGAAAGTGCTTGTGCTCGGGCCGCCACTACTGTC  
TGCGCGGGATGGACTACTACGCCAGCGGGGCGCCCTGCACCGACCGCCGCTGCGGCCTT  
GGTGCCGGGAGCAACAGACTATGTGTACCTCCCTACGAGCACCAGCCCGGAATGCAGTGT  
GCTGTTACAACCTCCCCGCGTCATACTACCCATATCCGAACCTTAGATGGGACACAAGT  
CACTTCAAGAAGTCTGGTGGTCCCCAGAGAAACAACATGTTATCCATCCTGAGTTTGTG  
TCTGAGACCTATCCCGACTATCGTTGCTGGTAGAGCCTGGGCTTCTGTGCAAGGGGGC  
AGAGCAATAAATAAACTCTTCAACCCCAAAAAAAAAAAAAAAAAAAAAA

**Restriction Sites:** Please inquire  
**ACCN:** NM\_032596



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_032596.3</a> , <a href="#">NP_115985.2</a>
<b>RefSeq Size:</b>	1184 bp
<b>RefSeq ORF:</b>	789 bp
<b>Locus ID:</b>	84688
<b>UniProt ID:</b>	<a href="#">Q8NCR6</a>
<b>Cytogenetics:</b>	9p13.3
<b>Gene Summary:</b>	<p>This gene encodes a nuclear- or perinuclear-localized protein with no predicted domains or similarity to other known proteins. Expression of this gene is induced during the differentiation of bronchial epithelial cells, and the encoded protein may play a role in ciliogenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>