

Product datasheet for **SC320237**

RTL8A (NM_001078172) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	RTL8A
Synonyms:	CXX1b; FAM127B; MAR8A; SIRH6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_001078172.1
 CCTTAGCGCTGACCCACGCAACCCAGCGAACTCCGCGGAGCGCGCGGCACGATGG
 ACGGTGCGGTGCAGCTGATGAAGGCCCTCCTGGCCGGGCCCTCCGGCCCGCGCGCGT
 GCTGGAGGAACCCGATTCCCTTTCCGAGACGTTTGACGGAGATACCGACCGACTCCCGG
 AGTTTCATCGTGACAGCTGCTCCTACATGTTCTGGACGAGAACACGTTCTCCAACGACG
 CCCTGAAGGTGACGTTCTCATACCCGCTCACGGGGCCAGCCCTGCAGTGGGTGATCC
 CCTACATCAGGAAGGAGAGCCCTGCTCAATGATTACCGGGGCTTCTGGCCGAGATGA
 AGCGGGTCTTTGGATGGGAGGAGGACGAGGACTTCTAGGCCGGGAGACCCTTGGGCCTGG
 GGGCGGTGCTCTGGGAAGAGTTGCTGTGCCAGTGGCCACCGCTAGGGTCTCCACAGGC
 GCCCTCCCTCCGCGCTCCCTCCCCCTCGAGCCGCGCATGTCCCTGCGCTCCTGTCC
 CCTCCCGCTAGTGCTTGCTTTGTTCCAGGAATAGCGCTCCAGGCTCCTGCTGCAGCTC
 CTGGGCCGCACTCTAGAGCGCTGCTCCGCCCTCTTTGCCAGCCAGCCCTCCATG
 CACATTTGGACGCTGCTGCACTTCAGGTGCAAGCTGGGCTCCTGTTACATACTGGACA
 GACCACCCACAGCGCTGCTGCCAACCCACCTCCACTCCTCCAGACTGCCAGCCAACCTAC
 ATCATTCTGCCACAGACCTACGCTACCACTGCCATCGCCATCCACCGCATCTACCAAC
 AGACTGTGGCTCCTAGTGATGTGGAATCACCTCGGAGGTATCCGAGCTGGACACAGCCCC
 TGGACAGTGATCCAGACAGCTGGCCGTCCCCAAGGAGCCTGTACCTTCAGCGAGACCC
 ATTTCTCCCCATCCCCAGAGACCTTTGTGTTCTGCCACATAGCTGCCAGGGCTTAAG
 TGTGCTGGCAACCAATCGAATCTCTATTTTCTCCTGTGGACGAGTTAGTTTGGCTA
 GAATCTGTTTCTCTGAATTTGCAGCTGTCTCTGATGGGTGCCTTTTGTCAACAC
 AGTAAGCCCTGCTCCCTCCCTGCTCTAATACTACCTGTACAAAGGTTTTTCTTAT
 TTTAATAAATGTCAGACACTATTAATAAGAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_001078172



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001078172.1</u> , <u>NP_001071640.1</u>
RefSeq Size:	1275 bp
RefSeq ORF:	342 bp
Locus ID:	26071
UniProt ID:	<u>Q9BWD3</u>
Cytogenetics:	Xq26.3