

Product datasheet for **SC320426**

RAB24 (NM_001031677) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAB24 (NM_001031677) Human Untagged Clone
Tag:	Tag Free
Symbol:	RAB24
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_001031677.2

```
GCGGGAGTCGCGCGCTGCGGGTAGGAGCCGGTTGCGGGAGACCCAGGTTGCGTTGGG
ATCCCCAGCCAGAACGGAGCTTAAGCCGGCAGCGAGCGAATGACGGAGTAGCGAGCTG
CACGGCGCGCTGCTGCGCTGTTGAGGACCTGTCCCGCGCTCCAGGCCGCCCGAGG
CTTGGGGTCTTCGAAGGATAATCGGCGCCCGGGGCCGAACAGCGGGGCACACGGGGCGC
TGCCGAAGTGCAAGGCCACGGCCAGAGCTCGAGCCGACGCGCTGTCTGGAGTCGTAGGT
TGGCGCCGTTTGGGGTCGGGGTCTGAGGCTTGGGCGCTGCCTGGGCCGAGCGGAGATCGG
GGTTTGCCTCCCGTCCCGCTCAGGACCCTGACGTGGCTGAAGCGGCCCGGGGAGCATGA
GCGGGCAGCGCTGGACGTCAAGGTGGTGTGCTGGGCAAGGAGTACGTGGGCAAGACTA
GCCTGGTGGAGCGCTACGTGCACGACCGCTTTCTGGTGGGGCCTTATCAGAACCACATCC
GGGCCCGCTTCGTGGCCAAGGTGATGTCGGTCCGAGACCGGACTGTGACATTAGGTATT
GGGACACAGCAGGCTCTGAGCGCTATGAGGCCATGAGTAGAATCTACTATCGGGGTGCCA
AGGCTGCCATCGTCTGCTATGACCTCACAGACAGCAGCAGCTTTGAGCGAGCAAAGTTCT
GGGTGAAGGAACTGCGCAGCCTAGAGGAGGGCTGCCAAATCTACTTATGTGGACCAAGA
GTGACCTGCTGGAAGAAGACCGGAGGCGTCGACGTGTGGACTTCCACGACGTCCAGGACT
ATGCAGACAATATCAAAGCTCAGCTCTTTGAAACATCCAGCAAGACAGGCCAGAGTGTGG
ACGAGCTCTTCCAGAAAGTGGCAGAGGATTACGTCAGTGTGGCTGCCTTCCAGGTGATGA
CAGAGGACAAGGGCGTGGATCTGGGCCAGAAGCCAAACCCTACTTCTACAGCTGTTGTC
ATCACTGAGTCAGCACTCACCTGGCCTGGGGGAATTAAGGAATTCCTCCGTAAGGGCTGG
ACCCAGCTCCTTTCTGGGCTTGGGTAGTCAAATGTCTGAGCTACCCAGGTCTCATGTC
AGCAGAGTGGCGCCTGCCTGTGCTGGCCATGGAACGGAGACAGCATTGGGCTGACTGTG
GGCATGAGGAGGGATAAGGCTGATTTGGACCCAGGCTTCTGCCCTGGACAGCACTTGTG
TCTGCAGATTATTTAAGTGGCTTTTGTCTGTAATAAAATCAGTGCAGTGTGAATCACA
CCCAGCCCTTTCCCTGCTGTGTGGATTAGGTGTCAAGACACCTAGTTCTTTCTGGGGCC
ACCCGGCTGGCCTCACTGCTTATATTAAGGCTCCTCCCACTCTCATTTTCTTTGGAAA
ACAAGACTTTTTTCCCCATGGTTACCGCTGAGATACTGGGGCTGTAGTAGTATAAAAGCT
CACAGTTCCTTCTGAGTGTGAAAAGAGTGCATGAGTTGCTTCGAAATAAAAGGGTCAAG
CATTAAAAAAAAAAAAAAAAAAAA
```



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001031677
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001031677.2</u> , <u>NP_001026847.1</u>
RefSeq Size:	1612 bp
RefSeq ORF:	612 bp
Locus ID:	53917
UniProt ID:	<u>Q969Q5</u>
Cytogenetics:	5q35.3
Protein Families:	Druggable Genome
Gene Summary:	<p>RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Olikkonen et al., 1993 [PubMed 8126105]).[supplied by OMIM, Aug 2009]</p> <p>Transcript Variant: This variant (1) encodes the functional protein. Both variants 1 and 2 encode the same protein.</p>