

Product datasheet for **SC320474**

C10orf84 (FAM204A) (NM_022063) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C10orf84 (FAM204A) (NM_022063) Human Untagged Clone
Tag:	Tag Free
Symbol:	C10orf84
Synonyms:	bA319I23.1; C10orf84
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_022063.1
 CGGGGTCTTAAGATGGCTGCTGGGCGACCACTTCTGACAGCTCCAGTTAATGGGACTGC
 CAGCAGGCTAGACATCCTTCGCCGGAATGTGCAGCGCTTAGTTCTGTGCCCGGAGGGG
 TTCTCTCTTTTTGGGATTGGAGTGTGGCATTGCCGCCAGAAAGCCTGCAGAATTCGGG
 TTGTGCAGCCAAGCAGGAACGGTGTACCTTTCTAGACCCATCCTCTTGAATCTTGGACC
 CTAAGAAAAGATGTGGAGTGGGCTGCTACCTCCTGGCCTAAATGAAAGTGACGCTGAGTC
 AAACCGGAAGATGAAGCTACGTTGGAGAACTCTGGACTTAACTTACAGGAAGATAAAGA
 GGATGAGAGCATCAGAAAAACAGAAATCATAGATTTCTCAACAGATGAACAAAACTGA
 AACAGAGTCAAAATGTAATGCCTATGAAGAGTGCCTTCTGGAATCCCATAGATATGTG
 GAATAAATTTCAAGAATGCATAAAAAACATTCTGAACAGAAAAGCACAACTCAAGATT
 CAGAGGAAAAGAAGAAAACGCTCCAGAAAAGATAAATTGAAGAATGAAAAGAATTACA
 TAGTGAACCGTCTCAAATGAAACCCAGTGGAAAGAGCTTACTCAGTATTTGGAGTCAA
 TGATAGATTTGACCCGCCTGTTAAAAGGAAAAAGTTGAGAAGTCAGGCCTTGAAAAGAG
 GATAGACCAGGCTGTGGAGGAGTGAATATTGAGAAGGCTGAGGAACTCAGCAACCAAGCT
 AGCTACTCGAGAGCTTGGTGTAAAAATTGCCAAAGCAGTTGCCTGCCACAACCTTTGTAA
 AGCCAAAAGGAGGTTGAAAATTCACAGGCTGCCCGAAAAAGAAAGAACTTGCATGGGG
 GTTTGAAGCAAAGAAGAGATGGGAAACCAAAAGCAACATGGGATACATGTAACCTGCCAG
 AGTGCTTCAAGACATTTGTAGACTCAAAATGCTCAATTTACTGAGAATGTTTTCTGCCTA
 TGTTAATATGTCAGAAAATTGATAGCACTAAAACAAAAATAAGCATAAAAATTTGGGATTT
 TGATTATCAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_022063



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022063.1 , NP_071346.1
RefSeq Size:	2746 bp
RefSeq ORF:	702 bp
Locus ID:	63877
UniProt ID:	Q9H8W3
Cytogenetics:	10q26.11