

Product datasheet for **SC320506**

Calneuron 1 (CALN1) (NM_031468) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calneuron 1 (CALN1) (NM_031468) Human Untagged Clone
Tag:	Tag Free
Symbol:	Calneuron 1
Synonyms:	CABP8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_031468.2
GACAAGGACAGGGAGGCACGCGCCAGATCGCGCTCCGAGTCGACGAGGGGGCTCGTGCA
GGTCCC CGCGCGCGCTGGGGAGCGGCCCGGGGGCGCCGAGGGAGCCGGCGCAGCGAGG
ATGCCGGCACGGGACGCGGCTGCCTGCTGGGGGTGTGAGCGAGAGGCTGAGCGCTGCC
GCTCGGGCATCCATTCCGCACCGCCCCCTCCCTGCGGGCCTCGGAGGAAGCCCCCGCG
CTGTGCGGAGGCGCTCGGCTGCGGGGTGCGGAGCCCCGGCCAGCAAGAGGAAAAGA
TGCCGTTCCACCATGTGACCGCGGCTTGTGTACAAGGGGAATTACCTCAACCGATCGC
TCTCTGCTGGCAGTGACAGCGAACAGCTGGCTAATATCTCCGTGGAGGAGCTCGATGAAA
TCCGAGAGGCCTTTCGGGTTCTGGACCGGGATGGGAACGGCTTCATCTCCAAGCAGGAGC
TGGGCATGGCCATGCGCTCTTTGGGTACATGCCAAGCGAGGTGGAGCTGGCCATCATCA
TGCAGCGCTTGGACATGGACGGGATGGCCAGGTGGATTTTGATGAATTCATGACCATTC
TTGGCCCCAAACTGGTGTCTTCAGAAGGTCGCGATGGTTTTCTTGGGAACACGATAGACA
GCATATTTCTGGCAGTTTGACATGCAAAGGATAACTCTGGAAGAGTTGAAGCACATTCTCT
ATCATGCCTTCCGAGACCACCTAACGATGAAGGACATTGAGAATCATTATCAATGAGG
AAGAGAGCCTGAATGAGACCTCGGGAACTGCCAAACAGAGTTTGAAGGAGTGCATTCCC
AGAAGCAGAACAGACAGACCTGCGTCCGGAAGAGCCTCATATGCCCTTTGCTATGGCCT
TCATCATCAGTGTGCTGATTGCAGCCAACAGATACTCCGGAGCGGCATGGAGTAGC
AGCCTCCCAGCCACACTGTGTTGCAGCTCACCGCGCATGTGCATGCCGTGCGGGCAG
ACTTTCTCCACACAGCAGACCGGACCTATGGACTATGGATGGATGCGGACGATGGAA
CGTTTCAGTAAAGAACCACAGGTTGAGTGCAACTGGGGTTGCACTTCAGCTTCATCTCCT
TGGCAGGGACACTAAAGGCTGTCTTCAGTGCTTAAATGGTTTTGTTTTCTAATGCAATA
AATCGCCAGGATTTCCAAAAAAAAAAAAAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_031468



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031468.2</u> , <u>NP_113656.1</u>
RefSeq Size:	5983 bp
RefSeq ORF:	660 bp
Locus ID:	83698
UniProt ID:	<u>Q9BXU9</u>
Cytogenetics:	7q11.22
Domains:	EFh
Protein Families:	Transmembrane
Gene Summary:	<p>This gene encodes a protein with high similarity to the calcium-binding proteins of the calmodulin family. The encoded protein contains two EF-hand domains and potential calcium-binding sites. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript.</p>