

Product datasheet for **SC320584**

MRPL28 (NM_006428) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPL28 (NM_006428) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRPL28
Synonyms:	MAAT1; p15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_006428.3

```

CGGTTCCGGTGGGCTCTGAACCTGAAAGGCTCGCGATGCCTCTACACAAGTATCCCGTG
TGGCTCTGGAAGCGGCTGCAGCTGCGGGAGGGCATCTGTTCCCGCCTGCCCGGCCACTAC
CTGCGCTCCCTGGAGGAGGAGCGGACGCCCACTCCCGTGCACTATAGGCCTCATGGGGCC
AAGTTCAAGATCAACCCAAGAACGGGAGCGGGAGCGTGTGGAGGACGTGCCCATTTCCC
ATCTACTTTCCCCCGAATCCCAGCGGGGTTGTGGGGCGGCGAGGGCTGGATCCTGGGC
CAAATATATGCCAACAACGACAAGCTCTCCAAGAGGCTGAAGAAAGTGTGGAAGCCACAG
CTGTTTGAGCGAGAGTTCTACAGTGAGATCCTGGACAAGAAGTTCACAGTGACTGTGACC
ATGCGGACCTGGACCTCATCGATGAGGCTTACGGGCTCGACTTTTACATCCTCAAGACC
CCGAAGGAGGACCTGTGCTCCAAGTTTGGGATGGACCTGAAGCGAGGGATGCTGCTGCGG
CTTGCCCGGCGAGGACCCAGCTGCACCCGAGGACCCGAGCGGCGGGCAGCCATCTAC
GACAAGTACAAGGAATTTGCCATCCCAGAGGAGGAGGCAGAGTGGGTGGGCCTCACGCTG
GAGGAGGCCATTGAGAAGCAGAGACTTTTGGAGGAGAAGGACCCCTGTACCCCTGTTCAAG
ATCTATGTGGCGGAGCTGATCCAGCAGCTGCAGCAGCAGGCACTGTCAGAGCCGGCGGTG
GTGCAGAAGAGAGCCAGTGGCCAGTGACCACACAGCTCCTCCATGCCGTGACCAACAGGCC
CAGCTTTCCCTGCCAGGCCCTTTGCACTGAGGACACAGATCCCGGGGAGCTGTGAGGGCC
ACCGGTGGGCGAGTGGGTGGATCCTGGTTTCGTGTGCTGCCCATGCACCTTCCAGCCCGGG
GCCAGCTTGGCAGGGATCCCCAGGAGGCTGGGCCGCCAGAGGCTCCTCTCAGGCTGGG
CCCCGACGTTTGGCGAGTGTTCTTGTGCCCGTGGGGCCGGGAGCGAGTAAAGTCTGGG
CCAGGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites:	Please inquire
ACCN:	NM_006428


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_006428.3 , NP_006419.2
RefSeq Size:	1127 bp
RefSeq ORF:	771 bp
Locus ID:	10573
UniProt ID:	Q13084
Cytogenetics:	16p13.3
Gene Summary:	<p>Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein, a part of which was originally isolated by its ability to recognize tyrosinase in an HLA-A24-restricted fashion. [provided by RefSeq, Jul 2008]</p>