

Product datasheet for **SC320618**

COX4NB (EMC8) (NM_006067) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	COX4NB (EMC8) (NM_006067) Human Untagged Clone
Tag:	Tag Free
Symbol:	COX4NB
Synonyms:	C16orf2; C16orf4; COX4NB; FAM158B; NOC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_006067.3
 CTTCCCAGAGGCGCGAGTCGCTTTTCCAGGAGCCGAGGCCAAACGGGGCCGGCCCTGAA
 GGGTCCCCCGCCAGCCCGCGGCCGCGCCATCAATCGCCGCGCCCTCGTCCCCTTCTCGG
 CTGAGGCGCCGCGCGGCCAGGCAGCGGGTCCAGGCCTCAGCCGCGCGCCAGGGGCTCC
 GGGGCCCTCCCGGGTCAAGTACCCGCGGTGAACTGACCACCCAGGCCTACTGCAAGAT
 GGTGCTGCACGGCGCAAGTACCCGCACTGCGCCGTCAACGGGCTCCTGGTGGCCGAGAA
 GCAGAAGCCGCGTAAGGAGCACCTCCCCCTGGGCGGCCCGCGGCCACACACCCTCTT
 CGTGGACTGCATCCCCCTTTCCAGGCACCCTGGCCCTCGCCCCATGCTGGAGGTGGC
 TCTCACCTGATTGATTCATGGTGCAAAGATCATAGCTACGTGATTGCTGGTTATTATCA
 AGCTAATGAGCGAGTAAAGGATGCCAGTCAAACCAGGTTGCAGAGAAGGTGGCCTCCAG
 AATCGCCGAGGGCTTCAGCGACTGCGCTCATCATGGTAGACAACACCAAGTTTACGAT
 GGACTGCGTAGCGCTACGATCCACGTGTACGAGCACCATGAGAACAGATGGCGGTGCAG
 AGACCCACACCATGACTACTGTGAAGACTGGCCAGAGGCACAGAGGATCTCAGCCTCGT
 CCTGGACAGCCGGTCTACGAGACGCTCGTGGATTTGATAACCACCTGGATGACATTTCG
 GAATGACTGGACAAACCCAGAGATCAATAAAGCTGTCCTACACTTGTGCTAGGCAGGCAC
 CGCTGTGACTGGGCTCCGGGCTTTCCACTACGTTGAAGAAGAAAACCTATTTTTAAATG
 TAAATAAAATATCTGGTAGCCTGTGTGGAAAGCTGACCGTTTTAAGAAGTGGCATGTGCC
 TTGAAAGGGGGCAGAATGTTTCAGTCGGTCGTGTTTTAACACAGAGTCTCTAGAAGAGGT
 GCAGACATCCCGTCTGACTGTCCCTGTGGACTCTCTCAGTTGTATGTTGTCTATAATCCTC
 CAAATCAAAGCTCTTCTGCTTGTGCAAGATTGTTCTATTAACAGTTTTAACTAACCT
 TTTATAATCTGGAGAGAATACTTTTAAAGGCTTTTGAACATCTTTTTTAAAAAAAAAAAA
 AAAAGCCACATCTTATACTCCAATGAGTGTGTTTGTGAGTTTTATCCTGTTTGAAATTT
 TGATGCTGTCCAGCGCTCCTGGGATCAGGGTAGACACGTGGAGTGTGGTGTGAGCTGT
 TTGAGAAGCTCTGTGTGGGTGTTTCTGTGGTCAGGTAGGCTGGGCTGTGAACAGCACTGC
 TCCACATGGTGTGCGCTCACAGAGCCACCGCAACTCTGCCAGCGGGAGTGACACCACTG
 ATGCCACAGAAGCCTTTATCGTTAGGGGCAAAGCAGACATCTCTGGTAGCAGGAGCTCA
 GTGTGGTGTCTTTCAGCCACGCAGAAAGCTCACGTGCCTGCAGGGCCAGGAAAGGGAAAT
 GATCGTTTCTTTGTATTAGGAAAAACAAGACTTAGACTTCCCGAGTTGTGCAGAGTGAC
 GTCTGTACTTTTTATACTTAGAACATTTGGTCAAGTTTTGTAACCAAACCTTAGGGAATT
 GAGAGGAAACTTTATGGTACGTTTACCTTTTTTCTTTTTATATGAATTATTTTGTACA
 AAGTTCAGAGAGCGGTACCTTTCTACTGTCTTTTCTTTGACTCCATGGAGGAATTTGTA
 ACGATCGGTGTGACCTGAAAGGGTTTGTGATACAGCGGTCACTATTTTCTTTAACACAT
 AGGAAAAAAGTTCAGACATCGCTACAATAAATGTGCTTTATTTGATCCAAAAA
 AAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_006067

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_006067.3](#), [NP_006058.1](#)

RefSeq Size: 1941 bp

RefSeq ORF: 633 bp

Locus ID: 10328

UniProt ID: [O43402](#)

Cytogenetics: 16q24.1

Domains: UPF0172

Gene Summary: Part of the endoplasmic reticulum membrane protein complex (EMC) that enables the energy-independent insertion into endoplasmic reticulum membranes of newly synthesized membrane proteins (PubMed:30415835, PubMed:29809151, PubMed:29242231, PubMed:32459176, PubMed:32439656). Preferentially accommodates proteins with transmembrane domains that are weakly hydrophobic or contain destabilizing features such as charged and aromatic residues (PubMed:30415835, PubMed:29809151, PubMed:29242231). Involved in the cotranslational insertion of multi-pass membrane proteins in which stop-transfer membrane-anchor sequences become ER membrane spanning helices (PubMed:30415835, PubMed:29809151). It is also required for the post-translational insertion of tail-anchored/TA proteins in endoplasmic reticulum membranes (PubMed:29809151, PubMed:29242231). By mediating the proper cotranslational insertion of N-terminal transmembrane domains in an N-exo topology, with translocated N-terminus in the lumen of the ER, controls the topology of multi-pass membrane proteins like the G protein-coupled receptors (PubMed:30415835). By regulating the insertion of various proteins in membranes, it is indirectly involved in many cellular processes (Probable).
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.