

Product datasheet for SC320657

LOC284912 (NM_203375) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC284912 (NM_203375) Human Untagged Clone
Tag:	Tag Free
Symbol:	LOC284912
Synonyms:	MGC3170
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_203375.1

```

AATTTGCTCTCTGCCTCACTGCCTGAGCTGGAACATTTGTCTTCTGCTCCTCAGACTGGGC
TTTACCCCATCTACTCCCAGGATACTCAGACCTTCAGACTCTGGCTGAATTATAAACTG
ACTTTCTGCGGTCTCCAGCTTGCAAATGGCAGGCTATGGAACCTTCAGCCACAAATCAC
CCTTCAGGGTAGACCGAGGAGAGCACCTCCCCTTCTGCTGAAGGGAGCCGATACACGC
TGGTGCCGGCTGGCCAAGAAGGAGGTGGACAAACCTGGCGACCGCTCCCTGGCACCACCT
CACCAGAGGACAGCTCACACCTCCACAGGGGTGAGAAGAAAGGAGCTGCCCTGGGGATGG
CTGTGTATCCCTACAACCTGCCAGGCACATGCCCTTGAACACCCGATGCTACGTGTCCCA
AAGGAACTGGTCTCCACCCACCCCGGCCGCTCCTGCTCCTGGGTACCCACCTTAGTA
AATGGCGCCACCCTGCGCGGTCACTCAGCGAGAACTCAACTCCTGGCAGCAGATGAC
GGGCACTCTGGTTAAATGACTCTCTCCAGCCTCCAAGTTCAACCTGCAGGGAAGCCCTGG
AAATCCTGTCTCCCTCTGCCCTGCCTCTCCTTCAGCTCCGAACCCACCCCACTCAGTAC
ATGGCCCCCAACAGGCTTTCAAATGTCTGCCTTCCCCAAACCAATCTGTGGTTTTAAAA
GAAGTAAGGTGAGGTACAGTAGCTCATGCCTGTAATCCCAGCACTTTGGGAGGCCGAGGC
AGGCGGATCACCTGAGGTGAGGAGTTTGAGACCAGTCTGGCCAACATGGCGAAACCCCGT
CTCTACTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites: Please inquire

ACCN: NM_203375

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_203375.1, NP_976309.1</u>
RefSeq Size:	899 bp
RefSeq ORF:	243 bp
Locus ID:	284912
Cytogenetics:	22q12.3