

Product datasheet for **SC320678**

Lano (LRRC1) (NM_018214) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lano (LRRC1) (NM_018214) Human Untagged Clone
Tag:	Tag Free
Symbol:	Lano
Synonyms:	dj523E19.1; LANO
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_018214.4
 GGAAGCTCCGCGCGGCGGGGCGGGGCGGCGACGGCGACTGGCGGGTGGGAGTGGAGGCACC
 GGCTGGCGGGGCGGGGTACAGGGACGGGGCAGGGGCTCCCGCTCCAGGTTCTTGAAGCA
 CTTCCGACCGGAAGCCCGGCGGAGAGCGAGCTAACCAAGAGCCAACAACGAGCGCG
 GAGAGGGCAGCGGACTGAGCGGAGCCGCGCCAGAGCGGGCTCGGAGCCCGGGTCTCCG
 CCGCTCGGACCCGGCTAGGCGGCGGGGCGGGCGGATGTTCCACTGCATCCCCCTGTG
 CGGTGCAACCGTCAATGTGGAGAGCATCGACAAGCGCCACTGCTCGTGGTCTACGTCCC
 CGAGGAGATCTACCGCTATGCCCGGAGCCTGGAGGAGCTGCTGCTGGACGCCAACAGCT
 CCGCGAGCTGCCCGAGCAATTTTTCCAGCTAGTCAAATTACGAAAAGCTTGGACTTAGTGA
 TAATGAAATTCAGCGGCTCCCTCCAGAAATAGCAAATTCATGCAGCTGGTGGAACTAGA
 TGTGTCTCGAAATGAGATTCTGAAATCCAGAAAGCATTTCATTCTGTAAGCACTGCA
 GGTAGCTGACTTCAGCGGAAACCACTGACTAGGTTGCCAGAAAGCTTTCCTGAATTACA
 GAATTTAATGTCTTTCTGTAATGACATCTCACTACAGTCTCTACCTGAAAATATTGG
 CAATCTTTATAACCTGGCTTCACTGGAAGTGAAGAGAAATCTTACATATCTTCTGTA
 CTCTCTTACCCAGCTGCGAAGACTAGAAGAACTTGATTTAGGAAACAATGAAATATATAA
 TTTGCCAGAATCAATTGGAGCCCTTACATCTAAAAGATCTCTGGTTGGATGGAATCA
 ACTGTCAGAATTACCTCAGGAAATAGGAAATCTGAAGAACCTGCTGTGTTAGATGTCTC
 TGAAAACAGGTTGAAAGACTTCTGAAGAAATCAGTGGCCTGACTTCATTAACGGATTT
 AGTCATTTCCAGAACTATTAGAAACGATTCCGGATGGCATTGGAAAATAAAGAACT
 GTCAATCTTGAAGGTGGATCAGAATAGACTCACACAGTTGCCTGAAGCAGTTGGGGAATG
 TGAAAGTCTCACTGAGTTAGTCTTACAGAAAATCAGCTCCTGACCCTGCCTAAAAGCAT
 AAAAGAGATCGGCGGGTCTGCAGCCTCACTGTGTTCTGTGTACGTGACAACAGACTAAC
 TCGGATACCTGCAGAGGTGTACAGGCAACAGAACTTCATGTCCTGGATGTGGCAGGAA
 CAGGTTGCTGCATCTACCTTTATCCCTGACTGCCTTGAAGTTGAAGGCTCTGTGGCTATC
 TGACAACCACTCCAGCCCTGCTTACATTCCAGACAGACACAGACTACCCACAGGAGA
 GAAGATTTAACCTGTGCTTACTTCTCAGCTGCCTTCTGAACCTACTTGTCAAGAGAA
 TCTGCCTCGTGTGGTGCAGTGGAGAACTTGGTAAATGATGTCTCTGATGAAGCCTGGAA
 CGAGCGTGTGTCAACAGAGTCACTGCGATCCGATTTGTGGAGGATGAGAAAGATGAAGA
 AGACAATGAGACGAGAACAATTAAGGCGAGCCACTCCACACCCAGGGGAGTTAAAGCA
 CATGAAAAGACAGTGGAGAATTTACGGAATGACATGAATGCTGCTAAAGGACTGGACTC
 AAACAAAACAGGTCATCATGCCATTGACCGAGTGACCACTTCTGTGTAGAGTTTCAC
 CTCCAAGTTTTACCTCCTGTGTCTTCTCTGCTGTGCGAGACGTTCTGTCTGCTTCCCGG
 GAGCCTCACGTGCTCCTTGTCTAACCCAGCCCCGCGCCATCTTCCCGTGGAGTGTGG
 GGAAGTGTGTCTCCAGGAAGTGCCTTACTCATCCCGCAACAGTCAAGCACCAGTGTG
 GTCTCCCGGTGTGATTTTTTTTTTTTTTAAATTTTCAGTTGTTTGAATAAGTAGAATACAC
 TACTGTCAACATAA

Restriction Sites: Please inquire

ACCN: NM_018214

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018214.4](#), [NP_060684.4](#)

RefSeq Size: 3180 bp

RefSeq ORF: 1575 bp

Locus ID: 55227

UniProt ID: [Q9BTT6](#)

Cytogenetics: 6p12.1

Domains: LRR, LRR_TYP, LRR_PS