

## Product datasheet for **SC320864**

### ATP5A (ATP5A1) (NM\_001001937) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ATP5A (ATP5A1) (NM_001001937) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATP5A
Synonyms:	ATP5A; ATP5A1; ATP5AL2; ATPM; COXPD22; hATP1; HEL-S-123m; MC5DN4; MOM2; OMR; ORM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >OriGene sequence for NM\_001001937 edited  
 GTCAGTCCGGAGGCTGCGGCTGCAGAAGTACCGCCTGCGGAGTAACTGCAAAGATGCTGT  
 CCGTGCGCGTGTGCTGCGGCCGTGGTCCGCGCCCTTCTCGGCGGGCCGGACTGGTCTCCA  
 GAAATGCTTTGGGTTTCATCTTTTCATTGCTGCAAGGAACTCCATGCCTCTAACACTCATC  
 TTCAAAAGACTGGGACTGCTGAGATGCCTCTATTCTTGAAGAGCGTATTCTTGGAGCTG  
 ATACCTCTGTTGATCTTGAAGAACTGGCGTGTCTTAAGTATTGGTGATGGTATTGCC  
 GCGTACATGGGCTGAGGAATGTTCAAGCAGAAGAAATGGTAGAGTTTTCTCAGGCTTAA  
 AGGGTATGTCCTTGAACCTGGAACCTGACAATGTTGGTGTGTCGTGTTGGAAATGATA  
 AACTAATTAAGGAAGGAGATATAGTGAAGAGGACAGGAGCCATTGTGGACGTTCCAGTTG  
 GTGAGGAGCTGTTGGTGTGATGTTGATGCCCTTGGTAAATGCTATTGATGGAAGGGTC  
 CAATTGGTTCCAAGACGCGTAGGCGAGTTGGTCTGAAAGCCCCGGTATCATTCTCGAA  
 TTTCACTGCGGGAACCAATGCAGACTGGCATTAAAGGCTGTGGATAGCTTGGTGCCAAATG  
 GTCGTGGTCAGCGTGAAGTATTGTTGACCGACAGACTGGGAAAACCTCAATTGCTA  
 TTGACACAATCATTAAACCAGAAACGTTTCAATGATGGATCTGATGAAAAGAAGAAGCTGT  
 ACTGTATTTATGTTGCTATTGGTCAAAAGAGATCCACTGTGCCAGTTGGTGAAGAGAC  
 TTACAGATGCAGATGCCATGAAGTACACCATTTGGTGTGCGGCTACGGCCTCGGATGCTG  
 CCCCACCTCAGTACCTGGCTCCTTACTCTGGCTGTTCCATGGGAGAGATTTTAGAGACA  
 ATGGCAAACATGCTTTGATCATCTATGACGACTTATCCAAACAGGCTGTTGCTTACCGTC  
 AGATGTCTCTGTTGCTCCGCCGACCCCTGGTGTGAGGCCTATCCTGGTGTGTTCT  
 ACCTACACTCCCAGTGTGGAGAGAGCAGCCAAAATGAACGATGCTTTTGGTGGTGGCT  
 CCTTGACTGCTTTGCCAGTCATAGAAACACAGGCTGGTGTGATGTCTGCTTACATTCCAA  
 CAAATGTCATTTCCATCACTGACGGACAGATCTTCTTGGAAACAGAATTGTTCTACAAG  
 GTATCCGCCCTGCAATTAACGTTGGTCTGTCTGTATCTCGTGTGCGGATCCGCTGCCAAA  
 CCAGGGCTATGAAGCAGGTAGCAGGTACCATGAAGCTGGAATTGGCTCAGTATCGTGAGG  
 TTGCTGCTTTTGGCCAGTTCGGTCTGACCTCGATGCTGCCACTCAACAATTTTGAAGT  
 GTGGCGTGCCTAACTGAGTTGCTGAAGCAAGGACAGTATTCTCCATGGCTATTGAAG  
 AACAAAGTGGCTGTATCTATGCGGGTGAAGGGGATATCTTGATAAACTGGAGCCAGCA  
 AGATTACAAAGTTTGAAGTGTCTTTCTTGTCTCATGTCGTGAGCCAGCACCAAGCCTTGT  
 TGGGCACTATCAGGGCTGATGGAAGATCTCAGAAACAATCAGATGCAAAGCTGAAAGAGA  
 TTGTAACAAATTTCTTGGCTGGATTTGAAGCTTAACTCCTGTGGATTACATCAATAC  
 CAGTTCAGTTTTGTCATTGTTCTAGTAAATTAGTCCATTTGTAAGGGTTACTCTCAT  
 ACTCCTTATGTACAGAAATCACATGAAAAATAAAGGTTCCATAATGCAAAAAAAAAAAAAA  
 AAAAAAAAAAAAAAAAAAAAAAAAAA

- Restriction Sites:** Please inquire
- ACCN:** NM\_001001937
- Insert Size:** 1900 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001001937.1</a> , <a href="#">NP_001001937.1</a>
<b>RefSeq Size:</b>	1950 bp
<b>RefSeq ORF:</b>	1662 bp
<b>Locus ID:</b>	498
<b>UniProt ID:</b>	<a href="#">P25705</a>
<b>Cytogenetics:</b>	18q21.1
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
<b>Gene Summary:</b>	<p>This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase catalyzes ATP synthesis, using an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F<sub>1</sub>, and the membrane-spanning component, F<sub>o</sub>, comprising the proton channel. The catalytic portion of mitochondrial ATP synthase consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled with a stoichiometry of 3 alpha, 3 beta, and a single representative of the other 3. The proton channel consists of three main subunits (a, b, c). This gene encodes the alpha subunit of the catalytic core. Alternatively spliced transcript variants encoding the different isoforms have been identified. Pseudogenes of this gene are located on chromosomes 9, 2, and 16. [provided by RefSeq, Mar 2012]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (a). Variants 1 and 2 both encode the same isoform (a).</p>