

Product datasheet for **SC320870**

COMMD9 (NM_014186) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: COMMD9 (NM_014186) Human Untagged Clone
Tag: Tag Free
Symbol: COMMD9
Synonyms: C11orf55; HSPC166; LINC00610
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_014186.1
 GGGGGACTTCGGCAAGATGGCTGCCCTGACAGCGGAGCATTTTGCAGCACTCCAGAGCCT
 GCTCAAGGCCTCCTCGAAAGATGTTGTCAGACAGCTGTGTCAGAAAGCTTTTCCAGTTC
 AGCCCTTGGCTTGAAAACTCTTGGATGTTACATGTTCCAGCTTGCTGTGACCCAGGA
 GGAGGCAGAGGAAGTCTCCAGGCTCTGCACCGCCTCACTAGGCTGGTGGCATTCCGTGA
 CCTGTCTCTGCCGAGGCAATTCTGGCTCTTTCCAGAAAAATTTCCACAAAACCTCAA
 AAACCTGTGACAAAGATCATCTAGAACATGTGTCTACTTGGAGAACCAGCCAGGC
 AAATCAGATCTCTGCCACGCCTGGTCGATCTGGACTGGAGAGTGGATATCAAACCTC
 CTCAGACAGCATCAGCCGCATGGCCGTCACCTGCCTGCTCCAGATGAAGATCCAAGA
 AGATCCCAGCCTATGCGGAGACAAACCCTCCATCTCAGCTGTACCGTGGAGCTGAGCAA
 AGAAACTGGACACCATGTTAGATGGCCTGGGCCGATCCGAGACCAACTCTCTGCCGT
 GGCCAGTAAATGATCCAGCCAGCTGCCAGGGCCACTGCCATGACCCAGCTGCTCATGAGT
 GATAAATGTCCTCCCATATGCAGGCTGCCCTTGCAGCTGCAGCTGACAACAGGCAGGATG
 GTGGGGACAGCAGGGGGCTACTGCCATCCAGAAGTTACAGTTGGATTGGGAAGAAGCAGC
 CAGATCCCCCGCTGTTCTCACTCATCTTTCTTTCTTTGAAGCTGGAGAGCAGAAGCC
 CCCATCTTTGAAAAGCTCCTGAGTGCACTTAATTACCACCATGGCAGGGTGAGGGAACA
 TTTGCATCGTCAGCTGCCTCTGCATAGCTGTTTGAGAAATTCAGGCCAAATCATGCAGC
 CTATCCAATAAGTAAGTTTATTTCCAACATTAGCTCTAATTAGTTCATTTCCAATCCCAG
 AACACATGGAGGGAATCGGACAGGTGATGCCAGCAGTTCCTGCTCCTCTGTCAGGGAAGC
 CAGGCAGAGCCACAGAGCATGGTCCATCCAGAGTGTTCCTGAGCCCTCCACCATAC
 TGAACCCCTCTTCAGTGTAGGAAGTCTGAAATGGGTGCTAATCCCTTCTTCATGAAAC
 CAGGGCCCTTCTTCTCATCTAATGCAGCCACTCCTAGGTGAAGAAGTGGGAATAATTGG
 AAATAACAACAGTTCTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire
ACCN: NM_014186



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_014186.1 , NP_054905.1
RefSeq Size:	1318 bp
RefSeq ORF:	597 bp
Locus ID:	29099
UniProt ID:	Q9P000
Cytogenetics:	11p13
Gene Summary:	<p>May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). May down-regulate activation of NF-kappa-B (PubMed:15799966). Modulates Na(+) transport in epithelial cells by regulation of apical cell surface expression of amiloride-sensitive sodium channel (ENaC) subunits (PubMed:23637203).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>