

Product datasheet for **SC320974**

Claudin 8 (CLDN8) (NM_199328) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Claudin 8 (CLDN8) (NM_199328) Human Untagged Clone
Tag:	Tag Free
Symbol:	Claudin 8
Synonyms:	HEL-S-79
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_199328.1
 GGGACAGAGAACCCTGCTTCAAAGCAGAAGTAGCAGTTCGGGAGTCCAGCTGGCTAAAAC
 TCATCCCAGAGGATAATGGCAACCCATGCCTTAGAAATCGCTGGGCTGTTTCTTGGTGGT
 GTTGAATGGTGGGCACAGTGGCTGTCACTGTGCCTCAGTGGAGAGTGTCCGGCCTTC
 ATTGAAAAACAATCGTGGTTTTTAAAACTTCTGGGAAGGACTGTGGATGAATTGCGTG
 AGGCAGGCTAACATCAGGATGCAGTGCAAAATCTATGATTCCTGCTGGCTCTTTCTCCG
 GACCTACAGGCAGCCAGAGGACTGATGTGTGCTGCTCCGTGATGCCTTCTTGGCTTTC
 ATGATGGCCATCCTTGGCATGAAATGCACCAGGTGCACGGGGACAATGAGAAGGTGAAG
 GCTCACATTCTGCTGACGGCTGGAATCATCTTCATCATCACGGGCATGGTGGTGCTCATC
 CCTGTGAGCTGGGTTGCCAATGCCATCATCAGAGATTTCTATAACTCAATAGTGAATGTT
 GCCCAAAAACGTGAGCTTGGAGAAGCTCTACTTAGGATGGACCACGGCACTGGTGCTG
 ATTGTTGGAGGAGCTCTGTTCTGCTGCGTTTTTGTGCAACGAAAAGAGCAGTAGCTAC
 AGATACTCGATACCTTCCATCGCACACCCAAAAAGTTATCACACCGGAAAGAAGTCA
 CCGAGCGTCTACTCCAGAAGTCAGTATGTGTAGTTGTGTATGTTTTTTAACTTTACTAT
 AAAGCCATGCAAATGACAAAAATCTATATTACTTTCTCAAATGGACCCAAAGAACTT
 TGATTTACTGTTCTTAACTGCCTAATCTTAATTACAGGAACTGTGCATCAGCTATTTATG
 ATTCTATAAGCTATTTACAGCAGAATGAGATATTAACCCAATGCTTTGATTGTTCTAGAA
 AGTATAGTAATTTGTTTTCTAAGGTGGTCAAGCATCTACTCTTTTTATCATTACTTTCA
 AAATGACATTGCTAAAGACTGCATTATTTTACTACTGTAATTTCTCCACGACATAGCATT
 ATGTACATAGATGAGTGTAAACATTTATATCTCACATAGAGACATGCTTATATGGTTTTAT
 TAAAAATGAAATGCCAGTCCATTACACTGAATAAATAGAAGTCAACTATTGCTTTTCAGG
 GAAATCATGGATAGGGTTGAAGAAGGTTACTATTAATTGTTTAAAAACAGCTTAGGGATT
 AATGTCCTCCATTTATAATGAAGATTAATGAAGGCTTAATCAGCATTGTAAAGGAAA
 TTGAATGGCTTTCTGATATGCTGTTTTTTAGCCTAGGAGTTAGAAATCCTAACTTTTA
 TCCTCTTCTCCCAGAGGCTTTTTTTTTCTGTGTATTAATTAACATTTTTAAAAAGCAG
 ATATTTGTCAAGGGGCTTTCATTCAAAGTCTTTTCCAGGGCTATACTCAGAAGAAAAG
 ATAAAAGTGTGATCTAAGAAAAAGTGTGGTTTTAGGAAAGTAAAAATTTTTGTTTTT
 GTATTTGAAGAAGATGATGCATTTTGACAAGAAATCATATATGTATGGATATATTTTAA
 TAAGTATTTGAGTACAGACTTTGAGGTTTCATCAATATAAATAAAAGAGCAGAAAAATAA
 AAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_199328

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_199328.1 , NP_955360.1
RefSeq Size:	1989 bp
RefSeq ORF:	678 bp
Locus ID:	9073
UniProt ID:	P56748
Cytogenetics:	21q22.11
Protein Families:	Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction
Gene Summary:	<p>This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. This protein plays important roles in the paracellular cation barrier of the distal renal tubule, and in the paracellular barrier to prevent sodium back-leakage in distal colon. Differential expression of this gene has been observed in colorectal carcinoma and renal cell tumors, and along with claudin-7, is an immunohistochemical marker for the differential diagnosis of chromophobe renal cell carcinoma and renal oncocytoma.[provided by RefSeq, May 2010]</p>