

Product datasheet for SC321066

CGR19 (CGRRF1) (NM_006568) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CGR19 (CGRRF1) (NM_006568) Human Untagged Clone
Tag:	Tag Free
Symbol:	CGR19
Synonyms:	CGR19; RNF197
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_006568.2 <pre> GGGGGGCTGGGCTGGGCTCCGCGGCTGGAGCCGGGCTCTACCCAGAGCAAGACCCTGATG GCTGCGGTGTTTCTGGTAACGCTTTATGAATACTCGCCGCTTTTCTACATCGCGGTGGTC TTTACCTGCTTCATCGTGACCACCGGCTGGTATTGGGATGGTTTGGTTGGGATGTTCCA GTAATTCTGAGAAATTCAGAAGAGACCCAGTTCAGCACAAAGATTTTCAAAAAGCAAATG AGACAAGTCAAGAATCCTTTTGGCTTAGAGATCACTAATCCATCTTCAGCTTCAATTACA ACTGGCATAACCTTGACAACAGATTGCCTTGAAGATAGCCTCCTTACATGCTACTGGGGG TGCAATGTTCAAAAATTATGAAGCTCTGCAGAAGCATGTTTATTGCTTCAGAATAAGC ACTCCCCAAGCATTAGAAGATGCTCTGTATAGTGAATATCTCTATCAGGAACAGTATTTT ATTAAGGATAGCAAAGAAGAAATATATTGCCAGTTACCAAGAGATACTAAAATTGAA GACTTTGGTACAGTACCCAGATCTCGCTATCCATTGGTAGCGCTATTGACCTTAGCTGAT GAGGATGACCGGAAATTTATGATATTATTTCCATGGTGTGAGTATTATTCCTGAT AGGACTTATAAACTATCCTGCAGAATATTGTATCAATATTTACTCTTGGCTCAAGGTCAA TTTCATGATCTTAAGCAACTTTTCATGTCTGCAAATAATAATTTCACTCCCTCCAACAAT TCCTCTTCAGAAGAAAAAACACAGACAGAAGTTTGTGGAAAAGGTGGGACTCTCTGAA AGTGAAGTTGAGCCATCGGAAGAGAACAGCAAGGACTGTGTTGTTTGCAGAAATGGGACT GTGAAGTGGTACTCTTACCATGCAGACACACATGCCTGTGTGATGGCTGTGTGAAGTAT TTTCAGCAGTGCCCAATGTGCAGGCAGTTTGTTCAGGAATCTTTGCACTTTGCAGTCAA AAAGAGCAAGATAAAGACAAACCGAAGACTCTTTGAAGACATCGTAACACTGAAAAGTAC ACTTTCTACTAAAGATGCAGAAATTGATGATCTTGGAAATTCATCATAACATGGAATCTAC AGTACTGACCATCAATGAAAATTATTTTAACTTCATATTTGTATGGTACTTGGATGAT AAAAATTAATTATTCCTTTCTGCTTAGTGAATGAATACTGGAATCCATCTGTGTTGATAC ATAAAATTCATTCAACTCCAAAAA</pre>
Restriction Sites:	Please inquire
ACCN:	NM_006568


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006568.2</u> , <u>NP_006559.1</u>
RefSeq Size:	1323 bp
RefSeq ORF:	999 bp
Locus ID:	10668
UniProt ID:	<u>Q99675</u>
Cytogenetics:	14q22.2
Domains:	RING
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	Able to inhibit growth in several cell lines.[UniProtKB/Swiss-Prot Function]