

Product datasheet for **SC321080**

MRPL22 (NM_014180) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPL22 (NM_014180) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRPL22
Synonyms:	HSPC158; L22mt; MRP-L22; MRP-L25; RPML25
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_014180.2
 GGGGAAAGATGGCGGCGGCAGTACTGGGACAGTTGGGTGCGTTATGGATACATAACCTGA
 GGAGCCGGGGAAGCTGGCCTTGGGTGTTTTACCTCAATCATATATCCACACAAGTGCTT
 CTCTTGACATTTCTCGAAAATGGGAGAAGAAGAAATAAAATTGTTTATCCTCCACAAGTGC
 CTGGAGAACCCTCGGAGACCAGCAGAAATCTACCACTGTGCAAGACAAATAAAATATAGCA
 AAGACAAGATGTGGTATTTGGCAAAATTGATACGAGGAATGTCTATTGACCAGGCTTTGG
 CTCAGTTGGAATTCAATGACAAAAAGGGGCCAAAATAATTAAAGAGGTTCTCTTAGAAG
 CACAAGATATGGCAGTGAGAGACCATAACGTGGAATTCAGGTCCAATTTATATAGCTG
 AGTCCACCTCAGGACGAGGCCAGTGCCTGAAACGCATCCGCTACCATGGCAGAGGTGCT
 TTGGGATCATGGAGAAGGTTTATTGCCATTATTTGTGAAGTTGGTGAAGGGCCCCAC
 CTCCACCTGAGCCACAAAGACGGCAGTTGCCCATGCCAAAGAGTATATTCAGCAGCTTC
 GCAGCCGGACCATCGTTCACACTCTATGATGAGGAGATTCACTCCACAGTGTATATAT
 TTTGCCATTTATTTCTAAAAATAAACAAAAATTGAAGGCTAAAAAAAAAAAAAAAAAAAA
 AAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_014180
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_014180.2, NP_054899.2</u>
RefSeq Size:	735 bp
RefSeq ORF:	621 bp
Locus ID:	29093
UniProt ID:	<u>Q9NWU5</u>
Cytogenetics:	5q33.2
Gene Summary:	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the L22 ribosomal protein family. A pseudogene corresponding to this gene is found on chromosome 4q. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a). Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.