

Product datasheet for **SC321165**

SEC61G (NM_014302) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SEC61G (NM_014302) Human Untagged Clone
Tag:	Tag Free
Symbol:	SEC61G
Synonyms:	SSS1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<p>>OriGene sequence for NM_014302.3</p> <pre>GGCACGAGGCTGGCATTGTTAGGTGTCGGTTGGGTAGGCAGTCATGGATCAGGTAATGCAG TTTGTGAGCCAAGTCGGCAGTTTGTAAAGGACTCCATTCGGCTGGTTAAAAGATGCACT AAACCTGATAGAAAAGAATTCCAGAAGATTGCCATGGCAACAGCAATAGGATTTGCTATA ATGGGATTCATTGGCTTCTTTGTGAAATTGATCCATATTCTATTAATAACATCATTGTT GGTGGCTGAATACATTTTGAAGAGAGTTTTTCATCTTAGAGATTGGTGAACAAGTGTA GGGTGTGAGAACTCACAGAATACAAATTTGCCTGTATGTTTTGTGGTTTTTTTTTTTCC TTTCAAGATGTTTTCTATTTCTAAATTAAGTAATTTCAAAGTAAAAAAAAAAAAAAAAAAAA AAAAAAAAAAAAAAAAAAAAAAAAAAAAA</pre>
Restriction Sites:	Please inquire
ACCN:	NM_014302
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014302.3](#), [NP_055117.1](#)

RefSeq Size: 484 bp

RefSeq ORF: 207 bp

Locus ID: 23480

UniProt ID: [P60059](#)

Cytogenetics: 7p11.2

Domains: SecE

Protein Families: Transmembrane

Protein Pathways: Vibrio cholerae infection

Gene Summary: The Sec61 complex is the central component of the protein translocation apparatus of the endoplasmic reticulum (ER) membrane. Oligomers of the Sec61 complex form a transmembrane channel where proteins are translocated across and integrated into the ER membrane. This complex consists of three membrane proteins- alpha, beta, and gamma. This gene encodes the gamma-subunit protein. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein.