

Product datasheet for **SC321283**

ARPM1 (ACTRT3) (NM_032487) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARPM1 (ACTRT3) (NM_032487) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARPM1
Synonyms:	ARP-T3; ARPM1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_032487.2
 GTTGAATCTTTTATGCGAAAGTCGGGTTTCTGATACTGGGATTCCGGGATTCCAGGT
 GTTGGGGTGGCCCAATTCCTGCGAGAAGCAATAGCGGGCGGTAAACATGAGGAGCACGGTG
 CGTCCAGCGAGTCCTTCCGCCTGGGGCCCTGCCGACCCCTGCCTGCGCCCCAGGACTC
 TGGCCTCACCCGGCCGTGCCGGGGCCTCTGTGACGCGGCGTTCAGGCACTCGGCCCCGG
 CCGAGCCCGTAGCTAGAGCGGCTCAGAGACAGGAGCGGCGGCAGCAGCGCGGCATGAA
 CCACTGCCAGCTACCGGTGGTGATCGACAACGGCTCGGAATGATCAAGCGGGCGTGCC
 TGGGTGCCGGGAGCCCCAGTTTATCTACCGAACATTATCGGCCCGCCCAAGGGCCAGAG
 CCGCGCGGCCAGGGCGGCTAGAACTCTGCGTGGGCGACCAAGCTCAGGACTGGAGGAG
 CTCGCTGTTTCATCAGTTACCCAGTGAGCGTGGTCTCATTACTTCATGCGAGGACATGGA
 GATCATGTGGAAGCATATCTATGACTATAACCTAAAGCTGAAGCCGTGTGATGGCCAGT
 CTTGATTACTGAGCCAGCGCTGAACCCACTGGCCAACCGCAACAGATCACGGAATGTT
 TTTTGAGCATCTGGGTGTTCTGCCTTCTATATGTCCATCCAGGCTGTGCTGGCTCTCTT
 TGCTGCTGGCTTCACTACTGGCCTTGCTGAATTCAGGTGCTGGGGTTACCCAGAGTGT
 GCCCATCTTTGAGGGTACTGTCTGCCTCATGGTGTGACGCAACTGGATCTGGCAGGCCT
 TGACCTCACCAACTACCTCATGGTGCTAATGAAGAACCATGGTATCATGTTGCTCAGTGC
 TTCAGACAGAAAGATTGTTGAAGACATCAAGGAGAGCTTTTGTATGTGGCAATGAACTA
 CGAAGAGGAAATGGCCAAGAAACCCGATTGTCTAGAGAAAGTTTACCAACTACCTGATGG
 GAAGGTATCCAGCTCCATGACCAGCTCTTTCTTGTCCAGAGGCCCTCTCTCTCCGTG
 TCATATGAACCTTGAGGCCCCCTGGCATTGATAAGATATGCTTCAGCAGCATAATGAAATG
 TGATACAGGCCTGAGGAATTCCTTTTCCAATATTATCCTTGCCGGGGGATCAACCTC
 TTTCCCTGGTTTAGACAAGCGGTAGTTAAGGATATAGCAAAGGTGGCTCCTGCCAACAC
 CGCTGTGCAAGTTATAGCTCCTCCAGAAAGGAAATATCAGTGTGGATGGGAGTTCTAT
 TCTTGATCCTTGTCTGCCTTCCAGGACATGTGGATCACTGCTGCAGAATTTAAAGAAAT
 TGGACCCAAACATAGTACACCAAGATGCTTCTGAAATACAGATAAAATGTTGGAAGAAA
 ATGTTTTGAGTATATGTGACAGAAAACCTTGGATATTATGTTTCTGGGAGAAGAGAAA
 ATACTTCACCTATTGGGATGCCAATATTCTGTTGATTCTATAATGGGTTTGGGGGAT
 AATAATGGTGAAGCTCAAGAACAGATGTCTATTGAGTAGAACCAAGTTAAATAAAAAA
 AAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_032487

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032487.2](#), [NP_115876.2](#)

RefSeq Size: 2002 bp

RefSeq ORF: 1119 bp

Locus ID: 84517

UniProt ID: [Q9BYD9](#)

Cytogenetics: 3q26.2