

Product datasheet for **SC321322**

CXX1 (FAM127A) (NM_001078171) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	CXX1
Synonyms:	CXX1; FAM127A; Mar8; MAR8C; Mart8; MART8C; SIRH5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_001078171.1
GAGGAGCAGCGATGGACGGTCGGGTGCAGCTGATAAAGGCCCTCCTGGCCTTGCCGATCC
GGCCTGCGACGCGTCGCTGGAGGAACCCGATTCCCTTTCCCGAGACGTTTGACGCGGATA
CCGACCGACTCCCGGAGTTCATCGTGCAGACGGGCTCCTACATGTTCTGGACGAGAACA
CGTTTCTCCAGCGACGCCCTGAAGGTGACGTTCTCATCACCCGCCCTCACAGGGCCCGCCC
TGCAGTGGGTGATCCCTACATCAAGAAGGAGAGCCCCCTCCTCAATGATTACCGGGGCT
TTCTGGCCGAGATGAAGCGAGTCTTTGGATGGGAGGAGACGAGGACTTCTAGGCCGGGA
GACCCTCGGGCCTGGGGCGGGTGCTCTGGGGAGGGTCCGCTGTGTTACTGGCCGCCGCC
AGGGTGCACACCGCGCCCTCCCTCCGCGCCTCCCTCCCCCTCGAGCCGCCGCGATGTCC
CCTGCGCTCCTGTTCCCTCCCGCGTAGTGCTTGCCTTTGTTCCAGGAATAGCGCTCCAGG
CTCCTGCTGCGCCCTGGGCTCACTCTGGAGCGAGCCGCCGCCCTCTCCTTCCAGCCA
GCCAGCCCCTCCCATGTACATTTGGACGCTGCTGCGCTCCAGCTGCAAGCTGGGCTCC
TGTTACACACTGGACAGACACCCACTGCCGCCGCTGCCAAGCCCTCTCCTCCCAACAG
ACTGCCAGACGACTACATCTGCCCCACAGACCTGCGCTGCCACAGCCATCGCCATCC
ATCGCATCCCACCGACAGACTGCTGCTCCTAGTGATCTGGACTCACCTCGGAGGTATCTG
GGCTGGCCACAGTCCCTGGACAGTATCCAGACAGCTGGCCGCCCCCAAGGGATCTGTC
ACCTTCAGCGAGACCTATTTCTCCCAACCCCAAGAACTCTTGTTCTTGCTAGGC
CCAGGTGTTCTGGCAGCCAAATCGAGTCTCTATTTCTCTTGTTGACAGTTAGTTTT
GCCATAACGCACTATTCTGAGTTTGCAACTGTCTCTGATGTGTCCTTTTGTTC AAC
ACAGTAACCCCTGCATTCTGCTCTGCTCTAATACTACCTGGAGAAAGTCTTTTCCTTA
TTTTCAATAAATGTCAGACATTATTGAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_001078171



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001078171.1 , NP_001071639.1
RefSeq Size:	1254 bp
RefSeq ORF:	342 bp
Locus ID:	8933
UniProt ID:	A6ZKI3
Cytogenetics:	Xq26.3