

Product datasheet for **SC321327**

ZSWIM1 (NM_080603) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZSWIM1 (NM_080603) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSWIM1
Synonyms:	C20orf162
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_080603.3
 GGCGACGGAGCCTGGCTGTGGGCCCATCTTTGGAAAAAGATCTGGGAATGATTGTCTAG
 CCTCCAGCCTCAACTTACTTGATGCTTGAGAGACTCAAAGCCCGTGGTCAGCTGCCCTG
 CAAAGAAAGTATTTTGACCTTGGCATTGGACAGCTCCCATCTCTCCATGGCCCTGACA
 ATGCTGAATGGGCTCCTGATTAAGGACTCAAGCCACCTATGCTGCTGCACCAGGTTAAC
 AAGACTGCCAGTTAGATACCTTCAACTACCAGAGCTGCTTTATGCAAAGTGCTTTTGAC
 CATTTCCCTGAGATCTTATTTATCCACCGGACCTATAACCCAAGGGTAAGGTCTTATAT
 ACCTTCCTGGTGGATGGACCTCGGGTGCAGCTGGAGGGTCATCTTGCCCGAGCAGTCTAC
 TTGCCATCCCTGCCAAGGAGGACACTGAAGGCCTGGCCAGATGTTCCAAGTATTCAAG
 AAGTTTAAATCCAGCATGGGAGAGAGTCTGTACCATCCTGGTGGATCCTCATTTCTTCCA
 CTGCCTATCCTAGCTATGGAGTTCCCCACAGCTGAGGTCTTCTCTCAGCCTTCCACATT
 TGTAAGTTCCTCCAGGCCAAGTTCTATCAGCTGTCCCTGAACGGCCCGTGGAAAGGCTG
 CTCCTGACCTCCCTGCAGAGCACAATGTGCTCAGCCACAGCAGGCAACCTGAGAAAGTTG
 TATACACTCCTGAGCAACTGCATCCCTCCAGCCAAGCTGCCGAGCTTCACTCACACTGG
 CTGCTCAACGACCGCATCTGGCTGGCTCACCCTGGAGAAGCCGAGCTGAGAGCAGCCAC
 TACTTCCAGAGCCTCGAGGTCAACACCCACATCCTCAGCCAGTTCTTTGGTACCACCCCA
 TCTGAGAAACAAGGTATGGCTTCTCTGTTCCGTTACATGCAGCAGAACTCTGCAGACAAG
 GCAAACCTCAACCAGGGCCTGTGTGCCCAAGCAATCATGCTCCCTCAGACACCATCCCC
 GAAAGCCCCAACTGGAGCAGCTGGTAGAATCCCACATCCAGCACTCCCTCAATGCCATC
 TGCACAGGGCCAGCAGCCCAACTGTGCCTGGGCGAGCTTGCTGTGGTCCAGAAATCCACA
 CACCTCATTGGCTCTGGCTCAGAAAAGATGAACATACAGATCCTGGAAGATACCCATAAG
 GTGCAGCCCCAGCCCCCTGCCAGCTGCAGCTGCTACTTTAACCAGGCCTTCCACCTGCCC
 TGCCGCCACATCCTAGCCATGCTCAGTGCCCGCCGCGAGGTGCTCCAGCCGACATGCTG
 CCGGCTCAGTGGACGGCAGGCTGTGCTACCAGTCTAGACAGCATCCTGGGCAGCAAGTGG
 AGTGAGACCTGGATAAGCACCTGGCAGTGACTCACCTACCGAGGAGGTGGGTGAGCTG
 TTGCAGCACTGCACCAAGGAGGAGTTTGAGCGGAGGTATAGCACCTGCGGGAAGTGGCC
 GACAGCTGGATTGGGCCTTATGAGCAGGTCCAACCTGATTATTCTCGATGCCAGAGAT
 GCTCATGCACCTGTGCACACTCACATCCACCCATACACACACACACACACACACACAC
 ACACACACACTCCCTTACACTGTTGTACTTCCGTGGGCCCTCCTTCCAGAACAAGGAC
 AACAAGGACAAGGTTGAAGGGTCTTCTCATCTACCATGGCCTGTACCTAGCATGTGTCT
 AGCTCAATGAGACAGGAGTCAGCAAATCTTAATCTGTTTACTTTACTCAGGTGGCCACAT
 ACAGTCTCTGTTGTATATTCTTGGTTTTGTTTTAATTTTTTTTTCTTTTTTTTTTTTTT
 TTTTTTTTTGAGAGGGAGTCTGTCTCCAGCCTGGGAGGGTGAGACTCCATCTAAAAAAA
 TAAATAAATGGCAACCCCTGGTCTAAGATAAGAGATAAAACATCAGGTGGTGGGTTGA
 GGTTTGGGGCTTGGTAGCAGTTGCCCCAGTCATGAGATGACTCACTTAACCCGCTCCTT
 TAAGTGAGCTGGGCTGGGAGGCTTCCTACAGGGGAAGAGGCCCTCTGGGAGCTGACTC
 AGCCAGGCTCCCTGAACTTTTTCTTGTCCCATCCTGGGGTCAATAAACTGAATGTTG
 CATATTCTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire
ACCN: NM_080603

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_080603.3](#), [NP_542170.3](#)

RefSeq Size: 2206 bp

RefSeq ORF: 1458 bp

Locus ID: 90204

UniProt ID: [Q9BR11](#)

Cytogenetics: 20q13.12

Domains: SWIM