

Product datasheet for **SC321426**

CYB5R3 (NM_000398) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CYB5R3 (NM_000398) Human Untagged Clone
Tag:	Tag Free
Symbol:	CYB5R3
Synonyms:	B5R; DIA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_000398.4
 GACAGAGCGAGCGCGCGGGCCACCATGGGGGCCAGCTCAGCACGTTGGGCCATAT
 GGTGCTCTTCCCAGTCTGGTTCCGTACAGTCTGCTCATGAAGCTGTTCCAGCGCTCCAC
 GCCAGCCATCACCTCGAGAGCCCGGACATCAAGTACCCGCTGCGGCTCATCGACCGGA
 GATCATCAGCCATGACACCCGGCGCTTCCGCTTGGCCGTCACCCAGCACATCCT
 GGGCCCTGTGCGCCAGCACATCTACCTCTCGGCTCGAATTGATGGAACCTGGTCTGT
 CCGGCCCTATACACCCATCTCCAGCGATGATGACAAGGGCTTCGTGGACCTGGTCATCAA
 GGTTTACTTCAAGGACACCCATCCCAAGTTTCCGCTGGAGGGAAGATGTCTCAGTACCT
 GGAGAGCATGCAGATTGGAGACACCATTGAGTTCGGGGCCCGAGTGGGCTGCTGGTCTA
 CCAGGGCAAAGGGAAGTTCGCCATCCGACCTGACAAAAAGTCCAACCCTATCATCAGGAC
 AGTGAAGTCTGTGGGCATGATCGCGGGAGGACAGGCATCACCCGATGCTGCAGGTGAT
 CCGCGCCATCATGAAGGACCTGATGACCACACTGTGTGCCACCTGCTCTTGGCAAACA
 GACCGAGAAGGACATCCTGCTGCGACCTGAGCTGGAGAACTCAGGAACAAACATTCTGC
 ACGCTTCAAGCTCTGGTACACGCTGGACAGAGCCCTGAAGCCTGGGACTACGGCCAGGG
 CTTTCGTGAATGAGGAGATGATCCGGGACCACCTCCACCCCCAGAGGAGGAGCCGCTGGT
 GCTGATGTGTGGCCCCACCCATGATCCAGTACGCTGCCTTCCAACCTGGACCACGT
 GGGCCACCCACGGAGCGCTGCTTCGTCTTCTGAGGGCCGGGACGGTACACGGCCACC
 CGCCCCGCGCACCCACGCCCTGTTACGCTCACCCAGTACCTCCCCACATCGCACACT
 GGGGCCCGGGTTCAGCCTGGCCTGCCCTGCCCTGGTGAATCACCTGGCTGAGCAGTTC
 CCCTGGAGCCCTTCGGGAGCAGGGCTGTGTCCAGATGGGCCACGGCTGAGCCTTCAGA
 GTACGTCTGCCTGGCACTTACTGGTCCCTACCAGAGACGCCAGCCCCATCCCTGTCTCT
 CATGACCCCTCGTCCACCCCCACACACTATAAGGCTGAGGGCTGCCAGCAGCCCGT
 CTGCCCCACCATTCGGGCGGTGGACCATAGTCGGGATGTCAGCAGACACACATGGGCAGC
 CCAAAGCTGCAGGTGCCAGGGCCACCCAGCCTCGCCTGTCACCCCACTCCCGCCTCA
 GGGCCAGGCCAGGCCTACCACCTGACCTGCATGAGACATTGACACCAGAAAGCCCTC
 TTGGGGCACTGCTCCCTACCCAGGGCCCTGGCCAGCCGGGAGCTTGGCTCTCCTCTGG
 CTAGAGTGGGAAGAGGGGGCTGGCCATGGGGCCCTCCCAGAACCTCAGCATTTCTTCCA
 GCCCATCCAAACTGAGGCAGCCTTGGGGAACCCGAGCTGGGGGTTGGCAGCCCACT
 GCACCGCCTCAGGGTTTTGGGGTCTGGGCTGGGGCCACCATCCCTGATGGCAGAACTCC
 CACAACCATGTATTTATCTCTGTCCTAAACCGTCCCCTCCTCCCTCACCCCCAGC
 ACAGGGGATTCTGAGCAGTGCCTCTTGTCTGAGGGACATATCAGTGACCTCGACGTTGC
 CTTTAGACTACAGTTGTGTTAGCCTCTTGCCTATTGGCTTTTTTTCAGAGTCATTTATGAGC
 AGAAAAAAAAAAGTAAACTTTGCTAATATTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
 AA

Restriction Sites: Please inquire

ACCN: NM_000398

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_000398.4](#), [NP_000389.1](#)

RefSeq Size: 1978 bp

RefSeq ORF: 906 bp

Locus ID: 1727

UniProt ID: [P00387](#)

Cytogenetics: 22q13.2

Domains: NAD_binding_1, FAD_binding_6

Protein Families: Druggable Genome

Protein Pathways: Amino sugar and nucleotide sugar metabolism

Gene Summary: This gene encodes cytochrome b5 reductase, which includes a membrane-bound form in somatic cells (anchored in the endoplasmic reticulum, mitochondrial and other membranes) and a soluble form in erythrocytes. The membrane-bound form exists mainly on the cytoplasmic side of the endoplasmic reticulum and functions in desaturation and elongation of fatty acids, in cholesterol biosynthesis, and in drug metabolism. The erythrocyte form is located in a soluble fraction of circulating erythrocytes and is involved in methemoglobin reduction. The membrane-bound form has both membrane-binding and catalytic domains, while the soluble form has only the catalytic domain. Alternate splicing results in multiple transcript variants. Mutations in this gene cause methemoglobinemias. [provided by RefSeq, Jan 2010]

Transcript Variant: This variant (1) encodes the membrane-bound form (isoform 1) of the enzyme. This isoform (1) is referred to in the literature as isoform m. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.