

Product datasheet for SC321435

FAM118A (NM_017911) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM118A (NM_017911) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM118A
Synonyms:	C22orf8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_017911.1

```

AGCGCGAGCTGGGACCAGCGTCTCGGAGGCGCCGAGAATTCACAGATGGATTCACTGGA
AAAGACAACAAATAGAAGTGAACAAAAATCCAGAAAGTTTTTAAAAAGCCTCATCCGCAA
ACAGCCCCAGGAAGTCTCTGGTTATCGGGACTGGCGTCAGCGCAGCAGTGGCCCCCGG
AATCCCTGCCCTTTGCTCGTGGAGAAGCTGCATCGAGGCCGTCATCGAGGCTGCAGAGCA
GCTGGAGGTGCTGCACCCCGGAGACGTCGCCGAGTTCCGGAGGAAAGTGACAAAGGACCG
GGACCTGTTGGTTGTCGCCCATGATCTGATCCGGAAGATGTCACCTCGCACAGGCGATGC
CAAGCCCAGCTTCTCCAGGACTGCCTGATGGAGGTGTTTGACGACCTGGAGCAGCACAT
CCGGAGTCTCTGGTGTGACGTCGATCCTCAGCCTGATGGACAGAGGCCCATGGTCTCT
GACCACCAACTATGACAACCTGCTGGAGGCCTTTGGCCGGCGGCAGAACAAGCCCATGGA
GTCCCTGGACTTGAAGGACAAGACCAAGTCTTGAATGGGCAAGAGGGCACATGAAGTA
CGGCGTCTCCACATTCACGGCCTCTACAGGGACCCCTGCGGGGTGGTGTGGACCCATC
GGGGTATAAAGACGCTCACTCAAGACGCGAGAAGTCATGGAAGTCCTCCAGAACTTATACCG
CACCAAGTCCTTTCTGTTTGTGGGCTGTGGGGAGACCCCTTCATGATCAGATATCCAGGC
CCTCTTTCTTTACTCCGTGCCGAATAAGGTGGATTGGAGCACTACATGCTTGTGCTGAA
GGAGAATGAAGACCATTCTTTAAGCATCAGGCAGATATGCTTCTGCACGGAATCAAAGT
TGTATCTACGGGACTGTTTTGACCACTTTCCAGGATATGTGCAAGACCTTGCCACTCA
GATCTGCAAAACAGCAAAGCCAGATGCTGATCGCGTGGACAGCACCATATTGGGTAA
TGATGCCAGGACTGTGCAAAGAGGAAGTTAGAAGAGAATGGAATTGAAGTTTCAAAAAA
ACGCACACAATCAGATACTGATGATGCTGGAGGGTCTTGAATCTTTACAGTAAACCTG
CAACTTGAAGAACTAGCCTTCTGTAACACAGTGCCCAACGAAGAGGAATGTATGGAGAA
CTCCACGTGGATCTCTGATTGCGAAACCGTCACATACACCAAGAGAGCCACATGGGCATG
TGGCCCTGAAGGCTGGGTGAGAGGGCTCCCTGTGTGTTGAAGTATGCAGGAGGGTGACG
CGGACACATTTCAAGTGGACTTTGCAAGGACTGATGGATAGCTACCTCAGGGACCAGAAT
CCGTGGGAAGGGATGGACCTGGTGTCCCGTCCCATCTGACAGGCTCTCTTTTGTCAAG
GTGGTATTTTTCGTAATAAAAGGGGAAGAGTAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAA
  
```


[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_017911
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017911.1</u> , <u>NP_060381.1</u>
RefSeq Size:	1707 bp
RefSeq ORF:	1074 bp
Locus ID:	55007
UniProt ID:	<u>Q9NWS6</u>
Cytogenetics:	22q13.31