

## Product datasheet for **SC321529**

### YIPF5 (NM\_030799) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	YIPF5 (NM_030799) Human Untagged Clone
Tag:	Tag Free
Symbol:	YIPF5
Synonyms:	FinGER5; SB140; SMAP-5; SMAP5; YIP1A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

<b>Fully Sequenced ORF:</b>	<p>&gt;OriGene sequence for NM_030799.6</p> <pre> GTAGCAACGGGGCTGGTTTCAGGGTCTGAAACAGAGTTTGGGGTTGTTTGGGATTAGTGA AGCTACTGCCTTTGCCGCCAGCGCAGCCTCAGAGTTTGATTATTTGCAATGTCAGGCTTT GAAAACCTAAACACGGATTTCTACCAGACAAGTTACAGCATCGATGATCAGTCACAGCAG TCCTATGATTATGGAGGAAGTGGAGGACCCTATAGCAAACAGTATGCTGGCTATGACTAT TCGCAGCAAGGCAGATTTGTCCCTCCAGACATGATGCAGCCACAACAGCCATACACCGGG CAGATTTACCAGCCAACCTCAGGCATATACTCCAGCTTCACCTCAGCCTTTCTATGGAAAC AACTTTGAGGATGAGCCACCTTTATTAGAAGAGTTAGGTATCAATTTTGACCACATCTGG CAAAAAACACTAACAGTATTACATCCGTTAAAAGTAGCAGATGGCAGCATCATGAATGAA ACTGATTTGGCAGGTCCAATGGTTTTTGCCTTGTCTTTGGAGCCACATTGCTACTGGCT GGCAAAATCCAGTTTGGCTATGTATACGGGATCAGTGCAATTGGATGTCTAGGAATGTTT TGTTTTATAAACTTAATGAGTATGACAGGTGTTTCATTTGGTTGTGTGGCAAGTGCCTT GGATATTGCTTCTGCCATGATCCTACTTTCCAGCTTTCAGTGATATTTTCTTTGCAA GGAATGGTAGGAATCATTCTCACTGCTGGGATTATTGGATGGTGTAGTTTTTCTGCTTCC AAAATATTTATTTCTGCATTAGCCATGGAAGGACAGCAACTTTTAGTAGCATATCCTTGC GCTTTGTTATATGGAGTCTTTGCCCTGATTTCCGTCTTTTGAAAATTTATCTGGGATGTG GACATCAGTGGGCCAGATGTACAAAAAGGACCTTGAACCTTAAATTGGACCAGCAAACCT GCTGCAGCGCAACTCTCATGCAGATTTACATTTGACTGTTGGAGCAATGAAAGTAAACGT GTATCTCTTGTTCAATTTTATAGAACTTTTGCATACTGTATTGGATTTACCTGCGGTGTG ACTAGCTTTAAATGTTTGTGTTTATACAGATAAGAAATGCTATTTCTTTCTGGTTCTCTGC AGCCATTGAAAAACCTTTTCTTGCAAATTATAATGTTTTTGATAGATTTTATCAACT GTGGGAAACCAACACAAAAGCTGATAACCTTTCTTAAAAACGACCCAGTCACAGTAAAGA AGACACAAGAAGGCCGGCGTGGTAGCTCACGCCTGTAATCCAGCACTTTGGGAGGCCG AGGCGGGCGGATCACAAGGGCAGGAGATCGAGACCATCCTGGTTAACACGGTGAACCCCC GACTCTACTAAAACACAAAAAAATTAGCTGGGCGTGGTGGCGGGCCCTGTAGTCCCA GCTACTCAGGAGGCTGAGGCAGGAGAAAAGTGTGAACCCAGGAGGCGGAGCTTGCAGTGA GCCGAGATCACACCACTGCACTCCACGCTGGGTGACAGGGTGTAGACTCTGTCTCAA AAAAAAAAAAAAAAAAA </pre>
<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_030799
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_030799.6</a></u> , <u><a href="#">NP_110426.4</a></u>
<b>RefSeq Size:</b>	3126 bp
<b>RefSeq ORF:</b>	774 bp
<b>Locus ID:</b>	81555
<b>UniProt ID:</b>	<u><a href="#">Q969M3</a></u>
<b>Cytogenetics:</b>	5q31.3
<b>Domains:</b>	Yip1
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	Plays a role in transport between endoplasmic reticulum and Golgi.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1 and 2 encode the same protein (isoform a).