

## Product datasheet for **SC321540**

### ASF1A (NM\_014034) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ASF1A (NM\_014034) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ASF1A  
**Synonyms:** CGI-98; CIA; HSPC146  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC (PS100020)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**Fully Sequenced ORF:** >OriGene sequence for NM\_014034.2  
AAACACTTGGAGTGCTCCCGTGTAAATAAAAAGAGGAAAAAGTTTCTCAAGTCGCCGCT  
GCACGACGTCTGGCCGGCGCTGGAGCGGGGTCTGCGCTCTCCCGAGCGCCGCGCGCTG  
GACTTTATTGTGCGCAACCAGCCCAAGTCCCAATTGTTGTGTTTTTTTCAAAAATATGG  
CAAAGGTTTCAGGTGAACAATGTAGTGGTCTGGATAAACCCTTCTCCTTTCTACAACCCGT  
TCCAGTTCGAGATCACCTTCGAGTGCATCGAGGACCTGTCTGAAGACTTGAATGAAAA  
TTATCTATGTGGGCTCTGCAGAAAGTGAAGAATACGATCAAGTTTTAGACTCTGTTTTAG  
TGGGTCTGTTCCCGCAGGAAGGCATATGTTTGTATTTTCAGGCTGATGCACCTAATCCAG  
GACTCATTCCAGATGCAGATGCAGTAGGCGTAAGTGTGCTAATTACTTGTACCTATC  
GAGGACAAGAATTTATTAGAGTTGGCTATTATGTAATAATGAATATACTGAGACAGAAT  
TAAGGGAAAATCCACCAGTAAAACCAGACTTTTCTAAGCTTCAAAGGAATATTTTGGCAT  
CTAATCCAGGTCACAAGATCCACATTAATTGGGAAGATAACACAGAAAACTGGAAG  
ATGCAGAGAGCAGTAATCCAAATCTACAGTCACTTCTTTCAACAGATGCATTACCTTCAG  
CATCAAAGGGATGGTCCACATCAGAAAACCTCACTAAATGTCATGTTAGAAATCCCACATGG  
ACTGCATGTGACCACCTACCATCCCTTTAGTACAAATTAAGCTATTAATAACACAGAA  
CTATTTCCCTGAAATCCGTAAGTACATAGTCAAAACACAATGTGAAGAAATTTGTTTAA  
AACATCCTGTAGAAAGTTATAAGAAAACCAGTATTTGAACAAATTTGGAATATAAATA  
CAACTATTTTTAAGTAAAAAAAAAAAAAAAAAAAA

**Restriction Sites:** Please inquire

**ACCN:** NM\_014034

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_014034.2</a> , <a href="#">NP_054753.1</a>
<b>RefSeq Size:</b>	2525 bp
<b>RefSeq ORF:</b>	615 bp
<b>Locus ID:</b>	25842
<b>UniProt ID:</b>	<a href="#">Q9Y294</a>
<b>Cytogenetics:</b>	6q22.31
<b>Domains:</b>	Anti-silence
<b>Gene Summary:</b>	This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair. [provided by RefSeq, Jul 2008]