

Product datasheet for **SC321541**

LIN37 (NM_019104) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIN37 (NM_019104) Human Untagged Clone
Tag:	Tag Free
Symbol:	LIN37
Synonyms:	F25965; lin-37; ZK418.4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_019104.1

```

CCCCGCCTTCTCCCGCCTTGACTTGTGACCCTAGGCCCTTTGGGGCGCCTCTGACCCAGC
TAGCCAGATCCCGGACCCAAACCATGTTCCCTGTGAAGGTGAAAGTGGAGAAATCAGAGC
TGGAGATGGCCAAAGCCCGGAACCAACTGGATGCTGTCTTGCAAGTGTCTGCTGGAGAAGA
GTCACATGGACAGGGAGCGTCTGGATGAGGAAGCTGGGAAAACACCCCTCAGACACCCACA
ATAAGGACTGCTCCATCGCAGCCACTGGCAAAAGGCCATCTGCCCGCTTCCCCACCCAGC
GGAGGAAGAAGAGGAGGGAGATGGATGATGGGCTGGCTGAGGGAGGGCCGAGCGATCCA
ACACATATGTGATCAAGCTGTTGACCGGAGCGTGGACTTGGCCAGTTTACGCGAGAACA
CGCCACTGTACCAATCTGCCGCGCCTGGATGCGCAACAGCCCTCTGTGCGCGAGCGTG
AATGCTCTCCAGCTCACCCCTGCCCGCCTGCCTGAGGATGAGGAGGGCTCAGAGGTA
CCAACAGCAAGAGTCGTGATGTGTACAAGCTGCCGCCACCCACACCCCGGGGCCACCCG
GAGATGCCTGCAGATCCCGCATCCCATCTCCACTGCAGCCTGAGATGCAGGGCACCCCTG
ACGATGAGCCCTCTGAGCCCGAGCCCTACCCCTCCACACTCATCTATCGAACATGCAGC
GCTGGAACGCATCCGCCAGAGGTGGAAGGAGGCCTCTCATCGGAACAGCTTTCGTTACT
CAGAAAGCATGAAGATCTACGAGAGATGTACGAACGACAGTGATGTTCCAGGTCCCCC
CACACCAGTAAACATCCCCAGCTCCACAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```

Restriction Sites: Please inquire

ACCN: NM_019104

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019104.1</u> , <u>NP_061977.1</u>
RefSeq Size:	929 bp
RefSeq ORF:	741 bp
Locus ID:	55957
UniProt ID:	<u>Q96GY3</u>
Cytogenetics:	19q13.12
Gene Summary:	This gene encodes a protein expressed in the eye. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) encodes the longer isoform (1).