

Product datasheet for **SC321662**

AP2 alpha (TFAP2A) (NM_001032280) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AP2 alpha (TFAP2A) (NM_001032280) Human Untagged Clone
Tag:	Tag Free
Symbol:	AP2 alpha
Synonyms:	AP-2; AP-2alpha; AP2TF; BOFS; TFAP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_001032280.2
CACGCACGCAGACACATTCTCCCTCTCTGCGCTCTCTCCTTCGGTCTCCCTCTCTGCCTC
TTTCTCTTTCACTTTTGCATGCCCTGCAACCTTTTAAATGTTGCCCTTCCTGTGATT
CGCCAGACGCCGCGCGGCCCGCCCGCGCTCGCCCGTCCCTCTCTCGCTCGTTTT
TGTCTCTCGCGTCCCTCTCCCACTCCGATTTGCTACACTGAGACTCCCGTCAATGGAC
TGCATTGAGACCGGCTCCGCGCGAGTTGCCTCTCCGCTTACGCTCGATTTCCAGGCA
TTCTTCCCTTATTAAGTATTCGTGTAATTAATAGTCATGAATATCTGCTATTAGGAGG
CTCCAGGAACGCTGCCAGCGCGTATTAGAAGCTCAAGCGAAGCCCGCTAAGAAAA
GAGGGGAGACACGGATTAAGGAACACGCACGCCACACACTCACACATACTTTTTCC
TTTTCTTTTTTGGGGTTTTTCATTTTTAAAGAACTTTGAATCATAACCAGTCGCGCA
GGATAGAGACCGTGGGTCGACCAGCTGAAGGCGCCGCGCAATCGGTGGTTCAAGTTTCG
GATGGATCCGAGCCGGCTCCCGCGCTCCGGGAGCGTCGGCGCCGTGACTGGAGTCTG
GGTGGCGCGGGCTCTCGCGCCTTTTGTGTGGGCGCGCTCGGGCGCGGGGACGCGG
GCGCTCCAGGCTTGTGTTTTTTCCACCTGACTCGTTACCCAGACTCTTCGCAGAT
GTTAGTTCACAGTTTTTACGCAATGGACCGTCACGACGGCACCAGCAACGGGACGGCAG
GTTGCCCAAGCTGGGCACTGTAGGTCAATCTCCCTACACGAGCGCCCGCGCTGTCCCA
CACCCCAATGCCGACTTCCAGCCCCATACTTCCCCCACCTACCAGCCTATCTACCC
CCAGTCGCAAGATCCTTACTCCCACGTCAACGACCCCTACAGCCTGAACCCCTGCACGC
CCAGCCGAGCCGAGCACCAGGCTGGCCCGCCAGAGGCAGAGCCAGGAGTCTGGGCT
CCTGCACACGCACCGGGGCTGCCTCACCAGCTGTCCGGCCTGGATCCTCGCAGGGACTA
CAGGGCGCAGGAGACCTCCTGCACGGCCACACGCGCTCAGCTCAGGACTCGGAGACCT
CTCGATCCACTCCTTACCTCACGCCATCGAGGAGTCCCGCATGTAGAAGACCCGGTAT
TAACATCCAGATCAAACGTGAATTAAGAAAGGCCCGTGTCCCTGTCCAAGTCCAACAG
CAATGCCGTCTCCGCCATCCCTATTAACAAGGACAACCTTTCGCGCGGTGGTGAACCC
CAACGAAGTCTTCTGTTCAAGTTCGGGTCGCTCTCGCTCCTCAGCTCCACCTCGAAGTA
CAAGGTACGGTGGCGAAGTGCAGCGCGGCTCTCACCACCGAGTGTCTCAACGCGTC
GCTGCTGGGCGGAGTGTCCGGAGGCGAAGTCTAAAAATGGAGGAAGATCTTTAAGAGA
AAAAGTGGACAAAATAGGATTAATCTGCCTGCAGGGAGACGTAAGCTGCCAACGTTAC
CCTGCTCACATCACTAGTAGAGGAGAAGCTGTCCACCTAGCCAGGGACTTTGGGTACGT
GTGCGAAACCGAATTTCTGCCAAAGCAGTAGCTGAATTTCTCAACCGACAACATCCGA
TCCCAATGAGCAAGTGACAAGAAAAACATGCTCCTGGCTACAAAACAGATATGCAAAGA
GTTACCGACCTGCTGGCTCAGGACCGATCTCCCTGGGAACTACGGCCCAACCCAT
CCTGGAGCCCGCATCCAGAGCTGCTTGACCCACTTCAACCTCATCTCCACGGTTCGG
CAGCCCCGCGGTGTGTCCGCGGTACGGCCCTGCAGAACTATCTCACCGAGGCCCTCAA
GGCCATGGACAAAATGTACCTCAGCAACAACCCCAATAGCCACACGGACAACAACGCCAA
AAGCAGTGACAAAGAGGAGAAGCACAGAAAAGTGAAGGCTCTCTCCCGCCCCGCCCTCCC
ACGCTCACCAGCCCCCGCGGCCACCCTCCGGCGGGTACAGCTCCGGGATCAGCAA
CCCTTCTGCTGCTACTGCTGCTGCTGCTGCCCGCCCGCCCGCCGCTGCCCTTG
GGTCCCCCGAGTCTCCGGGACTGCCCTCTCGACTGTCAGTGGGGCAGCCTCTCCGACTC
TGCACCCGCTCGACCTCCCACCCGCTCCCACACCCTGTGCCCTCATGTGGAGCCTAA
GAGAACAGAACAGCCGTGAAGCCAGCAGAGAAAAGTTCTGCCAAGTTTGTGAACCTTT
TTTTTTTAAACAAAACAACAAATCAACAACAGCAACAACAACAATAAAAACTT
TTTTCTAAAAAAAAGTAAAAATAAAAAAATATATGCGCTTATGGGACTGAGTCACC
ACCTTCCCTTACATACTTCAGTTGATGAGCCATACTTAAAAAAAAAAAAA

Restriction Sites: Please inquire
ACCN: NM_001032280

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001032280.2</u> , <u>NP_001027451.1</u>
RefSeq Size:	3808 bp
RefSeq ORF:	1296 bp
Locus ID:	7020
UniProt ID:	<u>P05549</u>
Cytogenetics:	6p24.3
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	<p>The protein encoded by this gene is a transcription factor that binds the consensus sequence 5'-GCCNNNGGC-3'. The encoded protein functions as either a homodimer or as a heterodimer with similar family members. This protein activates the transcription of some genes while inhibiting the transcription of others. Defects in this gene are a cause of branchiooculofacial syndrome (BOFS). Three transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Dec 2009]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (b) has a shorter and distinct N-terminus compared to isoform a.</p>