

## Product datasheet for **SC321728**

### ZCCHC9 (NM\_032280) Human Untagged Clone

#### Product data:

|                           |                                         |
|---------------------------|-----------------------------------------|
| Product Type:             | Expression Plasmids                     |
| Product Name:             | ZCCHC9 (NM_032280) Human Untagged Clone |
| Tag:                      | Tag Free                                |
| Symbol:                   | ZCCHC9                                  |
| Synonyms:                 | PPP1R41                                 |
| Mammalian Cell Selection: | Neomycin                                |
| Vector:                   | pCMV6-AC (PS100020)                     |
| E. coli Selection:        | Ampicillin (100 ug/mL)                  |
| Fully Sequenced ORF:      | >OriGene sequence for NM_032280.1       |

```

CCACGCGTCCGCCCGGCTCGCCAACTCGGCTGCTCTGGGGGATTCTGCGCGGTACCAC
TGATGAAATTATGACCAGGTGGGCGCGAGTTAGTACCACATATAACAAGAGAGCCTTGCC
TGCAACATCATGGGAGGACATGAAGAAGGATCCTTTGAGGGAACAAGCCAAAACCTACC
AAAGCGTAAACAACCTTGAAGCCAATAGGCTATCCCTCAAAAATGATGCACCCCAAGCAAA
ACATAAAAAGAACAAGAAAGAAAAAGAGTAAATGAAGATGTGAATGGATTCATGGA
ATACCTAAGACAGAATTCACAGATGGTTCACAATGGGCAAATTATAGCAACAGACAGTGA
GGAAGTAAGGGAAGAAATTGCAGTTGCTTTAAAGAAAGACAGTCGACGGGAAGGAAGAAG
ATTAAGGACAGACGCGGCAAGAAAAATGCAATGGTGTGTTTCCATTGTAGAAAACCTGG
TCATGGAATTGCAGATTGCCCGCTGCCCTTGAAAATCAAGACATGGGCACTGGGATATG
TTACAGGTGTGGGTCCACAGAGCACGAAATAACCAAGTGAAGGCTAAAGTAGACCCGGC
TCTTGGCGAATTTCTTTTGCAAAATGTTTTGTTGTGGAGAAATGGGGCACCTGTCTAG
ATCTTGTCTGATAATCCCAAAGGACTCTATGCTGATGGTGGCGGTTGCAAACTTTGTGG
CTCTGTGGAACATTTAAAGAAAGATTGCCCTGAAAGTCAGAAATTCAGAGCGAATGGTCAC
AGTTGGTTCGCTGGGCAAGGGAAATGAGTGCAGACTATGAAGAAATTTGGATGTACCTAA
ACCGCAAAAACCCAAAACAAAAATACCTAAAGTTGTTAATTTTGATAACAGCTAGCACT
ATCATGAGTTACTACCTATTGTTACTTTCTAAACCAGGCCCGCTTCACGAGTTAGAGTT
GAGCTCCCCTGTAGCCAGGACTATGCTGTAGATATCAGTATGATCTGGGTGTGGCCAAAA
ACAATTTTCTTTATTCTGTCTATCAAATAGTACTTCTACCACTGTTTGGAGAAAATTGAA
GAAAAGAATAAGATGATTAAATGAATTCTCTAAAGAACATATTTTAAGAGACAGAACTT
AGACATAACCAAGTAGTTGTATACCTGATTGTAACAATCATCTTTTATAAAAGCAAAATT
ATGCATAAATGTAACACCCCTAAACAAAGGGAACAAATGAATTTTGGTTATATTGTT
TTCTTGGCTCAAAGTATCAATTTATTATAATAGAAATTTATAAGTTGCTAGAAAGTC
TATATTTTAAAGTACGGAATAAAATATGTTTCTATTAAGCAAAAAAAAAAAAAA
  
```

Restriction Sites: Please inquire


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ACCN:</b>                  | NM_032280                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_032280.1, NP_115656.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq Size:</b>           | 4635 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 816 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 84240                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <u>Q8N567</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 5q14.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Domains:</b>               | zf-CCHC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Gene Summary:</b>          | <p>May down-regulate transcription mediated by NF-kappa-B and the serum response element. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) is the longest transcript. Variants 1-3 encode the same protein.</p>                                                                                                                                                                                                                                                                              |