

Product datasheet for SC321828

DNAJC25 (NM_001015882) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNAJC25 (NM_001015882) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNAJC25
Synonyms:	bA16L21.2.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_001015882.2

```

GTGGCAGAGCCGCCAGCGAGGCTGGGGATGGGGGCGCCGCTGCTCTCCCGGCTGGGGA
GCCGGGGCTGCCGGCCGGCGCTGGTGGATGCTGCTGGCGCCCTGCTGCCGGCGCTGCTG
CTGGTGCGGCCCCCGGGGGCCCTGGTGGAGGGGCTCTACTGCGGCACGCGGGACTGCTAC
GAGGTGCTGGGCGTGAGCCGCTCGGCGGGCAAGGCGGAGATCGCGCGGGCCTACGCCAG
CTGGCCCCGGCGCTACCACCCTGACCGCTACCGGCCCCAGCCCGGAGACGAGGGCCCCGGG
CGGACGCCGAGAGCGCCGAGGAGGCTTTCCTGCTGGTGGCAACCGCTACGAGACACTC
AAGGATGAAGAAACACGAAAAGATTATGATTACATGCTGGATCATCCAGAAGAGTACTAC
AGCCATTACTACCACTACTATAGCAGGCGCTTGGCCCTAAGGTGGATGTTAGAGTAGTG
ATTTTGGTCAGCGTGTGTGCTATTTTCGGTGTTCAGTTTTTCAGCTGGTGAATAGCTAC
AATAAGGCAATCAGCTACCTAGCCACAGTGCCCAAGTACCGTATCCAAGCTACAGAGATT
GCCAAGCAGCAGGACTGCTCAAAAAAGCCAAAGAAAAAGGCAAAAAACAAAAAGTCCAAA
GAAGAAATTCGTGACGAGGAGGAGAATCATATAAGAACATTATAAAAGTAAATAGAT
ATAAAGGGGGGCTATCAGAAACCCCAATCTGTGATCTTCTCCTGTTCAAATTATCTTA
GCTCCTTTTACCTATGCTCATATATAGTTTGGTATTGTGCGGTGGATCTATAATTTAAC
ATCAAAGGCAAGAATATGGAGAGGAAGAGAGATTATACATTATACGTAATCTATGAAG
ATGTCAAAGTCTCAATTTGATAGTCTAGAAGATCATCAGAAAGAACTTTTCTTAAACGA
GAGCTCTGGATCAAGGAGAATTATGAGGTCTACAAGCAAGAACAAGAGGAAGAATTAAAG
AAAAAGTTGGCAAATGACCCAGATGGAAGAGATACAGGAGATGGATGAAGAATGAAGGG
CCTGGGCGGTAAACATTTGTCGATGACTGAAGATTGATGGAATGCTACTATGCCAAACCT
TAATTGTGATATTATTTTCATACTGAATTATTTTAGAAATGTATCAATTGACTGCTGCT
CAGCAGTAACTAAAATTCCTCAAGTATTTGATTAACAGAATAATGTCAAAATTTAAACC
TTCCCTTAAACTTTTATACATAAGACATTTATGATTGTTCAATTTTATAATCTATTTGT
GGATTTTGTATAAAGATTTACATGAAGATTTATTAGTTGCCATTTAAAAATTTTATATG
TTTAGTTAAAGATTTGACATGAAAATTTATTAGATACCATTTAAAAATTTAAAAAAA
AAAAAA
  
```


[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001015882
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001015882.2, NP_001015882.2</u>
RefSeq Size:	2268 bp
RefSeq ORF:	1083 bp
Locus ID:	548645
UniProt ID:	<u>Q9H1X3</u>
Cytogenetics:	9q31.3
Protein Families:	Transmembrane