

Product datasheet for SC321961

SPIN2A (NM_019003) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPIN2A (NM_019003) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPIN2A
Synonyms:	dj323P24.1; DXF34; SPIN2; TDRD25
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC321961

```

GGAGCCCAGGAGGGGCCCAAGAAAAGAAAGGCAGATCTGCTCTGTGATTAGAATGGAAG
ACATCACTACTAAGAGGAACTAACACTGTAGTCTTGTGCAGAAAAAGGCAGAGACTAC
AGCCACTGCATCAGGGACTGGTGGCAGCCCAACGTGGCTGAGGGGTGGCCGCCACCAACA
GGGAGGTAGACGCCAGCGAGCGAAGCGCTGGGGACACCGCTGCAGTGCTCTTTGCCTGT
CTCGCCGCGCTCTACACACGCTCCTCATCTATGGATCCTACCCTCCCCGCTCTTGCAAG
CAGTCACTCAGCGAGTCGCTGTGACCCTCCTTCGCCACAGGCTCACAGCAGAGGCAGC
ATCGAGCAGGCATGAAGACCCCCAACGCACAGGAAGCCGAAGGGCAACAAACCAGGGCAG
CTGCAGGACGGGCCACTGGGTCTGCAACATGACAAAGAAAAAGTCTCCCAAAAGAAGC
AGAGAGGCAGACCTTCATCCCAACCCTGCAGGAACATCGTGGGCTGCAGAATTTCTCATG
GATGGAAGGAAGGAGATGAGCCCATCACGCAGTGGAAAGGAACCGTTCTGGATCAGGTGC
CTATAAATCCCTCTCTTTATCTGGTGAATATGATGGAATTGACTGTGTCTATGGACTGG
AACTTCACAGAGATGAAAGGGTTTTGTCTCTTAAATTTCTTTCTGACAGGGTGGCATCAT
CTCACATTAGTGATGCCAACCTTGCAAATACCATAATTGGCAAAGCAGTGGAAACACATGT
TTGAGGGAGAGCATGGTTCTAAGGATGAATGGAGGGGGATGGTCTTAGCTCAAGCACCTA
TCATGAAAGCCTGGTTTTATATTACCTATGAGAAAGATCCTGTCTTGATCATGTACCAGC
TTCTAGATGATTATAAGGAAGGTGACCTCCGCATCATGCCAGAATCCAGTGAGTCTCCTC
CAACAGAGAGGGAGCCAGGAGGAGTTGTAGATGGCCTAATAGGTAAGCATGTGGAATATA
CCAAAGAAGATGGCTCCAAAAGGATCGGCATGGTCATTACCAAGTGGAACCAAAACCCT
CTGTGTATTTTCATCAAGTTTGATGATGATTTCCATATCTATGTCTACGATTTGGTGAAAA
AGTCCTAACTGTTAGGGTAAAAATTTGGCACATGTGTGGAACAAATGTATAATTTGTAGA
CATGCAAAAAATGTTGCCTTTTCAGTGTACTGAAAGCTTATGGAATCCCTGATAACTAAAC
ATCTTTGCCAGCATTAAGTGTGTTTGTCTCTAAAAAATACAAATTTGTGAATACATGAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites: Please inquire


[View online »](#)

ACCN:	NM_019003
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019003.3</u> , <u>NP_061876.3</u>
RefSeq Size:	1351 bp
RefSeq ORF:	777 bp
Locus ID:	54466
UniProt ID:	<u>Q99865</u>
Cytogenetics:	Xp11.21
Domains:	Spin-Ssty
Gene Summary:	This gene encodes one of three members of the DXF34 gene family, located in a 100-kb region of chromosome Xp11.21. [provided by RefSeq, Nov 2009]