

Product datasheet for **SC321973**

ZNF784 (NM_203374) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF784 (NM_203374) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF784
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for SC321973
 GGCTTGGGCACAAGATGGCCGCTGCGCGCCAGAGGCCAGAGTCGGAGCTCACCAGCTC
 CGGAGTCGCGATCCCAGGAGCCACTGGACCTGGTCTGGTGCCTGATGACTGCCGGCCTG
 GCACACCCCCGAGTGACCTCATCGAGATCCAGGTGGTGAAGGTGACGGACACCACGCTGG
 TCCCTGAGCCCCCGGAGCCAGGTTCTTTCCACTGTGCCTTGTGCCCTGCTGCCTTCCGGC
 TGGTTTCCGAGCTGCTGTTCCACGAACATGGCCACTTGGCTGGGGCCGAGGGAGGCGGGC
 AGGGTGGGGACCCGAGCCGGTGCCACGTGTGCGGCCACAGCTGCCCGGGCCCCGCGAGCC
 TGCAGCGCGCACTACAGCTTGCACACGGGGAGCGGCCCTACCGCTGCGCGCTCTGCCCCC
 GCGCCTTCAAGGCCTTGGCGCCCCCTGCTCCGGCACCAGCACCGGCACGGGTGGAGCCGG
 GCACCTCTCGGAGGCCTCCGGACACAGCGGCCGTTGCAGAACAGAGGCCCGGGTGGCCC
 CGGAGAGGGCGGAGGTGGTATGGCGGCGGCGGGCGGCGCAGCGGTGGGAAGCCTT
 TTGCTGACAGTTCTGCGCCAAGCCCTCCGCCGCTCTCAGACATGCGAGACCACGAGC
 GCGTGACACTGGCGAGCGGCCATACATTGCGGCATCTGCGGAAGGGCTTACCCAGT
 CCTCGGTGCTTAGCGGCCACGCCGCATCCACTGGCGAGCGCCGTTCCGCTGCACGC
 TCTGCGACCGCACCTTCAACAACCTCCAACCTCCGCAAGCACCAGCGCACCCACTTCC
 ACGGGCCCGGGGCGGGGCTGGGAGACTCTGGAGGCCAGCTGGGCTCGTCCGCGGCTGAGG
 GGTCCGGGAGCGGGTGTGGGGTAGGGGACCTGCGGAGGAGGGGCGGGGGAGACCGCGA
 AGGTGAAGGTGGAGGCCGACAGTAGGCTGCAGACCAGGGACATGGGTAGTTAAGGCGG
 AGGTACGAGGCCGACCCGGGAGGGGACGGGTGGGGAACGAGGAGGGGACAGAAG
 AGGGGAGATTGGGGGAGAGATGACAGCTGCAGGGATGGTTGTGAGCCGCTAGTGATGGAG
 AGCAGAGGGAGAGGGCCAGGCTCCAGACTCCCACAGCCACACAGCACCTCTGCCAGGC
 CTAGGAGAGGACAGGTGCAGCTCTTGAGCTCTGCGGGTGTGCGGCCAAAGGCAAGGCC
 ACGGGCTGGATGTCACTTCCCCGACTGTCTTTGGTTGGCTTGTCTTGTGCAAGACCCA
 GCCTGTACGACAGAGCCTGGGCACTTCAGAGGAGGAGCCAGGATCGAATGGAAGGGGGG
 AATTGGGGTCCACCATAGTCTTCTGCTCTGGTCTCCACGGGTGGACCAGGATGGAGTCT
 CCTGCCTACCTCACTGCATTGCACTGGACCTGGGATGCCTATCCACCCTCTGGCAGAAGA
 CACTCACCAGGTTATCTGTGAAGAGACTCTGGGATCCCATCACCTCAAAGCCAGAGGGTC
 CCCAAGTCACCGCTGAGAGCACTTGAGCCTCAAGGATGTAAGCCTGACCATAGGATCTTG
 ACTCCAACAGCGCAACCCCCACCCCATTTGTGGTCCGTCCTTAACCCATCCACTCTTCT
 TCGGAGGCAACTGAGAACACATAAAGCAAGCAGCTACCTAGCATCCCCCTCTAAAGCTT
 TAGACTCAGAGCCAGGGTCCCCACAAGCCTCAAGGTAGCCTCAGGTTTCTCTAATTTT
 CTCCACTCCAGTTTGAAGCAAACAGCTTACTGCCTAGTCCCGCAATCCCAAGGGCGG
 GCTGGCTGATGGCAGCATGGTGGCTGGCTGGGTGTGGAGTGAAGAGTCACTGTGGTG
 GGGGCGAGAGGAGGACTTGGGAGCTGGAGGTGTGACACCTTCAGTTCTGTTCTATTAA
 GGACCTTCTGAAGGGAAA

Restriction Sites: Please inquire

ACCN: NM_203374

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_203374.1, NP_976308.1</u>
RefSeq Size:	2006 bp
RefSeq ORF:	972 bp
Locus ID:	147808
UniProt ID:	<u>Q8NCA9</u>
Cytogenetics:	19q13.42
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]