

Product datasheet for SC322015

ACTH (POMC) (NM_000939) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACTH (POMC) (NM_000939) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACTH
Synonyms:	ACTH; CLIP; LPH; MSH; NPP; OBAIRH; POC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for SC322015
 GCGTCCCCGCCCTCAGAGAGCAGCCTCCCGAGACAGGGTCCCACCAATCTTGTTTGCTT
 CTGCAGAGCCTCAGCCTGCCTGGAAGATGCCGAGATCGTGCTGCAGCCGCTCGGGGGCCC
 TGTTGCTGGCCTTGCTGCTTCAGGCCTCCATGGAAGTGCCTGGCTGGTGCCTGGAGAGCA
 GCCAGTGTGAGGACCTCACCCAGGAAAGCAACCTGCTGGAGTGCATCCGGGCCTGCAAGC
 CCGACCTCTCGGCCGAGACTCCCATGTTCCCGGAAATGGCGACGAGCAGCCTCTGACCG
 AGAACCCCGGAAGTACGTATGGGCCACTTCCGCTGGGACCGATTTCGGCCGCGCAACA
 GCAGCAGCAGCGGCAGCAGCGGCAGGGCAGAAGCGCGAGGACGTCTCAGCGGGCGAAG
 ACTGCGGCCCGCTGCCTGAGGGCGGCCCGAGCCCGCAGCGATGGTGCCAAGCCGGGCC
 CGCGCGAGGGCAAGCGCTCCTACTCCATGGAGCACTTCCGCTGGGCAAGCCGGTGGCA
 AGAAGCGCGCCAGTGAAGGTGTACCCTAACCGCGCCGAGGACGAGTCGGCCGAGGCCT
 TCCCCCTGGAGTTCAAGAGGGAGCTGACTGGCCAGCGACTCCGGGAGGGAGATGGCCCCG
 ACGGCCCTGCCGATGACGGCGCAGGGGCCAGGCCGACCTGGAGCACAGCCTGCTGGTGG
 CGGCCGAGAAGAAGGACGAGGGCCCTACAGGATGGAGCACTTCCGCTGGGGCAGCCCGC
 CCAAGGACAAGCGCTACGGCGGTTTCATGACCTCCGAGAAGAGCCAGACGCCCTGGTGA
 CGCTGTTCAAAAACGCCATCATCAAGAACGCCTACAAGAAGGGCGAGTGAGGGCACAGCG
 GGGCCCCAGGGTACCCTCCCCAGGAGGTCGACCCCAAAGCCCTTGCTCTCCCTGCC
 CTGCTGCCGCTCCAGCCTGGGGGTCGTGGCAGATAATCAGCCTCTTAAAGCTGCCTG
 TAGTTAGGAAATAAAACCTTTCAATTTCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_000939



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_000939.2</u> , <u>NP_000930.1</u>
RefSeq Size:	1245 bp
RefSeq ORF:	804 bp
Locus ID:	5443
UniProt ID:	<u>P01189</u>
Cytogenetics:	2p23.3
Protein Families:	Druggable Genome
Protein Pathways:	Adipocytokine signaling pathway, Melanogenesis

Gene Summary:

This gene encodes a preproprotein that undergoes extensive, tissue-specific, post-translational processing via cleavage by subtilisin-like enzymes known as prohormone convertases. There are eight potential cleavage sites within the preproprotein and, depending on tissue type and the available convertases, processing may yield as many as ten biologically active peptides involved in diverse cellular functions. The encoded protein is synthesized mainly in corticotroph cells of the anterior pituitary where four cleavage sites are used; adrenocorticotrophin, essential for normal steroidogenesis and the maintenance of normal adrenal weight, and lipotropin beta are the major end products. In other tissues, including the hypothalamus, placenta, and epithelium, all cleavage sites may be used, giving rise to peptides with roles in pain and energy homeostasis, melanocyte stimulation, and immune modulation. These include several distinct melanotropins, lipotropins, and endorphins that are contained within the adrenocorticotrophin and beta-lipotropin peptides. The antimicrobial melanotropin alpha peptide exhibits antibacterial and antifungal activity. Mutations in this gene have been associated with early onset obesity, adrenal insufficiency, and red hair pigmentation. Alternatively spliced transcript variants encoding the same protein have been described. [provided by RefSeq, Jan 2016]

Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 3. Variants 1-4 encode the same protein.