

Product datasheet for SC322030

ANKRD37 (NM_181726) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD37 (NM_181726) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD37
Synonyms:	Lrp2bp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322030 <pre> GGATTAGGCTGCTCGGGCGTAACCGGAGCTGGGGCGCGGCTGCCAAGGGCCGGCCCGGAA GTCCCAGCGGTCTTTAAATTCTCCCGTGCTAGGGCCAGCCTGCGCATTCTTACCTGTCGG GGTGCGGCGAGTGTCTCACCTCTCTGCACTTCCAAGGACTCTTGTCTCATCTGCCTTAGCGG GGAAATGCTGTTGCTGGATTGCAACCCCGAGGTGGATGGTCTGAAGCATTGCTGGAGAC AGGGGCGCTCGGTCAACGCACCCCGGATCCCTGCAAGCAGTCGCCTGTCCACTTAGCCGC AGGAAGCGGCTTGCTTGCTTTCTCTCTGGCAGCTGCAAACGGGCGCTGACCTCAACCA GCAGGATGTTTTAGGAGAAGCTCCACTACACAAGGCAGCAAAAGTTGGAAGCCTGGAGTG CCTAAGCCTGCTTGAGCCAGTGATGCCAAATTGATTTATGTAATAAGAACGGGCAAAC AGCTGAAGATCTCGCTTGGTCATGTGGATTTCCAGACTGTCCAAGTTTCTTACAACAAT TAAATGTATGCAGACAATAAAAGCAAGTGAACACCCTGACAGGAATGATTGTGTTGCCGT GCTCAGACAGAAACGGAGTCTCGGAAGTGTAGAAAATACCAAGTGGGAAAAGGAAGTGCTG ATGTCACGTGGGTATGAAGAAGTCTGAAGAAGCCTTCATTTTCATGCAAATCTATAAGC TCCTGCTTTTGGCTTTACCATATGTTGTGTCTAATCTCCTTCTGAGAAGGACGAAAACT TTCTTCCAAGTGAAGATCCATTTAAGAACACATGTATTTACATGCCTATAATATGCTGGT TGTGTATGCTTTGTCTTTAAGTTATTAAGGAACGTCTAAAAAAAAAAAAAAAAAAAAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_181726
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_181726.1, NP_859077.1</u>
RefSeq Size:	941 bp
RefSeq ORF:	477 bp
Locus ID:	353322
UniProt ID:	<u>Q7Z713</u>
Cytogenetics:	4q35.1
Protein Families:	Druggable Genome