

Product datasheet for **SC322053**

VWC2 (NM_198570) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VWC2 (NM_198570) Human Untagged Clone
Tag:	Tag Free
Symbol:	VWC2
Synonyms:	PSST739; UNQ739
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_198570, the custom clone sequence may differ by one or more nucleotides

```
ATGCCAGCTCCACTGCGATGGCAGTTGGCGCGCTCTCCAGTTCCTCCTGGTCACCTGCTGCCTGATGG
TGGCTCTGTGCAGTCCGAGCATCCCGCTGGAGAAGCTGGCCAGGCACCAGAGCAGCCGGGCCAGGAGAA
GCGTGAGCACGCCTCTCGGGACGGCCCGGGCGGGTGAACGAGCTCGGGCGCCCGGCGAGGGACGAGGGC
GGCAGCGGCCGGGACTGGAAGAGCAAGAGCGGCCGTGGGCTCGCCGGCCGTGAGCCGTGGAGCAAGCTGA
AGCAGGCCTGGGTCTCCAGGGCGGGGGCGCAAGGCCGGGGATCTGCAGGTCGGCCCCCGGGGGACAC
CCCGCAGGCGGAAGCCCTGGCCGACGCCCCAGGACGCGATTGGCCCGAACTCGGCCACGCCCGAG
CCACCCGAGGAGTACGTGTACCCGGACTACCGTGGCAAGGGCTGCGTGGACGAGAGCGGCTTCGTGTACG
CGATCGGGGAGAAGTTCGCGCCGGGCCCTCGGCCTGCCCGTGCCTGTGCACCGAGGAGGGCCCGCTGTG
CGCGCAGCCCGAGTGCCCGAGGCTGCACCCGCGCTGCATCCACGTGACACGAGCCAGTGCCTGCCCGCAG
TGCAAGGAGAGGAAGAATACTGCGAGTTCGGGGCAAGACCTATCAGACTTTGGAGGAGTTCGTGGTGT
CTCCATGCGAGAGGTGTCGCTGTGAAGCCAACGGTGAGGTGCTATGCACAGTGTGACGCTGTCCCGAGAC
GGAGTGTGTGGACCCTGTGTACGAGCCTGATCAGTGTGTCCCATCTGCAAAAATGGTCCAAACTGCTTT
GCAGAAACCGCGGTGATCCCTGCTGGCAGAGAAGTGAAGACTGACGAGTGCACCATATGCCACTGTACTT
ATGAGGAAGGCACATGGAGAATCGAGCGGCAGGCCATGTGCACGAGACATGAATGCAGGCAATGTAG
```

Restriction Sites: Please inquire

ACCN: NM_198570



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198570.1</u> , <u>NP_940972.1</u>
RefSeq Size:	1942 bp
RefSeq ORF:	978 bp
Locus ID:	375567
UniProt ID:	<u>Q2TAL6</u>
Cytogenetics:	7p12.2
Gene Summary:	This gene encodes a secreted bone morphogenic protein antagonist. The encoded protein is possibly involved in neural function and development and may have a role in cell adhesion. [provided by RefSeq, Oct 2009]